

U.S. ENVIRONMENTAL PROTECTION AGENCY

Office of Pesticide Programs
Antimicrobials Division (7510P)
1200 Pennsylvania Ave., N.W.
Washington, D.C. 20460

NOTICE OF PESTICIDE:

X Registration
Reregistration
(under FIFRA, as amended)

EΡΔ	Reg	Number:
LIA	IXCZ.	rumber.

Date of Issuance:

96706-1

12/11/2020

Term of Issuance:

Conditional

Name of Pesticide Product:

MBS MedTech Germicidal Disposable Wipes

Name and Address of Registrant (include ZIP Code):

Medical BioWaste Solutions, Inc. 25917 Pala, Suite 101 Mission Viejo, CA 92691

Note: Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Antimicrobials Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby registered under the Federal Insecticide, Fungicide and Rodenticide Act.

Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA section 3(c)(7)(A). You must comply with the following conditions:

1. Submit and/or cite all data required for registration/registration/registration review of your product under FIFRA when the Agency requires all registrants of similar products to submit such data.

Signature of Approving Official:

Date:

Date:

Date:

Date:

Date:

Date:

12/11/202020

Jacqueline Hardy
Product Manager, Team 34
Regulatory Management Branch 2
Antimicrobials Division (7510P)
Office of Pesticide Programs

EPA Form 8570-6

- 2. You are required to comply with the data requirements described in the DCI identified below:
 - a. Alkyl (68% C₁₂, 32% C₁₄) dimethyl ethylbenzyl ammonium chloride GDCI-30906-1748
 - b. Alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈) dimethyl ethylbenzyl ammonium chloride GDCI-30881-1739

You must comply with all of the data requirements within the established deadlines. If you have questions about the Generic DCI listed above, you may contact the Reevaluation Team Leader (Team 36): http://www2.epa.gov/pesticide-contacts/contacts-office-pesticide-programs-antimicrobial-division

- 3. Make the following label changes before you release the product for shipment:
 - Revise the EPA Registration Number to read, "EPA Reg. No. 96706-1."
- 4. Submit one copy of the final printed label for the record before you release the product for shipment.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) lists examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Assurance.

If you fail to satisfy these data requirements, EPA will consider appropriate regulatory action including, among other things, cancellation under FIFRA section 6(e). Your release for shipment of the product constitutes acceptance of these conditions. A stamped copy of the label is enclosed for your records. Please also note that the record for this product currently contains the following CSFs:

Basic CSF dated 9/14/2020

If you have any questions, please contact Terria Northern by phone at 703.347.0265, or via email at Northern. Terria@epa.gov.

Enclosure Stamped Label

MBS MedTech GERMICIDAL* DISPOSABLE WIPES

ACTIVE INGREDIENTS:

n-Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzyl ammonium chlorides	0.25%
n-Alkyl (60% C ₁₄ , 30% C ₁₆ , 5%C ₁₂ , 5%C ₁₈) dimethyl benzyl ammonium chlorides	0.25%
Isopropyl Alcohol	55.00%
OTHER INGREDIENTS	44.50%
TOTAL	100.00%

[Does not include the weight of the [cloth] [towelette] [or] [wipe]]

KEEP OUT OF REACH OF CHILDREN WARNING

See [back] [outer carton] [side] [container[s] inside] [panel] for additional precautionary statements

[Contains:] [NET CONTENTS] [and/or Net Wt.] _____oz (g) [or lbs. (kg)]

[# OF CLOTHS] [SIZE of each cloth]

EPA Reg. No.: 96706-R EPA Est. No.: 96706-TX-1

Manufactured by:
Medical BioWaste Solutions
Medical Technologies
25971 Pala, Suite 101, Mission Viejo, CA 92961

ACCEPTED

12/11/2020

Under the Federal Insecticide, Fungicide and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No. 00700

96706-1

First Aid:

Call a poison control center or doctor [or healthcare professional] for treatment advice. Have the product container or label with you when calling a poison control center or doctor [or healthcare professional], or going for treatment. [For additional information in case of emergency call toll free [phone number to be inserted]].

If in eyes:

- -Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin:
- -Take off contaminated clothing.
- -Rinse skin immediately with plenty of water for 15-20 minutes.

Precautionary Statements

Hazards to Humans and Domestic Animals

Warning. Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash hands thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using [the] [toilet] [restroom]. Remove and wash contaminated clothing before reuse.

Physical or chemical hazard

Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

To Dispense WIPES:

Remove lid and discard inner seal from container. Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the container lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep lid closed to prevent moisture loss.

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

[See directions for use on the container[s] inside.]

TO [CLEAN] [AND] [,] DISINFECT [AND DEODORIZE] HARD, NONPOROUS SURFACES: If present, use a wipe to remove visible soil prior to disinfecting. Unfold [a] [clean] wipe and thoroughly wet surface. Allow surface to remain wet for 2 minutes. Repeated use of wipe may be needed [Use enough wipes] to ensure that the surface remains wet for 2 minutes. Let air dry. A potable water rinse is required for food contact surfaces.

[The following set of directions are optional and for hospital/medical facility use only:

Ultrasound transducers and probes: For use on external surfaces of ultrasound transducers and probes only.

To disinfect: Preclean gel from external surface of ultrasound transducer or probe to be disinfected by using a clean washcloth. Then thoroughly wet the probe with a wipe. Allow surface to remain wet for 2 minutes. Repeated use of wipe may be needed [Use enough wipes] to ensure that the surface remains wet for 2 minutes. Let air dry.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal. [Statement not required for residential or household use products.]

Storage: Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss.

Towelette Disposal: Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling. If recycling is not available, put in trash collection.

Areas of Use:

Hospital and Healthcare Settings: Ambulances, Ambulatory Surgical Centers (ASC), Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctors' Offices, Donation Centers [blood][plasma][semen][milk][apheresis], Emergency Medical Settings, Emergency Rooms, E.R., Emergency Vehicles, Eye Surgical Centers, Healthcare Settings, Health Care Facilities, Home Health Care, Hospices, Hospitals, ICU, Isolation Areas, Laboratory, Laundry rooms, Long Term Care Centers, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating Rooms (OR), Ophthalmic Offices, Orthopedics, Out-Patient Surgical Centers (OPSC), Patient Care Areas, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physicians' offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory therapy, Skilled Nursing Facility, Surgery, Surgery Intensive Care Unit (SICU), Surgical Centers, Transport vehicles and X-Ray.

Critical Care Areas: CCU, Emergency Rooms, E.R., ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Unit (SICU).

and

Hospital, Healthcare, and Critical Care Use Sites: May be used on hard, nonporous surfaces of: ambulance equipment surfaces; bed railings; bladder scan; blood glucose meters; cabinets; cardiac gym equipment; carts; chairs; counters; defibrillators; dental chairs; dental countertops; dental unit instrument trays; dopplers; electronic thermometers; endodontic equipment such as apex locators; exam tables; gurneys; infant incubators [interior and exterior surfaces of]; infant warmers [interior and exterior surfaces of]; isolettes; IV poles; IV pumps; laboratory equipment and surfaces; loupes; operating room tables and lights; operatory light switches; Ophthalmoscopes; Otoscopes; oxygen hoods; physical therapy (PT) equipment surfaces; patient lift equipment; pulse oximeter; radiology equipment; slit lamps; spine back boards; stethoscopes; stretchers; stools; tables; telemetry equipment; telephones; toilet seats; tympanic thermometers; and hard, nonporous outside surfaces of: amalgamators and dental curing lights;

anesthesia machines and respiratory therapy equipment; apheresis machines, diagnostic equipment, dialysis machines; patient monitoring equipment, patient support and delivery equipment, pulp testers and motors; toilets; ultrasound monitors; and ultrasound transducers and probes.

Grocery Stores *or* **Shopping Centers** *or* **Supermarkets** *or* **Delicatessens:** Hard nonporous surfaces of grocery cart handles and grocery cart child seats.

Nursery or Daycare Centers: Hard nonporous surfaces of bassinets [interior and exterior surfaces of], changing tables, cribs, diaper changing stations, diaper pails, infant/child care equipment surfaces, trash cans, toys and other hard, nonporous surfaces listed on this label.

Bathroom[s] or Public Restrooms: Hard nonporous surfaces of Basins, [bath]tubs, bathroom fixtures, cabinets, countertops, doorknobs, faucet[s], handles, hampers, handrails, shower[s] [stalls], showers, sinks, glazed tiled walls, glazed tiles, toilet [handle] [rim] [seats] [tank] [and other exterior surfaces of], urinals [exterior surfaces of] and vanity tops, and other hard, nonporous surfaces listed on this label.

Health club facilities and Gyms: Hard nonporous surfaces of gym equipment, weight-lifting surfaces and other hard, nonporous surfaces listed on this label.

Schools or Classrooms or Churches or Office Buildings or Public Facilities or Manufacturing Facilities or Libraries or Public Telephones: Hard nonporous surfaces of cabinets, chairs, computers, computer monitors, computer mouses, counters, countertops, desks, desktops, doorknobs, fax machines, filing cabinets, filing doors, filing handles, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, telephones, toys, trash cans, work stations.

Transport vehicles, school busses, public busses, fire trucks, police cars: Hard, nonporous surfaces of vehicles including steering wheel, hard plastic or metal flooring, seats, and walls.

Restaurants: Hard nonporous, nonfood contact surfaces only including chairs, display cabinets, garbage cans, booster seats, floor, sinks, bathroom surfaces, glazed tiled walls, telephones, cash registers, and doorknobs.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to [sterilization] [or] [high level disinfection].

[Symbols including asterisks, bullets or footnotes and similar may be used to designate categories of organisms.]

[ATCC numbers are not required to be listed on the product label. Organisms may be listed in a table and in any order.]

Bacteria [(2 Minute Contact Time)]

Pseudomonas aeruginosa [ATCC 15442] Salmonella enterica [ATCC 10708] Staphylococcus aureus [ATCC 6538]

<u>Viruses* [(2 Minute Contact Time)]</u>

*Human Coronavirus [ATCC VR-740] [Strain 229E]

MBS Germicidal Wipe is virucidal against Human Coronavirus in the presence of organic soil (5% blood serum).

MBS Germicidal Wipes is a premoistened, nonwoven durable wipe containing a quaternary/alcohol based solution. For use in hospitals and other critical care areas where reduction of the hazards of cross-contamination between treated surfaces is required. Suitable for use on equipment requiring alcohol-based products. Use on hard, nonporous surfaces and equipment listed on this label made of stainless steel, plastic, Formica® and glass. [Formica® is a registered trademark of The Diller Corporation.]

<u>Optional Claims and Descriptors</u> [Note: Deodorizing directions must be used when deodorizing claims (claims involving odor control) are used]

X" x Y" [wipes]

pre-saturated towelettes X" x Y"

© [insert year] [insert company name]

[Absorbs] [Absorbent] [wipe]

After use, just toss wipe(s) away

Always ready to use

An effective method to deodorize and disinfect against odor-causing organisms on hard, nonporous surfaces listed on this label.

An effective way to disinfect hard, nonporous environmental surfaces listed on this label.

A nonwoven disposable cloth for use in hospitals and other critical care areas reduction of cross contamination is required.

Antibacterial

Antimicrobial

[Antibacterial] [Antimicrobial] health club equipment wipe

A Quaternary [ammonium chloride] solution impregnated in a wiping cloth

Are pre-measured to have the recommended amount of germicidal solution, wipe after wipe.

Assures proper dilution strength, product effectiveness and standardized technique.

[Bactericidal], [and] [Virucidal*] in [only] [just] [2 [two] minutes] Bleach-free [formula]

Can help reduce the risk of cross contamination on treated, hard, non-porous surfaces

Cellulosic nonwoven disposable cloth

Compatible for use on hard, nonporous surfaces commonly found in healthcare settings listed on this label

Convenient

Deodorize hard, nonporous surfaces [Use to]

Deodorizes, and kills common household germs⁺

Disinfect external surfaces of ultrasound transducers

Disinfect the precleaned surface

Disinfecting towelettes

Disinfects [and deodorizes] in 2 [Two] minutes

Disinfects grocery cart handles *or* shopping cart handles *or* grocery cart child seats *or* shopping cart child seats

Disinfects hard, nonporous environmental surfaces listed on this label [Use to]

Disinfects hard, nonporous surface in 2 minutes

[Disinfects] hard surfaces [that people touch]

Disinfects [health club equipment] [treadmill handles] [elliptical handles][exercise bike handles][weight benches]

Disinfects hard, non-porous surfaces [in [just] two [2] minutes]

Disinfects hard, non-porous surfaces on the go

Dispenses one wipe at a time

Does not contain bleach

Does not produce mist or overspray

Easy

Easy [and convenient] way to [touch-up] [bathroom] [office] everyday

[Easy] [convenient] [dispensing] lid [closure] [cap] [package]

Easy drop-in refill[s]

[Easy] pop-up-tub

Easy start lid [package]

Easy to use [easy to refill]

Easy to use package [canister] [container] keeps [your] wipes moist

Effective against [the following micro-organisms on hard, nonporous surfaces] [list organisms]

Eliminates odors

Eliminates potential for overspray

Gets rid of -or- eliminates odors

[Helps] [kill][s] [eliminate][s] odor-causing bacteria

Helps reduce cross contamination on treated surfaces

Hospital disinfectant [disinfection] [use only when hospital use sites on label]

Includes count refill [pack] [contains] wipes

[Just] [pull] wipe and toss

Kill germs⁺ and bacteria instead of germy paper towels, sponges and rags [Use this product to]

Kills bacteria in 2 minutes

Kills odor causing bacteria

Leaves [bathroom] [office] disinfected [and smelling fresh]

No measuring

No mixing

No sponges, no rags, no buckets, no mess

Non-scratching [wipe]

Nonwoven disposable cloth

Not for in-home use

Not for sale in [insert name of state(s)]

Odor [elimination] [neutralization]

One [wipe] at a time dispensing

On-the-go disinfectant

[Package] [Insert type of package] made from [x%] recycled materials

Perfect for disinfecting exterior toilet surfaces -or- around the exterior toilet surface

Picks up [particles] [dust] [dirt] [hair] [crumbs]

Pop-up [dispensing] [top]

Portable

[Premixed], [premeasured], [premoistened]

[Quickly] eliminates [lingering] [food] odors on surfaces

Ready to Use

Recyclable [package]

Reduces cross contamination on treated hard, nonporous surfaces

Requires no mixing or measuring and are ready-to-use

Resealable [stay-moist] pouch [package] [container] inside

Suitable for use on equipment requiring alcohol-based products

Tear at perforation

The 2-minute germicidal wipe [disinfectant wipe]

Thick wipes are [bactericidal,] [and] [virucidal*] in 2 minutes

Thick, large wipe disinfects hard, non-porous surfaces

Use on hard, nonporous surfaces only: [stainless steel] [Formica®] [glass] [glass tables] [glazed tile] [Test wipe on small inconspicuous area first.] [Formica® is a registered trademark of The Diller Corporation]

Wet wipes

Where control of the hazards of cross contamination between treated hard, non-porous surfaces is of prime importance

[Wipes out] [Kills] [eliminates] 99.9% of [list organisms] in [just] [2 minutes]

[Wiping] [soft] side for disinfecting

With a bleach-free formula

Master Label Reader Notes

[* All pesticidal "germ" claims must be linked to bacteria (minimum: Salmonella and Staph) and the Human Coronavirus.

[Graphics depicting use sites listed on this label can appear on product label and labeling.]

[The following graphics and statements are optional. They may be placed anywhere on label/container.]





NOT A BABY WIPE NOT A SKIN WIPE NOT FOR USE ON SKIN FOR USE ON HARD, NONPOROUS SURFACES ON





