



**US Environmental Protection Agency
Office of Pesticide Programs**

Office of Pesticide Programs

**Proposed Procedure: EPA MLB SOP MB-19: Growing a Biofilm using the
CDC Biofilm Reactor**

Date Revised: 07/20/16

DRAFT

Growing a Biofilm using the CDC Biofilm Reactor

I. Overview

- A. This document describes the methodology intended for use for growing a *Pseudomonas aeruginosa* or *Staphylococcus aureus* biofilm in the CDC biofilm reactor for the purpose of product efficacy testing.
- B. The method is based on ASTM E2562-12, however, it has been modified to support biofilm claims for antimicrobial products tested using the Single Tube Method.
- C. Growing the organism in the reactor is partitioned into 2 phases. The biofilm is established by operating the reactor in batch phase (i.e., no flow) for 24 hours followed by 24 hours with continuous flow of growth medium to form biofilm on coupons (CSTR phase). The coupons are harvested to conduct the efficacy test. The main components of this procedure are as follows:
1. Verify reactor operating volume (once per reactor) and periodically calibrate the pump
 2. Clean and screen coupons
 3. Prepare the reactor and steam-sterilize the assembled reactor
 4. Sterilize the 20 L carboy containing 19-20 L de-ionized water
 5. Prepare the inoculum
 6. Initiate batch phase
 7. Initiate CSTR phase
 8. Harvest coupons for efficacy testing (refer to “Single Tube Method for Determining the Efficacy of Disinfectants against Bacterial Biofilm,” dated 08/09/16)

II. Data Requirements

- A. None

III. Special Apparatus and Materials

- A. Test organisms.
1. *Pseudomonas aeruginosa* (ATCC No. 15442) obtained directly from ATCC.
 2. *Staphylococcus aureus* (ATCC No. 6538) obtained directly from ATCC.
- B. *Growth media.* Various concentrations of TSB for inoculum production, batch phase, and CSTR phase.

- 35 C. *Recovery media.*
- 36 1. R2A: for *P. aeruginosa*
- 37 2. Tryptic soy agar (TSA): for *S. aureus*
- 38 D. *Calibrated micropipettes.* For performing culture transfers.
- 39 E. *Ultrasonic water bath.* For cleaning the CDC Reactor coupons. Use any bath
- 40 capable of maintaining a homogeneous sound distribution of 45 ± 5 kHz and which
- 41 has a volume large enough to accommodate 50 mL or 250 mL conical tubes.
- 42 F. *Peristaltic pump.* Example: Masterflex L/S Computerized Drive Model 755-50
- 43 with Easy-Load II pump head, model 77201-60, or equivalent. The pump head is
- 44 capable of holding tubing with inner diameter (ID) 3.1 mm and outer diameter
- 45 (OD) 3.2 mm. Use Masterflex Norprene tubing (Cole Palmer 06404-16).
- 46 G. *Magnetic stir plate.* Top plate 10.16×10.16 cm, capable of rotating at 60-125
- 47 rpm.
- 48 H. *Silicon tubing.* Two sizes of tubing: one with ID 3.1 mm and OD 3.2 mm and the
- 49 other with ID 7.9 mm and OD 9.5 mm. Both sizes must withstand sterilization.
- 50 I. *Glass flow break.* Any that will connect with tubing of ID 3.1 mm and withstand
- 51 sterilization.
- 52 J. *Clamp.* Used to hold flow break, extension clamp with 0.5 cm minimum grip
- 53 size.
- 54 K. *Clamp stand.* Height no less than 76.2 cm, used with clamp to suspend glass flow
- 55 break vertically and stabilize tubing above reactor.
- 56 L. *Reactor Components.* For schematic illustrations of CDC reactor components,
- 57 see ref. VII.A. Reactor components available from BioSurface Technologies, Inc.
- 58 1. *Berzelius borosilicate glass tall beaker.* 1000 mL without pour spout,
- 59 9.5 ± 0.5 cm diameter. Barbed outlet spout added at 400 ± 50 mL mark.
- 60 Spout angled to $30\text{-}45^\circ$ to ensure drainage. Spout to accommodate
- 61 flexible tubing with an ID of 8-11 mm. NOTE: The rods and baffle
- 62 described in III.L.3 and III.L.6, respectively, will displace approximately
- 63 50 mL of liquid when the system is completely assembled. Therefore, an
- 64 outlet spout at the 400 mL mark will result in a reactor operating volume
- 65 of approximately 350 mL.
- 66 2. *Reactor top.* UHMW (ultra-high molecular weight) polyethylene top
- 67 (10.1 cm diameter tapering to 8.33 cm) equipped with 3 holes
- 68 accommodating 6-8 cm long pieces of stainless steel or other rigid
- 69 autoclavable tubing with OD of 5-8 mm for media inlet, air exchange and

- 70 inoculation port. Center hole, 1.27 cm diameter, to accommodate the
71 glass rod used to support the baffle assembly. Eight rod holes, 1.905 cm
72 diameter, notched to accommodate stainless steel rod alignment spike
73 (0.236 cm OD).
- 74 3. *Polypropylene rods.* Eight polypropylene rods, 21.08 cm long, machined
75 to hold three coupons at the immersed end. 316 stainless steel set screws
76 imbedded in side to hold coupons in place. Rods to fit into holes in
77 reactor top and lock into preformed notches.
- 78 4. *Coupons.* Twenty-four cylindrical coupons (i.e., borosilicate glass) with a
79 diameter of 1.27 ± 0.013 cm, thickness of approximately 3.0 mm.
- 80 5. *Small Allen wrench.* For loosening set screws.
- 81 6. *Stir blade assembly (baffled stir bar).* Teflon blade (5.61 cm) fitted into
82 cylindrical Teflon holder (8.13 cm) and held in place with a magnetic stir
83 bar (2.54 cm). Teflon holder fits onto a glass rod (15.8 cm), fitted into the
84 reactor top. The glass rod is held in place with a Swagelock fitting and
85 acts as a support for the moving blade assembly.
- 86 M. *Carboys.* Two 20 L autoclavable carboys, one used for waste and one for the
87 growth medium.
- 88 1. *Carboy lids.* Two carboy lids. One carboy lid with at least 3 barbed
89 fittings to accommodate tubing ID 3.1 mm (one for the growth medium,
90 one for bacterial air vent (filter), and one for injecting TSB medium
91 concentrate). One carboy lid with at least two 1 cm holes bored in the
92 same fashion (one for effluent waste and one for bacterial air vent).
- 93 N. *Bacterial air vent (filter).* Autoclavable 0.2 μm pore size, spliced into tubing on
94 waste carboy, carboy with growth medium, and reactor top; recommended
95 diameter 37 mm.
- 96 O. *Detergent.* Micro-90 Concentrated Cleaning Solution for Critical Cleaning;
97 International Products Corporation. For cleaning coupons and reactor parts.

98 **IV. Procedure and Analysis**

99 **A. Coupon preparation**

- 100 1. Coupons may be used repeatedly with proper cleaning and screening
101 between each use (see below). Check each coupon under 20X
102 magnification for scratches, chips, other damage or accumulated debris
103 before each use. Discard those with visible damage to surface topography.
- 104 i. After use in the reactor, place contaminated coupons in an

- 105 appropriate vessel, cover with liquid, and autoclave with the other
106 parts of the contaminated reactor system for 30 min.
- 107 ii. For glass coupons, sonicate coupons individually (e.g., in plastic
108 50 mL conical tubes or in a 24 well plate) for 5 min in a 1:100
109 dilution of detergent and tap water. Cover the coupons completely
110 with soapy water. Process coupons individually to minimize
111 damage to the coupons.
- 112 iii. Rinse coupons with reagent grade water and sonicate for
113 approximately 1 min in reagent grade water. Repeat rinsing and
114 sonication with reagent grade water until no soap is left on the
115 coupons.
- 116 iv. Once the coupons are clean, prevent oils and other residue from
117 contaminating the surface. Store screened and cleaned coupons in
118 a Petri dish.
- 119 B. Preparation of reactor
- 120 1. Place a cleaned and screened coupon into each hole in the reactor rods,
121 leaving the coupon flush with the inside rod surface. Tighten the set
122 screw.
- 123 2. Place rods loosely into reactor top (not yet fitted into notches).
- 124 3. Invert the reactor top and place baffled stir bar onto glass rod positioned in
125 the center of the reactor top.
- 126 4. Place assembled top into the reactor beaker.
- 127 5. Connect the bacterial air vent by fitting the vent to a small section of
128 appropriately sized tubing and attach to one of the rigid tubes on the
129 reactor top.
- 130 6. Splice the glass flow break into the growth medium tubing line near the
131 reactor top.
- 132 C. Sterilization of reactor system
- 133 1. Cover the end of the injection ports, the growth medium tubing that
134 connects to the growth medium carboy, and the overflow (waste) tubing
135 with aluminum foil. Cover any extra openings on the reactor top with
136 aluminum foil or plastic caps.
- 137 2. Steam-sterilize the empty reactor system for 20 min.
- 138 D. Culture preparation
- 139 1. Refer to Attachment 1 for stock culture generation and QC.

- 140 2. For *P. aeruginosa*, defrost a single cryovial and briefly vortex to mix.
141 Add 10 μ L of the thawed frozen stock (single use) to a tube containing 10
142 mL of TSB (300 mg/L), vortex, and incubate at $36\pm 1^\circ\text{C}$ for 24 ± 2 h.
- 143 3. For *S. aureus*, defrost a single cryovial and briefly vortex to mix. Add 10
144 μ L of the thawed frozen stock (single use) to a tube containing 10 mL of
145 TSB (30 g/L), vortex, and incubate at $36\pm 1^\circ\text{C}$ for 24 ± 2 h.
- 146 4. Inoculate an agar plate (e.g., TSA with 5% sheep blood, BAP) with a
147 loopful from the inoculated tube and streak for isolation. Incubate plate
148 and examine for purity.
- 149 5. If desired, verify that the concentration of the 24 ± 2 h culture is at least 10^7
150 CFU/mL.
- 151 E. Growth of biofilm in CDC reactor – batch phase
- 152 1. Make sure that the overflow (waste) line is clamped and aseptically add
153 500 mL of the batch culture medium to the cooled reactor (e.g., carefully
154 remove one rod, pour the medium into the reactor through the rod
155 opening, and re-insert the rod).
- 156 i. For *P. aeruginosa*, the batch culture medium is 300 mg/L TSB.
- 157 ii. For *S. aureus*, the batch culture medium is 3 g/L TSB.
- 158 2. Secure the rod alignment pins into the reactor top notches.
- 159 3. Place prepared reactor on a stir plate.
- 160 4. Clamp the flow break in an upright position.
- 161 5. Vortex the 10 mL tube of culture and use 1 mL to inoculate the reactor
162 (see I.D.2-3) through one of the available rigid stainless steel tubes in the
163 reactor top.
- 164 6. Turn on the magnetic stir plate.
- 165 i. For *P. aeruginosa*, the rotational speed of the baffle is 125 ± 5 rpm.
166 Run the reactor system in batch phase at room temperature (e.g.,
167 $21\pm 2^\circ\text{C}$) for 24 ± 2 h. Record ambient temperature.
- 168 ii. For *S. aureus*, the rotational speed of the baffle is 60 ± 5 rpm.
169 Incubate the reactor system in batch phase at $36\pm 1^\circ\text{C}$ for 24 ± 2 h.
- 170 F. CSTR phase for *P. aeruginosa*
- 171 1. For *P. aeruginosa* biofilm, run the CSTR phase at room temperature.

- 172 2. Add 50 mL of sterile 40 g/L TSB to 20 L sterile water (in carboy) to
173 achieve a final growth medium concentration of 100 mg/L TSB. Shake
174 the carboy to thoroughly mix the contents. Other concentrations of TSB
175 may be used to achieve a final concentration of 100 mg/L TSB.
- 176 3. Aseptically connect the tubing from the reactor to the carboy containing
177 the CSTR growth medium (100 mg/L TSB). Record ambient temperature.

178 G. Preparation of CSTR medium for *S. aureus*

- 179 1. For *S. aureus* biofilm, run the CSTR phase at $36\pm 1^\circ\text{C}$. Preheat the CSTR
180 phase growth medium prior to it entering the reactor in order to maintain
181 $36\pm 1^\circ\text{C}$ during CSTR phase.
- 182 2. Add 0.5 L of sterile 40 g/L TSB to 19 L sterile water (in carboy), then fill
183 to 20 L with sterile water to achieve a final growth medium concentration
184 of 1 g/L TSB. Shake the carboy to thoroughly mix the contents. Other
185 concentrations of TSB may be used to achieve a final concentration of 1
186 g/L TSB.
- 187 3. Aseptically connect the tubing from the reactor to the carboy containing
188 the CSTR growth medium. At least one day prior to CSTR phase, preheat
189 a portion of the CSTR growth medium from the carboy (assuming the
190 carboy is outside the incubator). To preheat, loosely coil ~12-15 feet of
191 tubing (ID 3.1 mm) inside the incubator with the reactor. Prime the tubing
192 with CSTR growth medium (1 g/L TSB).
- 193 4. Refer to Attachment 3 for the reactor/tubing setup.

194 H. Growth of biofilm in CDC reactor – CSTR mode

- 195 1. Pump a continuous flow of growth medium into the reactor to achieve a
196 30 ± 2 min residence time based on the reactor's operating volume (e.g.,
197 10.8 mL/min). Connect the end of the reactor drain to the waste carboy
198 and remove the clamp.
- 199 2. Allow the reactor to run in CSTR mode for 24 ± 2 h.
- 200 3. For conducting disinfectant efficacy or control coupon count evaluations,
201 proceed to "Single Tube Method for Determining the Efficacy of
202 Disinfectants against Bacterial Biofilm," dated 08/09/16.

203 V. Data Analysis and Calculations

- 204 A. None

205 VI. Attachments

- 206 A. Attachment 1: Procedures for Maintenance of Vegetative Bacterial Cultures –

- 207 Preparation of Frozen Stock Cultures
- 208 B. Attachment 2: Typical Growth Characteristics of strains of *P. aeruginosa* and *S.*
209 *aureus*
- 210 C. Attachment 3: Reactor/Tubing Setup for Generation of *S. aureus* Biofilm
- 211 **VII. References**
- 212 A. ASTM International, 2012. E-2562-12: Standard Test Method for Quantification
213 of *Pseudomonas aeruginosa* Biofilm Grown with High Shear and Continuous
214 Flow Using CDC Biofilm Reactor.
- 215 B. ASTM International, 2013. E-2871-13: Standard Test Method for Evaluating
216 Disinfectant Efficacy Against *Pseudomonas aeruginosa* Biofilm Grown in CDC
217 Biofilm Reactor Using Single Tube Method
- 218 C. Krieg, Noel R. and Holt, John G. 1984. Bergey's Manual of Systematic
219 Bacteriology Volume 1. Williams & Wilkins, Baltimore, MD. *P. aeruginosa* p.
220 164.
- 221 D. Sneath, P., Mair, N., Sharpe, M.E., and Holt, J. eds. 1986. Bergey's Manual of
222 Systematic Bacteriology Volume 2. Williams & Wilkins, Baltimore, MD. *S.*
223 *aureus* p. 1015.
- 224 E. Package Insert – Gram Stain Kit and Reagents. Becton, Dickinson and Company.
225 Part no. 882020191JAA. Revision 07/2011.
- 226 F. Package Insert – Oxidase Reagent Droppers. Becton, Dickinson and Company.
227 Part no. L001133. Revision 06/2010.
- 228 G. Package Insert – Catalase Reagent Droppers. Becton, Dickinson and Company.
229 Part no. L001237. Revision 06/2010.
- 230 H. Package Insert – Staphaurex Plus. Remel. Part no. R30950102. Revised
231 11/23/07.
- 232

233 Attachment 1

234 Procedures for Maintenance of Vegetative Bacterial Cultures – Preparation of Frozen Stock
235 Cultures

236 I. Preparation of Frozen Stock Cultures.

237 A. Initiate new stock cultures from lyophilized cultures of *Pseudomonas aeruginosa*
238 and *Staphylococcus aureus* from ATCC at least every 18 months.

239 B. Open ampule of freeze dried organism per manufacturer's instructions. Using a
240 tube containing 5-6 mL of TSB, aseptically withdraw 0.5 to 1.0 mL and rehydrate
241 the lyophilized culture. Aseptically transfer the entire rehydrated pellet back into
242 the original tube of broth. Mix thoroughly. Incubate broth culture at $36\pm 1^{\circ}\text{C}$ for
243 24 ± 2 h.

244 C. After incubation, streak a loopful of the suspension on TSA to obtain isolated
245 colonies. Incubate the plates for 18-24 h at $36\pm 1^{\circ}\text{C}$. Refer to section II for QC of
246 stock cultures.

247 D. Select 3-5 isolated colonies of the test organism and re-suspend in 1 mL of TSB.
248 For *S. aureus*, select only golden yellow colonies. Multiple phenotypes are
249 present for *P. aeruginosa* – the stock culture should be representative of all
250 phenotypes present on the streak isolation plate. Spread plate 0.1 mL of the
251 suspension on each of 6-10 TSA plates. Incubate the plates for 18-24 h at
252 $36\pm 1^{\circ}\text{C}$.

253 E. Following the incubation of the agar plates from I.D, place approximately 5 mL
254 sterile cryoprotectant solution on the surface of each plate. Re-suspend the
255 growth in the cryoprotectant solution using a sterile spreader without damaging
256 the agar surface. Aspirate the suspension from the plate with a pipette and place it
257 in a sterile vessel large enough to hold about 30 mL. Repeat the growth harvesting
258 procedure with the remaining plates and continue adding the suspension to the
259 vessel (more than 1 tube may be used if necessary). Mix the contents of the
260 vessel(s) thoroughly; if more than 1 vessel is used, pool the vessels prior to
261 aliquoting culture.

262 F. Immediately after mixing, dispense 0.5-1 mL aliquots of the harvested suspension
263 into cryovials; these represent the frozen stock cultures.

264 G. Store the cryovials at -70°C or lower for a maximum 18 months then reinitiate
265 with a new lyophilized culture.

266 II. QC of Stock Cultures.

267 A. Conduct a purity check (isolation streak) throughout each step of the frozen stock
268 culture generation process prior to freezing the culture.

269 1. For section I.B, conduct a streak isolation onto BAP from the rehydrated

- 270 lyophilized culture. In addition, streak a loopful onto both MSA and
271 Cetrimide.
- 272 2. For section I.C, conduct a streak isolation onto BAP from the TSA or NA
273 18-24 hour plate. In addition, streak a loopful onto both MSA and
274 Cetrimide.
- 275 3. For section I.D, conduct a streak isolation onto BAP from the resuspended
276 culture tube. In addition, streak a loopful onto both MSA and Cetrimide.
- 277 B. Conduct QC of the pooled culture concurrently with freezing (section I.E). Streak
278 a loopful on a plate of BAP. In addition, streak a loopful onto both MSA and
279 Cetrimide. Incubate all plates at $36\pm 1^{\circ}\text{C}$ for 24 ± 2 hours.
- 280 C. Following the incubation period, record the colony morphology as observed on
281 the BAPs and selective media plates (including the absence of growth) and Gram
282 stain. See Attachment 2 for details on cell and colony morphology, colony
283 characteristics on selective media, and stain reactions.
- 284 D. For each organism, perform a Gram stain (refer to VII.E) from growth taken from
285 the BAPs according to the manufacturer's instructions. Observe the Gram
286 reaction by using brightfield microscopy at 1000X magnification (oil immersion).
- 287 E. Conduct additional biochemical and antigenic analyses; refer to sections VII.G-
288 VII.H for *S. aureus* and VII.F for *P. aeruginosa* for instructions. Alternatively,
289 Vitek 2 Compact may be used.
- 290 F. Record all confirmation results.

291 Attachment 2

292 Typical Growth Characteristics of strains of *P. aeruginosa* (see ref. VII.C) and *S. aureus* (see ref.
293 VII.D).

294

	<i>P. aeruginosa</i> *	<i>S. aureus</i> *
Gram stain reaction	(-)	(+)
Typical Growth Characteristics on Solid Media		
Mannitol Salt	No Growth	circular, small, yellow colonies, agar turning fluorescent yellow
Cetrimide	circular, small, initially opaque, turning fluorescent green over time; agar fluorescent yellowish green	No Growth
Blood agar (BAP)	flat, opaque to off-white, round spreading (1), metallic sheen, slightly beta hemolytic	small, circular, yellow or white, glistening, beta hemolytic
Typical Microscopic Characteristics		
Cell dimensions	0.5-1.0 μm in diameter by 1.5-5.0 μm in length*	0.5-1.5 μm in diameter*
Cell appearance	straight or slightly curved rods, single polar flagella, rods formed in chains	spherical, occurring singly, in pairs and tetrads, sometimes forming irregular clusters

295 *After 24 \pm 2 hours

296 (1) *P. aeruginosa* may display three colony types: a) circular, undulate edge, convex, rough and opaque; b) circular,
297 entire edge, convex, smooth and translucent; c) irregular, undulate edge, convex, rough, spreading, and translucent.
298 Pyocyanin is not produced.

299 Attachment 3
300
301 Reactor/Tubing Setup for Generation of *S. aureus* Biofilm
302



303
304