

Methods for Ambient Speciated Mercury Monitoring

National Air Monitoring Conference

November 8, 2006
Las Vegas, Nevada

Matthew S. Landis

U.S. EPA Office of Research and Development, RTP, NC

Acknowledgements

Frank Schaedlich & Dan Schneeberger
Tekran, Inc.


Robert Stevens & Tom Atkeson
Florida Department of Environmental Protection

Eric Prestbo
Frontier Geosciences, Inc.

RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions





RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

WARNING

The Florida Department of Health and Rehabilitative Services has issued a health advisory urging limited consumption of largemouth bass and warmouth caught in certain portions of the Everglades due to excessive accumulation of the element mercury.

- Fish caught in Arthur R. Marshall Loxahatchee National Wildlife Refuge (Water Conservation Area 1) should not be eaten more than once per week by adults and not more than once per month by children under 15 and pregnant women.
- Fish caught in Water Conservation Areas 2a and 3 should not be eaten at all.

For additional information, contact the Florida Department of Health and Rehabilitative Services at (405) 355-3018.

SENT IN YOUR WALLET

THURSDAY, August 10, 1989

Sentinel

RECEIVED
AUG 17 1989
ENVIRONMENTAL HAZARDS
EPIDEMIOLOGY

Mercury pollutes 3 top fishing lakes

State officials: Mercury's effect on health unclear

One mistake was telling pregnant women to limit consumption of mercury-tainted fish. In fact, all women of childbearing age are supposed to limit the mercury-laden fish they eat in case they become pregnant.

Another problem with the early advice was they contemplated fish portions of half pound but did not say that. As many people felt safe in eating several pounds of fish a couple of times a week. That's what the problem is in a B...

High mercury levels taint S. Florida fish

BROWARD METRO
Sunday, September 15, 1991 \$1



Mercury in North Florida fish too high

By Carrie Teegardin
Democrat staff writer

High-mercury areas

MIAMI — The...
nio Noriega is tak...
tion into uncharted...
provide the federal...
broadest test yet of...
ment applies to gov...
against foreign citizen...
The criminal case is a...
opening statement b...
nt on Monday morning...
he historic drug-traf...
eteering trial has bee...
victors as just another...
constitutional scholars...
say the Noriega cas...
it routine...
say that if Noriega is co...
conviction is appeal...
eme Court and uphe...
tely to establish pre...
give U.S. milit...
l officials...

Ambient Measurements to Support Atmospheric Mercury Research

- Elucidate Atmospheric Chemistry
- Quantify Impacts from Specific Sources
- Evaluate Deterministic Models
- Provide Emission Reduction Accountability

Why Speciate Mercury?


➤ Species have Different Behaviors

- Elemental Mercury: Hg^0
- Reactive Mercury: Hg^{2+}
- Particulate Mercury: $\text{Hg}(\text{p})$

➤ Atmospheric Transport & Deposition Modeling

➤ Bioaccumulation, Exposure & Risk Assessment

Contemporary Atmospheric Mercury Conceptual Model




Global - Regional - Local

RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Spatial Distribution of Global Atmospheric Mercury Emissions




Source: UNEP Global Mercury Assessment, 2002, using J. Pacyna 1995 data, as presented by AMAP (1998)

RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Anthropogenic Air Emissions of Mercury: Distribution by Industrial Sector in 1995



Total: 2,382 metric tons

- Coal and fuel combustion is by far the largest source category
- Estimates are rough; most countries do not have Hg inventories
- We need to further develop reliable emissions inventories

Source: EPA estimates derived from UNEP Global Mercury Assessment, UNEP, Geneva, December 2002

RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Anthropogenic Hg Emissions Observed as Hg²⁺

Source	N	Hg ²⁺ /Total Hg(%) \pm Std. Dev.
Cement Production ^a	3	25 \pm 4
Incinerator (medical) ^a	3	95 \pm 5
Incinerator (municipal) ^a	8	78 \pm 8
Coal-fired Utility Boiler ^b	19	67 \pm 27

^a Stevens *et al.*, 1996

^b Prestbo and Bloom, 1995

RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions


Ambient RGM Methodologies

- **Impregnated Filter**
- **Refluxing Mist Chamber**
- **Tubular Denuder**

Challenges Measuring RGM

- Reject much larger *elemental* component
- Method must be 1-2 orders of magnitude more sensitive than total mercury methods
- Quantitatively pass *RGM* to the collector
- Exclude *particulate phase* mercury while avoiding filter artifact

Quartz Annular Denuder




* Landis et al. *Environ. Sci. Technol.*, **2002**, 36, 3000-3009

RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions


SEM Micrographs of Quartz Surfaces



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions


Laboratory Evaluation



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions



Collection Efficiency for HgCl_2



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Annular Denuder Performance




RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Importance of Properly Heating Sampling System


Barrow, Alaska (2001)



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Collocated Denuded & Undenuded Hg(p)



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

KCl Denuder Method Detection Limits

Field Blanks

N	Mean (pg)	Std Dev (pg)
66	2.22	1.24

10 LPM Flow Rate

Duration (hours)	Air Volume (m ³)	MDL (pg m ⁻³)
1	0.6	4.1
2	1.2	2.1
6	3.6	0.7
12	7.2	0.3

Observed RGM Concentrations

Site	N	Mean (pg m ⁻³)	Min (pg m ⁻³)	Max (pg m ⁻³)	Std Dev (pg m ⁻³)
Baltimore, MD	30	22.7	5.4	138.7	25.9
Everglades, FL	45	15.4	2.5	54.3	12.0
RTP, NC	26	15.5	3.5	51.0	12.0

Field Monitoring Setup



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Manual Denuder Sampling Box



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions


Analysis of Manual Denuder



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Flow Diagram of Tekran 2537A Mercury Analyzer




RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Tekran System Description

- **Model 2537A Vapor Phase Mercury Analyzer**
 - Semi-continuous Hg⁰ Measurements
 - Gold Traps for Hg⁰ Pre-concentration
 - Cold Vapor Atomic Fluorescence Spectrometer
- **Model 1130 Speciation Unit**
 - KCl Thermal Annular Denuder (Hg²⁺)
 - Zero Air Source & Pumping Module
- **Model 1135 Particulate Unit**
 - Thermal Quartz Filter & Pyrolyzer Column

Automated Speciation Instrumentation



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Typical External Installation



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Typical Internal Installation




RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Mauna Loa Hg Time Series

7/31 – 9/20, 2001



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Tekran System SOP

➤ Weekly

- Replace: denuder, soda & lime trap, impactor frit, filters
- Clean: inlet glassware, impactor, couplers
- Check: trap performance, flows, heaters
- Conduct perm tube calibration

➤ Monthly

- Manual injections (2537A, 1130 inlet)

➤ Bi-monthly

- Replace RPF, zero air filters

➤ Six Months

- Calibrate 2537A permeation tube

➤ Annual

- Replace: 1130 Carbon traps, pump brushes, 2537A lamp
- Calibrate 2505 syringe, flow controllers

QA/QC Inlet Injection



RESEARCH & DEVELOPMENT

Building a scientific foundation for sound environmental decisions

Speciation Methods Summary

- KCl-Coated annular denuders quantify RGM without known interference problems when following SOP
- RGM is thermally ($\sim 500^{\circ}\text{C}$) converted to Hg^0
- $\text{Hg}(\text{p})$ is thermally ($\sim 800^{\circ}\text{C}$) converted to Hg^0
- Low denuder MDL's allow high-resolution sampling
- Manual method is inexpensive, simple, and mobile
- Tekran provides automated high-resolution data

Questions ???

