

TECHNICAL SUPPORT SECTION EFFICACY REVIEW - II

Disinfectants Branch

EPA Reg. No. or File Symbol 46781-1

Date Division Received 07-01-87

Data Accession No.(s). 402521-01 & 402521-02

Product Manager No. PM-31 (Lee)

Product Name Metricide MX-1400 Activated Dialdehyde Solution

Company Name Metrex Research Corporation

202.0 Recommendations

202.1 Efficacy Supported by the Data

- a. The submitted sporicidal data support efficacy of of this product as a sterilant after reuse of the activated undiluted (2% Glutaraldehyde) solution for 14 days in (bucket or tray) manual systems at a contact time of 8 hours at 20°C when tested by the A.O.A.C. Sporicidal Test with 2 different batches.

203.0 Labeling

No revised labels were submitted.