

Technical Support Section Efficacy Review-II

Disinfectants Branch

EPA.Reg. No.or File Symbol 46781-2

Data Division Received 10-17-85

Data Accession No(s).

Product Manager No. PM 31 (Lee)

Product Name Metricide MX-2800

Company Name Metrex Research Corporation

202.0 Recommendations

202.1 Additional Data/Information

The submitted re-use sporicidal data report for Sample A fulfills the batch sample requirement to support sterilization claims for the re-use solution for a 10-day period.

The corrected copy of the Schedule of Bioburden Tasks for Sample B provides evidence that the testing was conducted in the same manner as Sample A and provided a cumulative bioburden load equivalent for a 10-day re-use period.

In view of the additional data/information submitted, an acceptable re-use claim would read: "This solution may be re-used for up to 10 days after activation for both disinfection and sterilization. - Do not use the activated stock solution beyond 28 days."