

SEP 24 1985

Metrex Research Corporation  
11270 S. Dransfeldt Road (Shipping Only)  
P.O. Box 965  
Parker, CO 80134

Attn: David Stump

Subject: Metricide Activated Dialdehyde Solution  
EPA Reg. No. 46581-1  
Submission dated July 31, 1985

The amendment referred to above, submitted in connection with registration under FIFRA to support subject product for disinfecting and sterilizing for a 14 day re-use period in a manual system is not acceptable. However, the submitted and reviewed data will support a re-use claim for seven (7) days.

Submit five (5) copies of labeling to reflect the supporting data, revise the re-use claim to read: "This solution may be re-used for up to seven (7) days after activation for both disinfection and sterilization. Do not use the activated stock solution beyond 14 days".

Refer to the attached Technical Support Section Efficacy Review II for a detailed analysis of the submitted data/information.

Sincerely,

John H. Lee  
Product Manager 31  
Disinfectants Branch  
Registration Division (TS-767C)

Attachment