

# RADIOL GERMICIDAL SOLUTION

504-6  
13

A solution to be used full strength for the chemical  
disinfection of surgical instruments and appliances

Item No.  
91500

## One Quart



**Active Ingredients 0.14 %**  
Cetyl dimethyl ethyl ammonium  
bromide 0.07%  
n-Alkyl (50% C<sub>12</sub>, 40% C<sub>14</sub>, 10% C<sub>16</sub>)  
dimethyl  
benzyl ammonium chloride 0.07%

**Inert Ingredients 99.86 %**  
Water, Sodium Nitrite, Dye

**CAUTION: Keep out  
of the reach of children**

**Harmful if swallowed**

Rinse empty container thor-  
oughly with water and dis-  
card it.

**Caution:** Avoid contact with  
skin and eyes. In case of con-  
tact, wash thoroughly with  
water. If irritation persists get  
medical attention. Avoid con-  
tamination of food.

**Note:** In the presence of spore  
organisms, mycobacterium tu-  
berculosis, or etiologic agent  
for viral hepatitis, quaternary  
ammonium compounds cannot  
be relied upon to produce  
thorough disinfection.

See accompanying circular for  
directions and cautions.

EPA Reg. No. 504 6

**NOVOCOL CHEMICAL MFG. CO., INC.**  
Brooklyn, New York 11207

Toronto 570, Ontario