

Cameraman's Note

Poor Copy

ALWAYS DILUTE ACCORDING TO  
"DILUTION TABLE" BEFORE USING

- DIRECTIONS FOR USE -

• Do Not Add to Soap or Detergent Solution

1. FOR ROUTINE USE: Swab, wipe or rinse surfaces (such as floors, walls, urinals, washstands, garbage pails) thoroughly, using one per cent solution (see "DILUTION TABLE"). Disinfect smaller objects by immersion in one per cent solution for at least 10 minutes. Permit all surfaces to dry before use.

2. FOR SPECIAL USE:

A. If known heavy contamination exists, use two per cent solution on surfaces. Swab, wipe or rinse surfaces with disinfectant solution at least twice, allowing solution to act for at least 10 minutes between applications. Disinfect smaller objects by immersion for at least 10 minutes. Permit all surfaces to dry before use.

(continued on far panel)

55-1-60-500

Pkd.1/60

Lot 634-2

ACCEPTED

AUG 30 1967

UNDER THE FEDERAL INSECTICIDE  
FUNGICIDE AND RODENTICIDE ACT  
FOR ECONOMIC POISON REGISTERED  
UNDER NO. 1691-81 SUBJECT  
TO ATTACHED COMMENTS.

1691-81

1 U.S. GALLON

Stk. No. GF6840-664-0541

DISINFECTANT

GERMICIDAL and FUNGICIDAL  
CONCENTRATE

FS O-D-406, Amend 2

USDA Reg. No. 1691-81

ACTIVE INGREDIENTS:

Soap  
Isopropyl Alcohol  
Mixed 4 and 6-chloro-2-phenylphenol.

22.54%

INERT INGREDIENTS:

Water

TOTAL

77.46%  
100.00%

CAUTION Keep Out of Reach of Children  
CHEMICAL COMPOUNDING CORP.

Jersey City, N.J.

3. P  
S  
L  
I  
S  
H  
  
Di  
  
One Per  
(Row  
Two Per  
(Row  
  
Rinse

**ACCEPTED**

**AUG 30 1967**

UNDER THE FEDERAL INSECTICIDE  
FUNGICIDE AND RODENTICIDE ACT  
FOR ECONOMIC POISON REGISTERED  
UNDER NO. 1691-81 SUBJECT  
TO ATTACHED COMMENTS.

1691-81

Cameraman's Note

Poor Copy

1 U.S. GALLON

Stk. No. GF6840-664-0541

**DISINFECTANT**

GERMICIDAL and FUNGICIDAL  
CONCENTRATE

FS O-D-406, Amend 2

**USDA Reg. No. 1691-81**

**ACTIVE INGREDIENTS:**

Soap	22.54%
Isopropyl Alcohol	
Mixed 4 and 6-chloro-2-phenylphenol.	

**INERT INGREDIENTS:**

Water	77.46%
<b>TOTAL</b>	<b>100.00%</b>

**CAUTION** Keep Out of Reach of Children  
CHEMICAL COMPOUNDING CORP.  
Jersey City, N. J.

(continued from first panel)

B. Cold disinfection of surgical instruments - sterilization by autoclave shall be performed wherever possible. Cold disinfection of surgical instruments with this material may be used as an emergency alternate. In the absence of heat, immerse in 2% solution for at least 10 minutes after removing excess soil by cleaning.

3. Precaution - To disinfect RUBBER and PLASTIC objects, swab surfaces if possible. Immerse only if required. LIMIT EXPOSURE TO 10 MINUTES. PROLONGED IMMERSION IN DISINFECTANT MAY AFFECT PLASTICS AND RUBBER SERIOUSLY. Do not get in eyes. Wash thoroughly after using. Harmful if swallowed.

DILUTION TABLE

Dilution	Ounce concen- trate per gallon water	Ml concen- trate per liter water	Teaspoons concentrate per gallon water
One Per Cent Solution (Routine Use)	1½	10	10
Two Per Cent Solution (Routine Use)	2½	20	20

Cont. No. NSCO-1-54336

Rinse empty container thoroughly with water, and  
discard it.