

Active Ingredients—N-Alkyl (C14, 50%; C12, 40%; C16, 10%) dimethyl Benzyl Ammonium Chloride	10%
inert Ingredients	90%
Phenol Coefficients (A.O.A.C. method 20 C.)	
<i>Salmonella typhosa</i>	35
<i>Micrococcus pyogenes</i> var. <i>aureus</i>	50
U.S. Registration No. 3837-5	

Avoid Contact with Children
Avoid Consumption of Food
Harmful If Swallowed. Avoid Contact with
Skin and Eyes
In Case of Contact, Flush with Plenty
of Water



Cold disinfection of surgical instruments:
Clean thoroughly; rinse, and immerse in a solution of 1 ounce of disinfectant per gallon of water for 10 minutes.

**A Product of
PHOENIX EMBALMING
CHEMICAL CO.
Concordia, Kansas**