

PM 31

6658-39

10/2

USE THE "MAGIC TOUCH"

Pro-Magic

formulated for professionals!

Disinfectant Deodorant

Disinfects Cleans Deodorizes

ACTIVE INGREDIENTS

n-Alkyl (60% C ₁₂ , 30% C ₁₄ , 5% C ₁₆ , 5% C ₁₈) dimethyl benzyl ammonium chlorides	0.1%
n-Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzyl ammonium chloride	0.1%
Isopropanol	53.0%
INERT INGREDIENTS	46.8%
Total	100.0%

EPA REG NO 6658-39 EPA EST 6658 MO-1

WARNING:

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMAN AND DOMESTIC ANIMALS

Causes eye and skin irritation. Do not get in eyes, skin or on clothing. Harmful if swallowed.

STATEMENT OF PRACTICAL TREATMENT

In case of contact, immediately flush eyes, or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink egg whites gelatin solution, or if these are not available, drink large quantities of water. Call a physician.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage.

STORAGE AND DISPOSAL

STORAGE: Do not puncture or incinerate container. Do not expose to heat or store at temperatures above 120°F.

DISPOSAL: Do not reuse empty container. Wrap container and put in trash collection.

DIRECTIONS FOR USE

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

SURFACE DISINFECTION AND DEODORIZING

Disinfectant Deodorant kills Staphylococcus aureus, Salmonella choleraesuis and Pseudomonas aeruginosa on inanimate environmental surfaces.

For use on telephone, door knobs, light switches, bed springs, diaper pails, toilet seats, sink tops, bed pans, garbage pails and desk tops.

For Hospitals, Institutions, Nursing Homes, Offices, Schools, Motels, and Hotels.

DISINFECTING: For treatment of previously cleaned hard non-porous surfaces, spray so as to completely cover and allow to remain for a contact time of ten minutes. Wipe off with clean cloth or sponge.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

ACCEPTED
MAY 17 1996

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 6658-39

1-17-96

KEEP OUT OF REACH OF CHILDREN

WARNING:

PHYSICAL OR CHEMICAL HAZARDS

Flammable. Contents Under Pressure.

Do not use near fire, sparks or flame. Do not puncture or incinerate container. Exposure to temperatures above 120°F may cause bursting. See back panel for additional precautionary statements

MANUFACTURED BY:

Midco Products Co., Inc.
Chesterfield, MO 63005-1005

NET WT. 15 OZ.

019618GA

#18