

Cidex

Activator Vial

Arbrook Inc.

Directions: See Back Label for Use
FORMULA

DANGER: Contact with eyes causes
irritation.

Sterilizing and Disinfecting Solution

Sporicidal
Virucidal
Fungicidal

Tuberculocidal
Bactericidal
Pseudomonacidal

aqueous formula
Cidex7*

long-life activated dialdehyde

Product
CF-750

Active ingredient
Glutaraldehyde ... 2.0%
Inert Ingredients ... 98.0%
Total ... 100.0%

Contains Sodium Nitrite
as a corrosion inhibitor.

Note — Attached vial
contents *must* be added
to solution before this
product is effective.
See directions.


Rinse empty container
thoroughly with water
prior to discarding.

DANGER: See Back Label for Precautions

* TRADEMARK

† See Directions for Use

Arbrook*



Sterilizing and Disinfecting Solution

Sporicidal
Virucidal †
Fungicidal

Tuberculocidal
Bactericidal
Pseudomonacidal

Cidex aqueous formula **7** *

long-life activated dialdehyde

Product

CF-750

Active Ingredient	
Glutaraldehyde	2.0%
Inert Ingredients	98.0%
Total	100.0%

Contains Sodium Nitrite
as a corrosion inhibitor

Note — Attached vial
contents *must* be added
to solution before this
product is effective
See directions

Rinse empty container
thoroughly with water
prior to discarding.

DANGER: See Back Label for Precautions

Arbrook.

Directions for Use:

DANGER:

Directions for Activating:

KEEP IN A COOL PLACE

When activated, solution changes color to green.

Do not use this solution beyond 28 days after mixing.

NOTE DATES HERE

MIXING DATE _____ EXPIRATION DATE _____

E. P. A. Reg. No. 7078-4

E. P. A. Est. No. 7078-TX-1

CONTENT ONE GAL. (128 FLUID OZ.)

MADE IN U.S.A.

PAT. 3016328

CONTROL NUMBER ON BOTTOM

© Arbrook, Inc. 1975

Arbrook inc.
Arlington, Texas 76010