

No. 5AL97250

50 PACKETS
2/3 FL. OZ. EACH



**BENZALKONIUM
CHLORIDE**
Tincture Solution
1:750

CAT. #L97250

50 PACKETS
2/3 FL. OZ. EACH



**BENZALKONIUM
CHLORIDE**
Tincture Solution
1:750

CAT. #L97250

**CAUTION. Keep
out of reach of
children. Harmful
if swallowed.**

ACTIVE INGREDIENTS: alkyl*
dimethyl benzyl ammonium
chloride0.13%
Isopropyl alcohol37%
INERT INGREDIENTS:
Water62.87%
*C₁₂-67%, C₁₄-25%, C₁₆-7%, C₈,C₁₀
and C₁₈-1%.

E.P.A. REGISTRATION NO. 9480-2
E.P.A. Est.NO. 9480-NY-1

E.P.A. Est. NO. 9480-IN-1

**CAUTION. Keep
out of reach of
children. Harmful
if swallowed.**

ACTIVE INGREDIENTS: alkyl*
dimethyl benzyl ammonium
chloride0.13%
Isopropyl alcohol37%
INERT INGREDIENTS:
Water62.87%
*C₁₂-67%, C₁₄-25%, C₁₆-7%, C₈,C₁₀
and C₁₈-1%.

E.P.A. REGISTRATIONNO. 9480-2
E.P.A. Est.NO. 9480-NY-1

E.P.A. Est.NO. 9480-IN-1

For use in thermometer holders. Immerse precleaned Thermometer in undiluted solution in holder for 10 minutes to disinfect, then store thermometer in holder at bedside until needed. Change solution in holder every 7 days or with each new patient.

Do not use, pour, spill or store near heat or open flame.

Supplied by
 PROFESSIONAL DISPOSABLES INC.
MT. VERNON, N.Y. 10550


SIDE PANEL

For use in thermometer holders. Immerse precleaned Thermometer in undiluted solution in holder for 10 minutes to disinfect, then store thermometer in holder at bedside until needed. Change solution in holder every 7 days or with each new patient.

Do not use, pour, spill or store near heat or open flame.

Supplied by
 PROFESSIONAL DISPOSABLES INC.
MT. VERNON, N.Y. 10550

FRONT PANEL



PACKING BOX