

ALLERGAN SOAKARE® SOAKING SOLUTION - 5 cc. Label J.B. #2067

(front)

1/6 fl. oz. Sterile

ALLERGAN
SOAKARE®
CONTACT
LENS
SOAKING
SOLUTION

(logo)

(back)

Allergan Soakare Soaking Solution is a germicidal solution specially designed to provide maximal conditioning and protection for contact lenses when not in use.

Active Ingredients: disodium edetate 0.25%, benzalkonium chloride 0.01%.

Inert Ingredients: 99.74%.

Directions: The lens storage container and/or carrying case should be kept amply supplied with Allergan Soakare Soaking Solution at all times. The lenses should be completely immersed for at least 30 minutes. Always keep lenses in the storage container between wearings. Just prior to use, rinse lenses and treat with Allergan LIQUIFILM® Wetting Solution.

Allergan Pharmaceuticals
Irvine, California 92664 U.S.A.

ACCEPTED WITH COMMENTS

USDA Reg. No. 9992-1
Control No.

SOAKARE[®] SOAKING SOLUTION - 120 cc. Unit Carton J.B. #1792

(front)

Sterile

ALLERGAN
SOAKARE[®]
CONTACT
LENS
SOAKING
SOLUTION

4 fl. oz.

(logo)

(right)

Directions: The lens storage container and/or carrying case should be kept amply supplied with Allergan Soakare Solution at all times. The lenses should be completely immersed for at least 30 minutes. Always keep lenses in the storage container between wearings. Just prior to use, rinse lenses and treat with Allergan LIQUIFILM[®] Wetting Solution.

Allergan Pharmaceuticals
Irvine, California 92664 U.S.A.

(left)

Allergan Soakare Soaking Solution is a germicidal solution specially designed to provide maximal conditioning and protection for contact lenses when not in use.

Active Ingredients: disodium edetate 0.25%, benzalkonium chloride 0.01%.

Inert Ingredients: 99.74%.

USDA Reg. No. 9992-1

(back)

Sterile

ALLERGAN
SOAKARE®
CONTACT
LENS
SOAKING
SOLUTION

4 fl. oz.

(logo)

(top)

No. 016

ALLERGAN
SOAKARE®
CONTACT LENS
SOAKING
SOLUTION

4 fl. oz.

(bottom)

Control No.

(front)

½ fl. oz. Sterile

ALLERGAN
SOAKARE[®]
CONTACT
LENS
SOAKING
SOLUTION

(logo)

(back)

Allergan Soakare Soaking Solution is a germicidal solution specially designed to provide maximal conditioning and protection for contact lenses when not in use.

Active Ingredients: disodium edetate 0.25%, benzalkonium chloride 0.01%.

Inert Ingredients: 99.74%.

Directions: The lens storage container and/or carrying case should be kept amply supplied with Allergan Soakare Soaking Solution at all times. The lenses should be completely immersed for at least 30 minutes. Always keep lenses in the storage container between wearings. Just prior to use, rinse lenses and treat with Allergan LIQUIFILM[®] Wetting Solution.

Allergan Pharmaceuticals
Irvine, California 92664 U.S.A.

USDA Reg. No. 9992-1
Control No.

ALLERGAN SOAKARE[®] SOAKING SOLUTION - 120 cc. Label

(front)

Sterile

ALLERGAN
SOAKARE[®]
CONTACT
LENS
SOAKING
SOLUTION

4 fl. oz.

(logo)

(back)

Allergan SOAKARE[®] Soaking Solution is a germicidal solution specially designed to provide maximal conditioning and protection for contact lenses when not in use.

Active Ingredients: disodium edetate 0.25%, benzalkonium chloride 0.01%.

Inert Ingredients: 99.74%.

Directions: The lens storage container and/or carrying case should be kept amply supplied with Allergan Soakare Soaking Solution at all times. The lenses should be completely immersed for at least 30 minutes. Always keep lenses in the storage container between wearings. Just prior to use, rinse lenses and treat with Allergan LIQUIFILM[®] Wetting Solution.

Allergan Pharmaceuticals
Irvine, California 92664 U.S.A.

USDA Reg. No. 9992-1
Control No.