



TEMPCARE
Thermometer
germicidal
solution

CAUTION: KEEP OUT OF REACH
OF CHILDREN. (SEE DIRECTIONS FOR
USE ON REVERSE PANEL.)

U.S.D.A. REG. No. 10197-1
4 FL. OZ.

Front

WHITE
COPY

DIRECTIONS:

Fill Tempcare holder to mid-point—immerse thermometer. Wipe dry before using, after use, clean thermometer thoroughly in water, rinse and immerse in solution. Change solution daily. Remove solution when storing holder between illnesses.

CONTAINS:

Active ingredients:	
n-alkyl (50% C14, 40% C12, 10% C16)	
dimethyl benzyl ammonium chloride	0.135%
Inert ingredients:	<u>99.865%</u>
	100.000%

dist. by
EISELE & COMPANY, INC.
NASHVILLE, TENNESSEE 37208
A subsidiary of SHULTON, INC.

MADE AND PTD. IN U.S.A. E8-1

Back