

EXHIBIT 2
COPY FOR CARION

AMERICAN HOSPITAL SUPPLY CORPORATION
DIVISION OF AMERICAN HOSPITAL SUPPLY CORPORATION
ILLINOIS 60073

TOMAC SENTRY™ THERMOMETER KIT

ACCEPTED
5 MAY 1980
RECEIVED
MAY 1980
MAY 1980

DIRECTIONS FOR USE:

- Remove contents from box
- Pull off plastic cover
- Remove thermometer
- Snap off center post on cap
- Insert thermometer
- Replace plastic cover
- Clean and dry thermometer after each use

PATIENT NAME _____

ROOM NO _____

CATALOG NUMBER
TOMAC SENTRY™ THERMOMETER KIT ORAL

- CONTENTS:**
- Thermometer
 - Prefilled Thermometer Vial™
 - Holder and Cover

• 70% by weight Isopropanol
Storage and Disposal:
Rinse vial with water before releasing to patient.

CATALOG NUMBER
TOMAC SENTRY™ THERMOMETER KIT ORAL

FOR INDIVIDUAL PATIENT USE ONLY
Rinse vial with water before releasing to patient.
Thermometer solution contains:
70% by weight Isopropanol solution for disinfecting clinical thermometer.
ACTIVE INGREDIENT:
Isopropanol 70%
INERT INGREDIENT 30%
EPA Reg No 10287-7 EPA Est 33695-11-1
WARNING - KEEP OUT OF REACH OF CHILDREN
Net Contents 0.67 Fl. Oz.

DISTRIBUTED BY
AMERICAN HOSPITAL SUPPLY
DIVISION OF AMERICAN HOSPITAL SUPPLY
CORPORATION • MCGRAW PARK • ILLINOIS 60073

TOMAC SENTRY™ THERMOMETER KIT

This thermometer conforms to all the requirements established in Voluntary Product Standard PS 39-70 Published by the National Bureau of Standards

BEST DOCUMENT AVAILABLE

