

# ETHYLENE OXIDE AND CARBON DIOXIDE STERILIZING GAS

**FOR INDUSTRIAL AND MEDICAL USE**

**ONLY IN COMMERCIAL GAS STERILIZERS**

## ETHYLENE OXIDE AND CARBON DIOXIDE MIXTURE

ACTIVE INGREDIENT:  
ETHYLENE OXIDE 10% BY WEIGHT  
INERT INGREDIENT:  
CARBON DIOXIDE 90% BY WEIGHT  
TOTAL CONTENTS \_\_\_\_\_ LBS. NET  
E.P.A. REG. NO. 11491-1  
EXP. 2/78


### WARNING: CONTENTS UNDER PRESSURE!

1. Do not puncture. 2. Do not use or store near heat or open flame. Exposures above 130°F may cause bursting. Never throw containers into fire or incinerator. 3. Do not use oil or grease on valve, regulator, fittings. 4. Do not recharge without written consent of Liquid Carbonic.

### WARNING: DO NOT IGNORE!

#### KEEP OUT OF REACH OF CHILDREN

1. Avoid contact with eyes, skin, or clothing. Avoid breathing vapor. In case of contact, remove all contaminated clothes immediately and flush skin or eyes with plenty of water for at least 15 minutes; for eyes, get medical attention. 2. Use in accordance with sterilizer manufacturer's operating instructions. 3. Ethylene oxide is subject to polymerization; use within 2 months after receipt. 4. Instructions for use and other information are in Industrial Gas Bulletin IG-2013.



**LIQUID CARBONIC**

Subsidiary of Houston Natural Gas Corporation  
135 South LaSalle Street • Chicago, Illinois 60603

Form #3121 Rev. 4/74

**A RED DIAMOND<sup>®</sup> PRODUCT (Non Life Supporting)**

**ACCEPTED**