

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

- 2 -

- d. Revise the statent: "Wash clothing before reuse" to read:
"Remove contaminated clothing and wash before reuse."
- e. Add the following additional precautionary statement:
Wash thoroughly with soap and water after handling.
- f. Revise the Environmental Hazards statement to comply with
PR Notive 93-10.

3. Make the labeling changes listed below before you release
the product for shipment: (Tag Label)

- a. Include the EPA Registration No.
- b. Add a referral statement to read: See outer label on
cylinder for additional precationary statements and Directions
for Use.
- c. Include the full company name on the tag.

If these conditions are not complied with, the registration
will be subject to cancellation in accordance with FIFRA sec.
6(e). Your release for shipment of the product constitutes
acceptance of these conditions.

Submit one (1) copy of the final printed label prior to
release of the product for shipment.

A stamped copy of the label is enclosed for your records.

Sincerely,

Marion Johnson
Product Manager (31)
Antimicrobial Program Branch
Registration Division (7505C)

Enclosures

CONCURRENCES							
SYMBOL							
SURNAME							
DATE							

307
ACTIVE INGREDIENT:
ETHYLENE OXIDE..... 8.5%

INERT INGREDIENTS
CARBON DIOXIDE..... 20.8%
HELIUM..... 0.7%

TOTAL..... 100.0%

KEEP OUT OF REACH OF CHILDREN

DANGER! PELIGRO!

PRECAUTION AL USARIO: SI USTED NO LEE INGLES, NO USE ESTE PRODUCTO
HASTA QUE LA ETIQUETA LE HAYA SIDO EXPLICADO AMPLIAMENTE

STATEMENT OF PRACTICAL TREATMENT

IF INHALED OR SUSPICION OF INHALATION: REMOVE PERSON TO FRESH AIR, KEEP WARM
IF BREATHING HAS STOPPED, GIVE ARTIFICIAL RESPIRATION. IF BREATHING IS
DIFFICULT, GIVE OXYGEN. CALL A PHYSICIAN EVEN IF NO SYMPTOMS ARE PRESENT
KEEP UNDER MEDICAL OBSERVATION. SYMPTOMS MAY BE DELAYED

IN CASE OF CONTACT WITH EYES SKIN OR CLOTHING IMMEDIATELY FLUSH EYES OR
SKIN WITH PLENTY OF WATER FOR AT LEAST 15 MINUTES WHILE REMOVING
CONTAMINATED CLOTHING AND SHOES. CALL A PHYSICIAN. WASH CLOTHING BEFORE
REUSE. DISCARD CONTAMINATED SHOES. NOTE TO PHYSICIAN: SEE ATTACHED
CYLINDER TAG.

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS:

DANGER. NON-FLAMMABLE LIQUID AND GAS UNDER PRESSURE. CAN CAUSE EYE AND
SKIN BURNS DO NOT GET LIQUID IN EYES, ON SKIN, OR ON CLOTHING.
MAY BE HARMFUL OR FATAL IF INHALED. DO NOT BREATHE GAS. CAN CAUSE IRRITATION
OF THE RESPIRATORY TRACT. SYMPTOMS OF EXPOSURE MAY BE DELAYED

DANGER. CONTAINS ETHYLENE OXIDE A CANCER AND REPRODUCTIVE
HAZARD.

ENVIRONMENTAL HAZARDS:

DO NOT DISCHARGE EFFLUENT CONTAINING THIS PRODUCT INTO LAKES, STREAMS,
PONDS, ESTUARIES, OCEANS OR PUBLIC WATERS UNLESS THIS IS SPECIFICALLY
IDENTIFIED AND ADDRESSED IN AN N.P.D.E.S. PERMIT. DO NOT DISCHARGE EFFLUENT
CONTAINING THIS PRODUCT TO SEWER SYSTEMS WITHOUT PREVIOUSLY NOTIFYING THE
SEWAGE TREATMENT PLANT AUTHORITY. FOR GUIDANCE, CONTACT YOUR STATE WATER
BOARD OR REGIONAL OFFICE OF THE EPA.

PHYSICAL AND CHEMICAL HAZARDS:

CONTENTS UNDER PRESSURE. BACK FLOW INTO CONTAINER COULD CAUSE CYLINDER
RUPTURE ALWAYS USE A BACK FLOW PREVENTING DEVICE IN THE CONTAINER DISCHARGE
PIPING. CONTAINER TEMPERATURE SHOULD NOT EXCEED 125 °F (52°C). SECURE THE
CYLINDER IN AN UPRIGHT POSITION DURING STORAGE AND USE. KEEP PROTECTIVE
VALVE CAP IN PLACE WHEN NOT IN USE. USE ONLY IN CLOSED SYSTEMS WHICH ARE
RATED FOR THE CYLINDER PRESSURE (2000 PSIG). CLOSE VALVE AFTER EACH USE AND
WHEN EMPTY. CYLINDER MUST NOT BE REFILLED EXCEPT BY OR WITH THE CONSENT OF
LIQUID CARBONIC.

LIQUEFIED GASES N.O.S.,
UN 1956
(Ethylene Oxide, Carbon Dioxide)

NON
FLAMM
GAS
2

DIRECTIONS FOR USE: IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRO
WITH ITS LABELING. USE IN ACCORDANCE WITH STERILIZER MANUFACTURERS OF
PRODUCT SHOULD ONLY BE USED BY PERSONS EXPERIENCED IN THE USE OF ETH
MIXTURES OR UNDER THE SUPERVISION OF SUCH AN EXPERIENCED PERSON. SEE
FURTHER DIRECTIONS FOR USE. USE IN ACCORDANCE WITH THE MATERIAL SAFETY

STORAGE AND DISPOSAL: DO NOT CONTAMINATE WATER, FOOD OR FEED BY STOR
STORAGE: DO NOT STORE NEAR SOURCES OF HEAT, OPEN FLAMES, OR IN DIREC
TEMPERATURES MUST NOT EXCEED 125°F (52°C). TO MINIMIZE POLYMERIZATION
TEMPERATURE IS RECOMMENDED. DO NOT PUNCTURE OR DROP CYLINDERS.
PESTICIDE DISPOSAL: THE CONTENTS OF THIS CYLINDER MUST EITHER BE USE
RETURNED IN THE ORIGINAL CYLINDER TO LIQUID CARBONIC FOR HANDLING. DO
PRODUCT WITHOUT THE EXPRESS DIRECTION OF LIQUID CARBONIC

CONTAINER DISPOSAL: RETURN EMPTY AND OUT-DATED CYLINDERS FOR RE-USE.
VALVES FOR LEAKS IF LEAKS ARE DETECTED, CONTACT SUPPLIER

EPA REG. NO:

EPA ESTABLISHMENT NUMBER:

11491-CA-1 11491-FL-1 11491-
11491-CA-2 11491-TX-1

NET CONTENTS _____ LBS
EXPIRATION DATE _____

FOR INDUSTRIAL AND MEDICAL USE IN COMMERCIAL GAS
DESIGNED FOR USE WITH ETHYLENE OXIDE
**EMERGENCY RESPONSE PHONE # (504) 677-
CHEMTREC (800) 424-9300**



LIQUID CARBONIC
INDUSTRIES CORPORATION
810 JORIE BLVD
OAK BROOK, IL 60521-2216

MEDIFUME PLUS

LIQUID

LIQUID CARBONIC

810 JORIE BLVD
OAK BROOK, IL
60521-2216

EMERGENCY RESPONSE PHONE # (504) 673-8831
CHEMTREC (800) 424-9300

DO NOT REMOVE THIS TAG

DANGER! NON FLAMMABLE LIQUID AND GAS UNDER PRESSURE.
HARMFUL IF INHALED.

THIS TAG CONTAINS ADDITIONAL INFORMATION WHICH IS CONSIDERED AN
INTEGRAL PART OF THE LABEL FOR THIS PRODUCT.

NOTE TO PHYSICIAN: THIS PRODUCT CONTAINS ETHYLENE OXIDE AS
THE ACTIVE INGREDIENT. EXPOSURE OF EYES OR SKIN MAY RESULT IN
IMMEDIATE OR DELAYED SKIN IRRITATION AND BLISTERING WITH
POSSIBLE ALLERGIC SKIN REACTION. THE AFFECTED AREAS SHOULD BE
THOROUGHLY CLEANED AND TREATED SYMPTOMATICALLY. INHALATION MAY
PRODUCE A NUMBER OF EFFECTS, THE MOST IMMEDIATELY SERIOUS OF
WHICH IS PULMONARY EDEMA, A RESULT OF ETHYLENE OXIDE'S
RESPIRATORY TRACT IRRITANCY AND WHICH MAY BE DELAYED IN
ONSET. IF THIS OCCURS TO A SIGNIFICANT DEGREE IT CAN BE
FATAL. MANAGEMENT OF INSULATION EXPOSURES SHOULD INCLUDE
EXTENDED OBSERVATION AND RESPIRATORY SUPPORT AS INDICATED BY
THE PATIENTS PROGRESS AND PRESENTATION. DUE TO THE UNDERLYING
RESPIRATORY IRRITATION IN CASES OF PULMONARY EDEMA THE USE OF
HYPERBARIC OXYGEN MAY, HOWEVER BE CONTRAINDICATED AND ITS USE
MUST BE CAREFULLY CONSIDERED. OTHER SEQUELAE OF INHALATION MAY
INCLUDE TIGHTNESS OF THE CHEST, HEADACHE, NAUSEA, VOMITING,
DIARRHEA, VERTIGO, WEAKNESS, DROWSINESS, CYANOSIS,
INCOORDINATION, CONVULSION AND/OR COMA DEPENDANT ON THE DEGREE
OF EXPOSURE. ANY OF THESE MAY BE DELAYED IN ONSET. MANAGEMENT IS
SUPPORTIVE AND SYMPTOMATIC.

ACCEPTED
VERIFICATION
IN DEPT. OF HEALTH
OCT 12 1995

LIQUID CARBONIC: Medifume Plus™ VOL I

11491-10

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING. THIS PRODUCT IS FOR INDUSTRIAL USE ONLY IN COMMERCIAL GAS STERILIZERS DESIGNED FOR USE WITH ETHYLENE OXIDE. THIS PRODUCT SHOULD BE USED ONLY BY PERSONS EXPERIENCED IN STERILIZATION OR FUMIGATION USING THIS MIXTURE, OR ARE UNDER THE DIRECT SUPERVISION OF PERSONS WHO ARE EXPERIENCED IN STERILIZATION OR FUMIGATION USING THIS MIXTURE. USE THIS PRODUCT IN ACCORDANCE WITH THE STERILIZER MANUFACTURER'S OPERATING INSTRUCTIONS FOR THE PARTICULAR STERILIZER AND ONLY AFTER READING LIQUID CARBONIC'S MATERIAL SAFETY DATA SHEET FOR MEDIFUME PLUSTM. ALWAYS CHECK FOR LEAKS OF THE CYLINDER VALVE BEFORE MOVING THE CYLINDER INTO THE STERILIZING FACILITY. DO NOT USE ANY OIL OR GREASE ON VALVES, REGULATOR OR FITTINGS FOR DELIVERY OF THIS GAS INTO THE STERILIZER. BACK FLOW INTO THE CONTAINER COULD CAUSE CONTAINER RUPTURE. ALWAYS USE A BACK FLOW PREVENTING DEVICE IN THE CONTAINER DISCHARGE PIPING. CHECK FOR LEAKS AFTER CONNECTING CYLINDER TO THE STERILIZER INLET PIPING. CONTAINER TEMPERATURE SHOULD NOT EXCEED 125°F (52°C) DURING STORAGE OR USE. KEEP PROTECTIVE CAP ON WHEN NOT IN USE. SECURE THE CYLINDER IN AN UPRIGHT POSITION DURING STORAGE AND USE. USE ONLY IN CLOSED SYSTEMS WHICH ARE RATED FOR THE CYLINDER PRESSURE (2000 PSIG). CLOSE VALVE AFTER EACH USE AND WHEN EMPTY.

ITEMS TO BE STERILIZED SHOULD BE THOROUGHLY CLEANED OF SOIL BEFORE BEING STERILIZED. THIS PRODUCT, WHEN USED IN ACCORDANCE WITH THE STERILIZER MANUFACTURER'S OPERATING PROCEDURES, WILL STERILIZE AND/OR REDUCE MICROBIAL LOAD ON MEDICAL AND LABORATORY ITEMS, PHARMACEUTICALS, ASEPTIC PACKAGING, WHOLE AND GROUND SPICES, AND OTHER PROCESSED NATURAL SEASONING MATERIALS EXCEPT MIXTURES TO WHICH SALT HAS BEEN ADDED (SEE 40 CFR 185.2850) AND COSMETICS.

LIQUID CARBONIC: Medifume PlusTM VOL IACCEPTED
VOL. 10, 1995
FEDERAL BUREAU OF INVESTIGATION

OCT 12 1995

11441-10