

## BACK PANEL

### BENZALKONIUM CHLORIDE U.S.P.

TINCTURE  
SOLUTION 1:750  
FOR USE IN HOSPITALS

Thermometer Kits: To disinfect thermometer, use undiluted. Allow thermometer to remain immersed at least ten minutes. Rinse thermometer with potable water before use.

**WARNING — FLAMMABLE! KEEP AWAY FROM HEAT OR OPEN FLAME**

#### CAUTION

**HARMFUL IF  
SWALLOWED**

Mixing benzalkonium chloride solution with ordinary soaps and with anionic detergents may decrease or destroy the bactericidal activity of the solution. Discard used package in a safe place.

## FRONT PANEL

TEAR HERE

### BENZALKONIUM CHLORIDE U.S.P.

TINCTURE  
SOLUTION 1:750

#### CAUTION

KEEP OUT OF REACH OF CHILDREN. SEE OTHER PRECAUTIONS ON THE BACK PANEL.

$\frac{2}{3}$  Fluid Oz.

#### ACTIVE INGREDIENTS:

ALKYL (50% C<sub>12</sub>, 30% C<sub>14</sub>,  
17% C<sub>16</sub>, 3% C<sub>18</sub>) DIMETHYL  
BENZYL AMMONIUM  
CHLORIDE ..... 0.13%  
ISOPROPANOL ..... 40.00%  
INERT INGREDIENT: .... 59.87%

EPA Reg. No.

11563-2

Manufactured by

**Liquipak Corporation**

ALMA, MICHIGAN 48801

**ACCEPTED**

MAY 7 1975

UNDER THE FEDERAL INSECTICIDE  
FUNGICIDE AND RODENTICIDE ACT  
FOR ECONOMIC POISON REGISTERED  
UNDER NO. 11563-2