


733 W. BROAD STREET
EMMAUS, PA. 18049


STERILIZING GAS MIXTURE
FOR USE IN COMMERCIAL GAS STERILIZERS

ACTIVE INGREDIENT
Ethylene Oxide 10 wt %

INERT INGREDIENT
Carbon Dioxide 90 wt %

NET CONTENTS 1 lb.
EPA REG. NO. 28472-6



 *McPherson and Chemicals*
SPECIALTY GAS DEPARTMENT
Homestead Facility
P.O. Box 351 RD 2, Tompaugh, Pa. 18252

WARNING!
KEEP OUT OF REACH OF CHILDREN
Contents under pressure. Do not use or store near heat or open flame. Exposure to temperatures above 125°F may cause bursting. Breathing of vapor is harmful. Exposure to eyes, skin and clothing may cause severe irritation. Avoid contact. In case of contact, immediately remove all contaminated clothing and flush with plenty of water for at least fifteen minutes. For eyes get medical attention. Instructions for use and other information are in Sterilizing Gas Brochure.

ACCEPTED
FEB 21 1974
EPA REG. NO. 28472-6


733 W. BROAD STREET
EMMAUS, PA. 18049

STERILIZING GAS MIXTURE
FOR USE IN COMMERCIAL GAS STERILIZERS

ACTIVE INGREDIENT
Ethylene Oxide 10 wt %

INERT INGREDIENT
Carbon Dioxide 90 wt %

NET CONTENTS 1 lb.
EPA REG. NO. 28472-6



Star Products and Chemicals
SPECIALTY GAS DEPARTMENT
Hometown Facility
P.O. Box 351 RD 2 Tamaqua Pa 18252

WARNING!
KEEP OUT OF REACH OF CHILDREN
Contents under pressure. Do not use or store near heat or open flame. Exposure to temperatures above 125°F may cause bursting. Breathing of vapor is harmful. Exposure to eyes, skin and clothing may cause severe irritation. Avoid contact. In case of contact, immediately remove all contaminated clothing and flush with plenty of water for at least fifteen minutes. For eyes get medical attention. Instructions for use and other information are in Sterilizing Gas Brochure.

[Faint, illegible handwritten text]

DO NOT REMOVE THIS TAG

WARNING

ACCEPTED
FEB 01 1974
28472 6

STERILIZING GAS MIXTURE

FOR USE IN COMMERCIAL
GAS STERILIZERS

Gross Wt _____ lb
Tare Wt _____ lb
Net Wt _____ lb

Batch No. _____

Ethylene Oxide Is Subject to Polymerization

Filling Date _____
Use Before _____

BEFORE USING, READ STERILIZING
GAS BROCHURE _____



Specialty Gas Department
Hometown Facility
P. O. Box 351
RD 2, Tamaqua, Pennsylvania 18252
Telephone: (717) 467-2981

DO NOT REMOVE THIS TAG

WARNING

GAS AND LIQUID UNDER PRESSURE

Protect cylinder against damage. Secure cylinder upright when discharging. Cylinder valve opens by turning to the left. Always open the valve wide. Close valve when not in use. Do not use or store near heat or open flame. Use and store in a cool, dry, well ventilated area. Never allow cylinder temperature to exceed 125°F. Do not empty cylinder completely. Leave a positive pressure remaining in the cylinder. Never recharge cylinder.

 *Air Products and Chemicals*
INC.