

TEAR AT NOTCH

CLINIPAD®

**Benzalkonium
Chloride**

U.S.P.

Aqueous Solution 1:750

FOR HOSPITAL USE

CAUTION

Keep Out of Reach of Children

.67 Fluid Oz.

**E.P.A. REG. NO.
31973-1**

Sold by
**THE CLINIPAD CORPORATION
STAMFORD, CONNECTICUT 06907 U.S.A.**

TEAR AT NOTCH

CLINIPAD®

**Benzalkonium
Chloride
U.S.P.**

ACTIVE INGREDIENTS:

Alkyl* dimethyl benzyl ammonium
chlorides 0.13%

INERT INGREDIENT:

micropore filtered water 99.87%

* (C12 67%, C14 25%, C16 7%, C8-C10-C18
1%) dimethyl benzyl ammonium chloride.

DIRECTIONS FOR USE: Product to
be used undiluted in disinfection
and storage of precleaned thermo-
meters and instruments. (A sub-
mersion time of 10 min. should be
allowed)

DISCARD IN A SAFE PLACE

CAUTION: mixing benzalkonium
chloride solution with ordinary
soaps or with anionic detergents
may decrease or destroy the bac-
tericidal activity of the solution.

ACCEPTED WITH COMMENTS