

ETHYLENE OXIDE AND REFRIGERANT STERILIZING GAS

FOR USE IN COMMERCIAL

GAS STERILIZERS

**ETHYLENE OXIDE AND
DICHLORODIFLUOROMETHANE
MIXTURE**

ACTIVE INGREDIENT:
ETHYLENE OXIDE 12% BY WEIGHT
INERT INGREDIENT:
DICHLORODIFLUOROMETHANE
88% BY WT.
TOTAL CONTENTS _____ LBS. NET
E. P. A. REG. NO. 32462-1

WARNING! CONTENTS UNDER PRESSURE!

1. Do not puncture. 2. Do not use or store near heat or open flame. Exposures above 130°F may cause bursting. Never throw containers into fire or incinerator. 3. Do not use oil or grease on valve, regulator, fittings. 4. Do not recharge without written consent of Liquid Air Inc.

WARNING! DO NOT IGNORE!

KEEP OUT OF REACH OF CHILDREN

1. Vapor harmful. Irritating to nose, throat, eyes and skin. Avoid breathing vapors. Avoid contact with eyes, skin and clothing. In case of contact immediately remove all contaminated clothing, and flush skin or eyes with plenty of water for at least 15 minutes; for eyes, get medical attention. Wash contaminated clothing before reuse. 2. Intended for gas sterilization only. 3. Ethylene oxide is subject to polymerization; use within 2 months after receipt. 4. Instructions for use and other information are in Gas Bulletin SG10.



LIQUID AIR INC.

Divisions:
AMERICAN CRYOGENICS
INDUSTRIAL AIR PRODUCTS

(Non Life Supporting)

ETHYLENE OXIDE AND REFRIGERANT STERILIZING GAS

FOR USE IN COMMERCIAL

GAS STERILIZERS

**ETHYLENE OXIDE AND
DICHLORODIFLUOROMETHANE
MIXTURE**

**ACTIVE INGREDIENT:
ETHYLENE OXIDE 12% BY WEIGHT
INERT INGREDIENT:
DICHLORODIFLUOROMETHANE
88% BY WT.
TOTAL CONTENTS _____ LBS. NET
E. P. A. REG. NO. 32462-1**


WARNING! CONTENTS UNDER PRESSURE!

1. Do not puncture. 2. Do not use or store near heat or open flame. Exposures above 130°F may cause bursting. Never throw containers into fire or incinerator. 3. Do not use oil or grease on valve, regulator, fittings. 4. Do not recharge without written consent of Liquid Air Inc.

WARNING! DO NOT IGNORE!

KEEP OUT OF REACH OF CHILDREN

1. Vapor harmful. Irritating to nose, throat, eyes and skin. Avoid breathing vapors. Avoid contact with eyes, skin and clothing. In case of contact immediately remove all contaminated clothing, and flush skin or eyes with plenty of water for at least 15 minutes; for eyes, get medical attention. Wash contaminated clothing before reuse. 2. Intended for gas sterilization only. 3. Ethylene oxide is subject to polymerization; use within 2 months after receipt. 4. Instructions for use and other information are in Gas Bulletin SG10.



LIQUID AIR INC.

Divisions:
**AMERICAN CRYOGENICS
INDUSTRIAL AIR PRODUCTS**

(Non Life Supporting)