



SUPER WEX-CEL

CONCENTRATED GERMICIDAL DETERGENT

CLEANS ——— DISINFECTS ——— DEODORIZES

GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL,
FUNGICIDAL and VIRUCIDAL*

in soft or hard water. (A.O.A.C. 400 ppm Hardness as Calcium Carbonate)

ACTIVE INGREDIENTS:

Sodium ortho-Phenylphenate 13.27%, Sodium ortho-Benzyl-para-Chlorophenate
13.02%, Sodium Dodecyl Benzene Sulfonate 3.06%, Tetrasodium Ethylenediamine
Tetraacetate 2.42%.

INERT INGREDIENTS

68.23%. Includes biodegradable detergents and no phosphates.

EPA REGISTRATION NO. 34810-7

DANGER: KEEP OUT OF REACH OF CHILDREN.

See rear panel for additional precautionary and first aid statements.

Wexford Labs, Inc.

1035 South Vandeventer

St. Louis, Missouri 63110

ACCEPTED

SEP 11 1975

UNDER THE FEDERAL INSECTICIDE
FUNGICIDE AND RODENTICIDE ACT
FOR ECONOMIC POISON REGISTER-
ED UNDER NO. 34810-7

FOR INDUSTRIAL AND INSTITUTIONAL USE
CONTENTS ONE U. S. GALLON

SUPER WEX-CEL

CONCENTRATED GERMICIDAL DETERGENT

SUPER WEX-CEL is used for cleaning, disinfecting and deodorizing floors, walls and other similar hard surfaces. Especially suited for use in Operating Rooms, Intensive Care, Nurseries, Emergency Areas, Patient Rooms and Rest Rooms.

When used as directed SUPER WEX-CEL is effective against Staphylococcus aureus, Salmonella choleraesuis, Pseudomonas aeruginosa, pathogenic fungi, and Lipophilic viruses.

DIRECTIONS

For all disinfection, including hospital and health care facilities, use SUPER WEX-CEL at a dilution of 1/2 ounce per gallon of water. *When used at a dilution of 1/2 ounce per gallon of water on environmental inanimate surfaces, SUPER WEX-CEL is effective against Influnza A₂, Herpes Simplex, Adenovirus Type 5, and Vaccinia viruses.

SUPER WEX-CEL may be applied with a cloth, mop, or brush, or by immersing equipment to be disinfected.

When disinfecting surgical instruments with SUPER WEX-CEL, these instruments must be precleaned to remove blood and serious exudates prior to immersion in, or treatment with SUPER WEX-CEL.

Rinse empty container thoroughly with water and detergent before discarding or returning drum to drum conditioner.

SUPER WEX-CEL is a synthetic germicidal detergent with broad spectrum kill of both gram positive and gram negative microorganisms. It is an excellent cleaner which dilutes clear (1:256) in soft or hard water up to 400 ppm water hardness. Contains biodegradable detergents and is phosphate free. Do not mix with other cleaners.

GERMICIDAL QUALIFICATIONS

(AOAC Use dilution Confirmation) 1:256

TUBERCULOCIDAL ACTIVITY TEST

(AOAC, M. Tuberculosis) 1:256

FUNGICIDAL TEST

(AOAC, T. Interdigitale) 1:256

DANGER

KEEP OUT OF REACH OF CHILDREN Corrosive Causes eye damage and severe skin irritation. Avoid breathing spray mist. Do not get in eyes, on skin, or on clothing. Protect eyes and skin when handling concentrate. Harmful or fatal if swallowed. Avoid contamination of food.

FIRST AID

In case of contact, immediately flush eyes or skin with plenty of water for a least 15 minutes. For eyes call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol and call a physician.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

Wexford Labs, Inc.

EPA ESTABLISHMENT NUMBER 34810-MO-1

ACCEPTED

SEP 11 1975

UNDER THE FEDERAL INSECTICIDE
FUNGICIDE AND RODENTICIDE ACT