

FOR INDUSTRIAL AND MEDICAL USE ONLY IN COMMERCIAL GAS
Cylinders in Accordance with the Manufacturer's Recommendations

NET WT. 135 LBS.

ACTIVE INGREDIENT:

ETHYLENE OXIDE

12

INERT INGREDIENTS:

DICHLORODIFLUOROMETHANE 88%

TOTAL CONTENTS 135 LBS. NET

WARNING

KEEP OUT OF
REACH OF CHILDREN

EPA Reg. No. 35084-1

(ACTIVE INGREDIENT)

ETHYLENE OXIDE

12%

(INERT INGREDIENT)

DICHLORODIFLUOROMETHANE 88%

ACCEPTED

TOTAL CONTENTS 135 LBS. NET

001 NO 1973

UNDER THE FEDERAL INSECTICIDE
FUNGICIDE AND RODENTICIDE ACT
FOR ECONOMIC POISON REGISTRY
REGISTRY NO. 35084-1

WARNING

**KEEP OUT OF
REACH OF CHILDREN**

1. Exposure to eyes, skin or clothing may cause severe irritation. Avoid contact. In case of contact, remove all contaminated clothes immediately and flush skin or eyes with plenty of water for at least 15 minutes; for eyes, get medical attention.
2. Use in accordance with sterilizer manufacturer's operating instructions.
3. Ethylene oxide is subject to polymerization; use within 2 months after receipt.
4. Close valve when not in use. Handle cylinder with care. Contents under pressure.
5. Use and store in cool, dry, well ventilated area. Keep cylinder temperature from exceeding 125°F.

EPA EST. No. 35084-NJ-1

**MG SCIENTIFIC GASES
DIV. OF MG TECHNICAL PRODUCTS, INC
1100 HARRISON AVE., KEARNY, N. J. 07029**