

**METRICIDE 28™
LONG-LIFE
ACTIVATED
DIALDEHYDE**

Activator Vial

For use with
**METRICIDE 28™
LONG-LIFE
ACTIVATED
DIALDEHYDE
FLUID**

Directions: Add contents of this vial to five gallons of fluid. Shake. Activated solution will be green.

Danger: Keep out of reach of children. Contact with eyes causes damage. In case of contact with eyes, flush with water immediately and get medical attention. For skin contact, flush thoroughly with water.

CONTAINER DISPOSAL


Rinse empty container thoroughly with water and discard it. Do not reuse empty container. Wrap container and put into trash collection.

Net weight: 1.0000 lb.

METREX RESEARCH CORP.

Littleton, Colorado

303-741-1111



BEST DOCUMENT AVAILABLE