

10031

52252-1

124

renalin™

DIALYZER REPROCESSING CONCENTRATE

COLD STERILANT FOR DIALYZER REPROCESSING

READ RENALIN INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.

These contents must be diluted with AAMI water prior to use.

Contains two containers. Each of the two containers has two liters of Renalin™ Concentrate

STORAGE

- | | |
|---|---|
| <ol style="list-style-type: none"> 1. Renalin™ is a strong oxidizer. 2. Store in shipping carton. 3. Do not expose to direct sunlight. 4. Maintain temperature below 75°F (24°C). 5. Avoid contact with skin and eyes. Wear rubber gloves and goggles. 6. Avoid contact with combustible materials. 7. Avoid contamination from any source, including metals, dust, etc. | <p>Such contamination may cause rapid decomposition, generation of large quantities of oxygen gas and high pressures.</p> <ol style="list-style-type: none"> 8. Store in original closed container—NEVER TAMPER WITH VENT. 9. If spilled, flush away with large quantities of water. 10. Use AAMI quality water in all dialyzer reprocessing steps when using Renalin.™ |
|---|---|

DISPOSAL

1. Do not contaminate water, food or feed by storage or disposal.
2. Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.
3. Triple rinse empty container with water. Then offer container for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities, dispose of by burning. If burned, stay clear of the smoke.

Active ingredients

Hydrogen peroxide 24.0%
 Peroxyacetic acid 4.0%

Inert ingredients 72.0%

Total 100.0%

NET CONTENTS: 1.06 U.S. gallons (4 liters)


LOT NO. _____

CONCENTRATE
 EXPIRATION DATE _____

For chemical emergency, spill, leak, fire, exposure and accident, call Chemtrec, day or night

(800) 424-9300.

In District of Columbia or Outside Continental U.S.A. Call
 (202) 483-7816



	<p>RENAL SYSTEMS, INC. 14905 28th Avenue North Minneapolis, Minnesota 55441 U.S.A. Telephone (612) 553-3300 • (800) 328-3340 Telex 29 0825 RSIMPLS PLOH</p>
--	--

ACCEPTED
SEP 25 1985
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 50252-1

Label

Renalin™

Concentrate

2.5 gallon container

renalin™

DIALYZER REPROCESSING CONCENTRATE

COLD STERILANT FOR DIALYZER REPROCESSING

READ RENALIN INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.

These contents must be diluted with AAMI water prior to use.
 Contains two containers. Each of the two containers has two liters of Renalin™ Concentrate

STORAGE

1. Renalin™ is a strong oxidizer.
2. Store in shipping carton.
3. Do not expose to direct sunlight.
4. Maintain temperature below 75°F (24°C).
5. Avoid contact with skin and eyes. Wear rubber gloves and goggles.
6. Avoid contact with combustible materials.
7. Avoid contamination from any source, including metals, dust, etc.
8. Store in original closed container—**NEVER TAMPER WITH VENT.**
9. If spilled, flush away with large quantities of water.
10. Use AAMI quality water in all dialyzer reprocessing steps when using Renalin.™

DISPOSAL

1. Do not contaminate water, food or feed by storage or disposal.
2. Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.
3. Triple rinse empty container with water. Then offer container for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities, dispose of by burning. If burned, stay clear of the smoke.

ACCEPTED
 SEP 25 1985
 Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No. 222-000000

Active ingredients	
Hydrogen peroxide	24.0%
Peracetic acid	4.0%
Inert ingredients	72.0%
Total	100.0%

NET CONTENTS: 1.06 U.S. gallons (4 liters)

LOT NO. _____

CONCENTRATE

EXPIRATION DATE _____

BEST AVAILABLE COPY

For chemical emergency, spill, leak, fire, exposure and accident, call Chemtrec, day or night

(800) 424-9300

COLD STERILANT FOR DIALYZER REPROCESSING

READ RENALIN INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.

These contents must be diluted with AAMI water prior to use.

Contains two containers. Each of the two containers has two liters of Renalin™ Concentrate

STORAGE


1. Renalin™ is a strong oxidizer.
2. Store in shipping carton.
3. Do not expose to direct sunlight.
4. Maintain temperature below 75°F (24°C).
5. Avoid contact with skin and eyes. Wear rubber gloves and goggles.
6. Avoid contact with combustible materials.
7. Avoid contamination from any source, including metals, dust, etc.

Such contamination may cause rapid decomposition, generation of large quantities of oxygen gas and high pressures.

8. Store in original closed container—**NEVER TAMPER WITH VENT.**
9. If spilled, flush away with large quantities of water.
10. Use AAMI quality water in all dialyzer reprocessing steps when using Renalin.™

DISPOSAL

1. Do not contaminate water, food or feed by storage or disposal.
2. Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.
3. Triple rinse empty container with water. Then offer container for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities, dispose of by burning. If burned, stay clear of the smoke.



Active ingredients

Hydrogen peroxide 24.0%

Peracetic acid 4.0%

Other ingredients 72.0%

Total 100.0%

NET CONTENTS: 1.06 U.S. gallons (4 liters)

LOT NO. _____

**CONCENTRATE
EXPIRATION DATE** _____

BEST AVAILABLE COPY

For chemical emergency, spill, leak, fire, exposure and accident, call Chemtrec, day or night

(800) 424-9300.

In District of Columbia or Outside Continental U.S.A. Call:
(202) 483-7616


Shipping Carton Label

Renalin™

Concentrate

Twin Pack

(2 - 2.5 gallon containers)



WARNINGS

DANGER—KEEP OUT OF REACH OF CHILDREN.

DANGER—Corrosive. Causes eye damage and skin irritation. Do not get in eyes, on skin or on clothing. Wear safety glasses and rubber gloves when handling. Harmful if swallowed. Wash thoroughly after handling.

INSTRUCTIONS FOR USE

1. It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
2. Read Renalin™ Instructions for Use and Technical Notes before using this product.
3. Use only with the Renatron® Dialyzer Reprocessing System, which has been modified for use with Renalin™.
4. Dilute Renalin™ Concentrate by filling original container to the shoulder with AAMI quality water.
5. Connect to Renatron® uptake hose and follow machine instructions.

renalin™

DIALYZER REPROCESSING CONCENTRATE

COLD STERILANT FOR DIALYZER REPROCESSING

READ RENALIN™ INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.
This product must be diluted with AAMI water prior to use.

Active Ingredients	
Hydrogen peroxide	24.0%
Peroxyacetic acid	4.0%
Inert Ingredients	72.0%
Total	100.0%

Concentrate
Expiration
Date

DILUTED EXPIRATION DATE

NET CONTENTS:
0.53 U.S. gallons (2 liters)
EPA Reg. No.: 52252-1
EPA Est. No.: 53263-1A-01
LOT NO.

DANGER—KEEP OUT OF REACH OF CHILDREN.

In case of contact with skin or eyes, flush with large amounts of water for at least fifteen (15) minutes. For eyes, get prompt medical attention.

If swallowed, drink water immediately to dilute, do not attempt to cause vomiting. Call physician immediately.

See Renalin™ Instructions for Use and Technical Notes (PN 50058-000) as well as side panel for additional precautionary statements.

 **renal
systems
inc.**

RENAL SYSTEMS INC.
145 N. 28th Avenue, North
Minneapolis, Minnesota 55412 U.S.A.
Telephone: 612-553-3300 • FAX: 612-553-3311
Telex: 29 0625 DS MP: S PLOH

67013-0000

BEST AVAILABLE COPY

ACCEPTED

SEP 25 1985

Under the Federal Insecticide,
Fungicide, and Herbicide Act,
as amended, for the pesticide
registered under
EPA Reg. No. **52252-1**

renalin™

DIALYZER REPROCESSING CONCENTRATE

COLD STERILANT FOR DIALYZER REPROCESSING

READ RENALIN™ INSTRUCTIONS FOR USE
AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.

This product must be diluted with AAMI water prior to use.

Ingredients
Hydrogen peroxide . . . 24.0%
Acetic acid . . . 4.0%
Other ingredients . . . 72.0%
Total . . . 100.0%

DILUTED EXPIRATION DATE

NET CONTENTS
0.53 U.S. gallons (2 liters)
EPA Reg. No.: 50262-1
EPA Prod. No.: 50263-1A-01

LOT NO.

DANGER—KEEP OUT OF REACH OF CHILDREN.

In case of contact with skin or eyes, flush with large amounts of water for at least fifteen (15) minutes. For eyes, get prompt medical attention.

If swallowed, drink water immediately to dilute, do not attempt to cause vomiting. Call physician immediately.

See Renalin™ Instructions for Use and Technical Notes (PN 50058-000) as well as side panel for additional precautionary statements.

STORAGE

1. Store in shipping carton.
2. Do not expose to direct sunlight.
3. Maintain temperature below 75°F (24°C).
4. Avoid contact with combustible materials.
5. Avoid contamination from any source, including metals, dust, etc. Such contamination may cause rapid decomposition, generation of large quantities of oxygen gas and high pressures.
6. Store in original closed container—**NEVER TAMPER WITH VENT.**
7. Use AAMI quality water in all dialyzer reprocessing steps when using Renalin™.
8. After dilution, the resultant solution must be used within seven (7) days. Fill in the expiration date of diluted solution at time of dilution in the space provided on this label. The expiration date of the diluted Renalin™ is seven days after dilution.
9. For chemical emergency, spill, leak, fire, exposure and accident, call Chemtrec, day or night (800) 424-9300, (202) 483-7616.

DISPOSAL

1. Do not contaminate water, food or feed by storage or disposal.
2. Wastes resulting from the use of this product may be disposed of on-site or at an approved waste disposal facility.
3. Triple rinse empty container with water. Then offer for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities dispose of by burning. If burned, stay clear of the smoke.

thoroughly after handling.

DIALYZER REPROCESSING CONCENTRATE

INSTRUCTIONS FOR USE

1. It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
2. Read Renalin™ Instructions for Use and Technical Notes before using this product.
3. Use only with the Renatron® Dialyzer Reprocessing System, which has been modified for use with Renalin™.
4. Dilute Renalin™ Concentrate by filling original container to the shoulder with AAMI quality water.
5. Connect to Renatron® uptake hose and follow machine instructions.

COLD STERILANT FOR DIALYZER REPROCESSING

READ RENALIN™ INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.
This product must be diluted with AAMI water prior to use.

Active Ingredients	
Hydrogen peroxide	24.0%
Peroxyacetic acid	4.0%
Inert Ingredients	
	72.0%
Total	100.0%

DILUTED EXPIRATION DATE

SEP 25 1985

NET CONTENTS
 0.53 U.S. gallons (2 lit)
 EPA Reg. No.: 50052-1
 EPA Est. No.: 50053-1

Concentrate Expiration Date

DANGER—KEEP OUT OF REACH OF CHILDREN.

In case of contact with skin or eyes, flush with large amounts of water for at least fifteen (15) minutes. For eyes, get prompt medical attention.

If swallowed, drink water immediately to dilute, do not attempt to cause vomiting. Call physician immediately.

See Renalin™ Instructions for Use and Technical Notes (PN 50058-000) as well as side panel for additional precautionary statements.

ACCEPTED

SEP 25 1985

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 50052-1



67013-000D

RENAL SYSTEMS, INC.
 14905 28th Avenue North
 Minneapolis, Minnesota 55441 U.S.A.
 Telephone (612) 553 3300 • (800) 328-3333
 Telex 29 0825 RSIMPLS PLOH

Label

Renalin™

Concentrate

for use with HPR™ Machine

32 oz. Container, Containing 6 fl. oz.

ACCEPTED
SEP 25 1985
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 52252-1

))))
))))
))))
))))
))))
))))
))))
))))
))))
))))

10.7.19

ACCEPTED
SEP 25 1985
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under 52252-1
EPA Reg. No.

FILL TO THIS LINE

renalin™

DIALYZER REPROCESSING CONCENTRATE

COLD STERILANT FOR DIALYZER REPROCESSING
FOR USE WITH HPR™
DIALYZER REPROCESSING SYSTEM ONLY.

READ RENALIN™ INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.
This product must be diluted with AAMI water prior to use.

WARNINGS

DANGER—KEEP OUT OF REACH OF CHILDREN.

DANGER—Corrosive. Causes eye damage and skin irritation. Do not get in eyes, or skin or on clothing. Wear safety glasses and rubber gloves when handling. Harmful if swallowed. Wash thoroughly after handling.

INSTRUCTIONS FOR USE

1. It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
2. If concentrated Renalin™ in unopened bottle is not within 1/16 inch of "Concentrate" line shown below, DO NOT USE.
3. Read Renalin™ Instructions for Use and Technical Notes before using this product.
4. Use only with the HPR™ Dialyzer Reprocessing System, which has been modified for use with Renalin™.
5. Dilute Renalin™ Concentrate by filling this container to the "Fill" line with AAMI quality water.
6. Connect to HPR™ uptake hose and follow machine instructions

Active ingredients	
Hydrogen peroxide	24.0%
Peroxyacetic acid	4.0%
Inert ingredients	72.0%
Total	100.0%

DILUTED EXPIRATION DATE

NET CONTENTS:
6 U.S. fluid oz.
(177.4 milliliters)
EPA Reg. No.: 52252-1
EPA Est. No.: 53263-IA-01

**Concentrate
Expiration
Date**

LOT NO.


DANGER—KEEP OUT OF REACH OF CHILDREN.
In case of contact with skin or eyes, flush with large amounts of water for at least fifteen (15) minutes. For eyes, get prompt medical attention.
If swallowed, drink water immediately to dilute, do not attempt to cause vomiting. Call physician immediately.
See Renalin™ Instructions for Use and Technical Notes (PN 50058-000) as well as side panel for additional precautionary statements.

CONCENTRATE



renalin™

DIALYZER REPROCESSING CONCENTRATE



FILL TO THIS LINE

COLD STERILANT FOR DIALYZER REPROCESSING

FOR USE WITH HPR™
DIALYZER REPROCESSING SYSTEM ONLY.

READ RENALIN™ INSTRUCTIONS FOR USE AND
TECHNICAL NOTES BEFORE USING THIS PRODUCT.
This product must be diluted with AAMI water prior to use.

STORAGE

1. Store in shipping carton.
 2. Do not expose to direct sunlight.
 3. Maintain temperature below 75°F (24°C).
 4. Avoid contact with combustible materials.
 5. Avoid contamination from any source, including metals, dust, etc. Such contamination may cause rapid decomposition, generation of large quantities of oxygen gas at high pressures.
 6. Store in original closed container—**NEVER TAMPER WITH VENT.** Use AAMI quality water in all dialyzer reprocessing steps when using Renalin™.
 7. After dilution, the resultant solution must be used within seven (7) days. Fill in the expiration date of diluted solution at time of dilution in the space provided on this label. The expiration date of the diluted Renalin™ is seven days after dilution.
- For chemical emergency, spill, leak, fire, exposure and accident, call Chemtrec, day or night (800) 424-9300, (202) 483-7616.

DISPOSAL

1. Do not contaminate water, food or feed by storage or disposal.
2. Wastes resulting from the use of this product may be disposed of on-site or at an approved waste disposal facility.
3. Triple rinse empty container with water. Then offer for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities dispose of by burning. If burned, stay clear of the smoke.

WARNINGS

DANGER—KEEP OUT OF REACH OF CHILDREN.

DANGER—Corrosive. Causes eye damage and skin irritation. Do not get in eyes, or skin or on clothing. Wear safety glasses and rubber gloves when handling. Harmful if swallowed. Wash thoroughly after handling.

INSTRUCTIONS FOR USE

1. It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
2. If concentrated Renalin™ in unopened bottle is not within 1/16 inch of "Concentrate" line shown below, DO NOT USE.
3. Read Renalin™ Instructions for Use and Technical Notes before using this product.
4. Use only with the HPR™ Dialyzer Reprocessing System, which has been modified for use with Renalin™.
5. Dilute Renalin™ Concentrate by filling this container to the "Fill" line with AAMI quality water.
6. Connect to HPR™ uptake hose and follow machine instructions.

Active ingredients

Hydrogen peroxide	24.0%
Peroxyacetic acid	4.0%
Inert ingredients	72.0%
Total	100.0%

DILUTED EXPIRATION DATE

NET CONTENTS:
1.05 U.S. fluid oz.
(177.4 milliliters)
EPA Reg. No.: 52252-1
EPA Est. No.: 53263-IA-01

ACCEPTED
SEP 25 1985
Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 52252-1

Concentrate
Expiration
Date

DANGER—KEEP OUT OF REACH OF CHILDREN.

In case of contact with skin or eyes, flush with large amounts of water for at least fifteen (15) minutes. For eyes, get prompt medical attention.

If swallowed, drink water immediately to dilute, do not attempt to induce vomiting. Call physician immediately.

See Renalin™ Instructions for Use and Technical Notes (PN 50057-000) as well as side panel for additional precautionary statements.

CONCENTRATE

renal systems inc.

RENAL SYSTEMS, INC.
13205 28th Avenue North
Minneapolis, Minnesota 55417, U.S.A.
Telephone: 612-551-1100 • Telex: 254311
Telex: 254311 RENAL SYSTEMS, INC.


Shipping Carton

Label

Renalin™

Concentrate

For use with HPR™ Machine
12 - 32 oz. containers,
Containing 6 oz. each



2814

renalin™

DIALYZER REPROCESSING CONCENTRATE

COLD STERILANT FOR DIALYZER REPROCESSING

FOR USE WITH HPR™ DIALYZER REPROCESSING SYSTEM ONLY.

READ RENALIN INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.

Case contains 12 bottles. Each bottle contains 177.4 milliliters of Renalin™ concentrate.

These contents must be diluted with AAMI water prior to use.

STORAGE

1. Renalin™ is a strong oxidizer.
2. Store in shipping carton.
3. Do not expose to direct sunlight.
4. Maintain temperature below 75°F (24°C).
5. Avoid contact with skin and eyes. Wear rubber gloves and goggles.
6. Avoid contact with combustible materials.
7. Avoid contamination from any source, including metals, dust, etc. Such contamination may cause rapid decomposition, generation of large quantities of oxygen gas and high pressures.
8. Store in original closed container—**NEVER TAMPER WITH VENT.**
9. If spilled, flush away with large quantities of water.
10. Use AAMI quality water in all dialyzer reprocessing steps when using Renalin.

DISPOSAL

1. Do not contaminate water, food or feed by storage or disposal.
2. Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.
3. Triple rinse empty container with water. Then offer container for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities, dispose of by burning. If burned, stay clear of the smoke.

DANGER—KEEP OUT OF REACH OF CHILDREN

Active ingredients

Hydrogen peroxide 24.0%

Peroxyacetic acid 4.0%

Inert ingredients 72.0%

Total 100.0%

NET CONTENTS: 72 U.S. fluid oz. (2.13 liters)

LOT NO. _____

CONCENTRATE
EXPIRATION DATE _____

For chemical emergency, spill, leak, fire, exposure
and accident, call Chemtrac, day or night

(800) 424-9300.


In District of Columbia or Outside Continental U.S.A. Call:
(202) 483-7818

 **renal
systems
inc.**

67018-0008

RENAL SYSTEMS, INC.

14905 28th Avenue North
Minneapolis, Minnesota 55441 U.S.A.
Telephone (612) 553-3300 • (800) 328-3340
Telex 29-0825 RSIMPLS PLOH



renalin™

DIALYZER REPROCESSING CONCENTRATE

COLD STERILANT FOR DIALYZER REPROCESSING

FOR USE WITH HPR™ DIALYZER REPROCESSING SYSTEM ONLY.

READ RENALIN INSTRUCTIONS FOR USE AND TECHNICAL NOTES BEFORE USING THIS PRODUCT.

Case contains 12 bottles. Each bottle contains 177.4 milliliters of Renalin™ concentrate.
These contents must be diluted with AAMI water prior to use.

STORAGE

1. Renalin™ is a strong oxidizer.
2. Store in shipping carton.
3. Do not expose to direct sunlight.
4. Maintain temperature below 75°F (24°C).
5. Avoid contact with skin and eyes. Wear rubber gloves and goggles.
6. Avoid contact with combustible materials.
7. Avoid contamination from any source, including metals, dust, etc. Such contamination may cause rapid decomposition, generation of large quantities of oxygen gas and high pressures.
8. Store in original closed container—**NEVER TAMPER WITH VENT.**
9. If spilled, flush away with large quantities of water.
10. Use AAMI quality water in all dialyzer reprocessing steps when using Renalin.™

DISPOSAL

1. Do not contaminate water, food or feed by storage or disposal.
2. Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.
3. Triple rinse empty container with water. Then offer container for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities, dispose of by burning. If burned, stay clear of the smoke.

DANGER—KEEP OUT OF REACH OF CHILDREN

Active ingredients

Hydrogen peroxide 24.0%

Peroxyacetic acid 4.0%

Inert ingredients 72.0%

Total 100.0%

NET CONTENTS: 72 U.S. fluid oz. (2.13 liters)

LOT NO. _____

CONCENTRATE
EXPIRATION DATE _____

For chemical emergency, spill, leak, fire exposure and accident, call Chemtrec, day or night

(800) 424-9300

In District of Columbia or Outside Continental U.S.A. Call
(202) 483-7615



RENAL SYSTEMS, INC.

14905 28th Avenue North
 Minneapolis, Minnesota 55441 U.S.A.
 Telephone (612) 553-3300 • (800) 328-3341
 Telex 29-0825 RSIMPLS PLOH

ACCEPTED
 SEP 25 1985
 For the Federal Insecticide,
 Fungicide and Rodenticide Act
 amended, for the pesticide
 registered under
 A Reg No. 5225-1