

Kills HIV-1 on precleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling on inanimate surfaces/objects with blood or other body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

**Special Instructions for Cleaning and Decontamination
Against HIV of Surfaces or Objects Soiled with Blood or
Body Fluids.**

Personal Protection: When handling items soiled with blood or other body fluids, disposable latex gloves, gowns, surgical mask and eye protection are recommended.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

Contact Time: Cide-Swipes are effective as a virucide against HIV-1 (AIDS Virus) for a three minute contact time.

Disposal of Infectious Material: Blood or other body fluids should be autoclaved and disposed according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.