

Diacide Dialyzer is a single use, disposable dialyzer. It is designed for use with standard dialysis machines (Peritoneal type) and hemodialysis machines (Type A). It is used to dialyze patients with acute renal failure and chronic renal failure with acute physiological status. To assure complete dialysis, use Diacide Dialyzer with sterile physiological saline. To assure complete dialysis, use Diacide Dialyzer with a standard dialysis machine and standard dialysis solution.

DIACIDE Solution has been demonstrated to be compatible with hollow fiber membranes having Cuprophane, regenerated cellulose or cellulose acetate membranes and has been shown to be safe for dialyzer storage.

This is a single use product.

Diacide Dialyzer

disinfecting solution

Active ingredient
Gluconolactide ... 29%
Inert ingredients ... 71%
Total 100%

NOTE: This product is a concentrate and must be activated and diluted immediately prior to use. Activator must be added to fluid before this product is effective. See Directions for Activation and Dilution.

TM

MANUFACTURED BY:
GULFSTREAM MEDICAL INC.
550 Pleasant St.
Winthrop, MA 02152

Keep out of reach of children.
DANGER: See right side panel for precautions.

GULFSTREAM

CONTENTS 28 FL. OZ. (0.83 LITERS)
MADE IN U.S.A.
CONTROL NUMBER ON BOTTOM

NOTE DATES HERE.
MIXING DATE _____
EXPIRATION DATE _____

GULFSTREAM

BEST AVAILABLE COPY

ACCEPTED

MAR 27 1991

Registered under the Federal Insecticide,
Fungicide, and Rodenticide Act.
This product is registered under the pesticide
law under EPA Reg. No. 63720-2