

65787-3

06/25/2009

1/4



U.S. ENVIRONMENTAL PROTECTION AGENCY  
Office of Pesticide Programs  
Antimicrobials Division (7510P)  
1200 Pennsylvania Avenue NW  
Washington, D.C. 20460

EPA Reg.

Number:

65787-3

Date of

Issuance:

JUN 25 2009

Term of Issuance:

Conditional

Name of Pesticide Product:

Bleach-Wipe 1:50

## NOTICE OF PESTICIDE:

☒ Registration  
☐ Reregistration

(under FIFRA, as amended)

**FILE COPY**

Name and Address of Registrant (include ZIP Code):

Alcavis HDC, LLC  
8322 Helgerman Court  
Gaithersburg, MD 20877

Note: Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby registered/reregistered under the Federal Insecticide, Fungicide and Rodenticide Act.

Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product (OPP Decision No. D-405633) is registered in accordance with FIFRA sec 3(c)(7)(A) provided that you:

1. Submit and/or cite all data required for registration of your product under FIFRA sec. 3(c)(5) when the Agency requires all registrants of similar products to submit such data; and submit acceptable responses required for re-registration of your product under FIFRA section 4.

2. Make the labeling changes listed below before you release the product for shipment:

a. Revise the "EPA File Symbol to read, "EPA Reg. No. 65787-3."

Signature of Approving Official:

Emily Mitchell

Product Manager Team-32

Regulatory Management Branch II

Antimicrobials Division (7510P)

Date:

Page 1 of 2

b. The comparative claim "The right and ready solution" is unacceptable and has been removed from label.

c. Change "*Streptococcus pneumonia*" to read "*Streptococcus pneumoniae*."

d. These towelettes/wipes appear to be individually wrapped. However, since the towellettes/wipes are not individually sold the **Container Disposal** language must be revised as follows:

Dispose of towelette in trash after use. Do not flush in toilet. Non-refillable container. Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available discard container in trash.

e. A "Mat" is generally considered as a porous surface. Therefore, the Mat use site has been removed from the Hard, Nonporous Areas of Use.

3. Submit three (3) copies of your final printed labeling before distributing or selling the product bearing the revised labeling.

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records.

Sincerely,



Emily Mitchell  
Product Manager 32  
Regulatory Branch II  
Antimicrobials Division (7510P)

Enclosures: (Stamped Label)

[bracketed text may be added or deleted depending on market]

# ALCAVIS

## BLEACH-WIPE 1:50

For wiping hard non-porous environmental  
surfaces [and patient care equipment]

**KEEP OUT OF REACH OF CHILDREN**

**CAUTION**

[See back panel for additional precautionary statements]

Manufactured by:  
Alcavis HDC, LLC  
8322 Helgerman Court  
Gaithersburg, MD 20877  
www.AlcavisHDC.com  
1-800-726-2308

EPA Reg. No.: 65787-  
EPA Est. No.: 65787-MD-1

ACCEPTED  
with COMMENTS  
EPA Letter Dated:

JUN 25 2009

8" x 10" Towelette

Net Contents: [# towels, x oz. each]

Under the Federal Insecticide,  
Fungicide, and Rodenticide Act as  
amended, for the pesticide,  
registered under EPA Reg. No. 65787-3

### Active Ingredient:

Sodium Hypochlorite: ..... 00.10%  
Other Ingredients: ..... 99.90%  
Total: ..... 100.00%

1:50 Dilution of 5.25% Sodium Hypochlorite Bleach  
1,000 ppm Available Chlorine

### First Aid

If in  
Eyes:

- Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing.
- Call a poison control center or doctor for treatment advice.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment. Contact the Poison Control Center at 1-800-222-1222 for 24 hour emergency medical treatment information.

## BACTERICIDAL

Effective in the presence of a 5% organic soil load against:

*Staphylococcus aureus* Methicillin Resistant [MRSA] [ATCC 33592]  
*Enterococcus faecalis* Vancomycin Resistant [VRE] [ATCC 51299]  
*Streptococcus pneumonia* [ATCC 6303]  
*Streptococcus pyogenes*, [ATCC 12344]  
*Staphylococcus aureus* [ATCC 6538]  
*Pseudomonas aeruginosa* [ATCC 15442]  
*Salmonella enterica* [ATCC 10708]

**BLEACH-WIPE 1:50**

Wipe has 1,000 ppm available chlorine

3/4

## Precautionary Statements

### Hazards to humans and domestic animals

CAUTION: Causes minimal eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

### Areas of Use

For use on hard, nonporous surfaces including stainless steel, chrome, metal, glass, vinyl, plastic, laminate, linoleum, glazed ceramic tile, glazed porcelain, and sealed or painted wood. Areas of use include: [Medical, [Dental and Laboratory] Counters], [Exam Tables], [Carts], [Dialysis Machines], [Point of Care Equipment], [Telephones], [Sink Tops], Sinks, Counter Tops, Floors, [Bathrooms], [Tubs], [Handrails], [Door] [Knobs], [Lockers], and [Desks].

### Directions for Use

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

- [Use personal protective equipment (gloves) when disinfecting biohazardous surfaces] *[Hospital, Medical and Dental uses only]*
- For heavily soiled surfaces, a preliminary cleaning is required
- Open Wipe packet
- Remove premoistened towelette
- Apply towelette and wipe desired surface to be disinfected
- Treated surface must remain visibly wet for 10 minutes
- Use additional wipe(s) if needed to assure continuous 10 minute wet contact time
- Allow surface area to dry and discard used towelette

[This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.] *[Hospital, Medical and Dental uses only]*

### Storage and Disposal

Do not contaminate water, food, or feed by storage and disposal.

**Storage:** Do not store near heat or open flame.

**Product and Container Disposal:** Do not reuse towelette. Dispose of used towelette and packet in trash. Do not flush in toilet.

ACCEPTED  
with COMMENTS  
EPA Letter Dated:

JUN 25 2009

Under the Federal Insecticide,  
Fungicide, and Rodenticide Act as  
amended, for the pesticide,  
registered under EPA Reg. No. 65787-3