



U.S. ENVIRONMENTAL PROTECTION AGENCY  
 Office of Pesticide Programs  
 Antimicrobial Division (7510C)  
 401 "M" St., S.W.  
 Washington, D.C. 20460

EPA Reg. Number:

69251-1

Date of Issuance:

JAN 11 1999

NOTICE OF PESTICIDE:

Registration  
 Reregistration

(under FIFRA, as amended)

Term of Issuance:

Conditional

Name of Pesticide Product:

ViraVac In-Line  
 Disinfectant System

Name and Address of Registrant (include ZIP Code):

Viatro Corporation  
 6779 Engle Road  
 Cleveland, Ohio 44130

**Note:** Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby registered/reregistered under the Federal Insecticide, Fungicide and Rodenticide Act.

Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA sec. 3(c)(7)(A) provided that you:

1. Submit and/or cite all data required for registration/reregistration of your product under FIFRA sec. 3(c)(5) when the Agency requires all registrants of similar products to submit such data; and submit acceptable responses required for reregistration of your product under FIFRA section 4.

2. Make the following label change:

Revise the statement "For the Safe Disposal of Suction Canister Contents" to read "For the Disposal of Suction Canister Contents".

3. Submit two copies of the revised final printed label for the record.

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records.

Signature of Approving Official:

Robert S. Brennis, PM 32

Date:

JAN 11 1999

Sanitizes Up To  
**1500**  
CC

**PRECAUTIONARY STATEMENTS:** Hazards to Humans and Domestic Animals. **CORROSIVE:** Causes irreversible eye damage and

skin irritation. Harmful if swallowed. Do not get in eyes, on skin or on clothing. Avoid contamination of food. Avoid breathing fumes or vapors. Wash hands thoroughly with soap and water after handling. Remove and wash contaminated clothing before reuse. **Environmental Hazards:** This pesticide is toxic to fish and aquatic organisms.

**STATEMENT OF PRACTICAL TREATMENT: IF SWALLOWED:** Call a doctor or get medical attention. Do not induce vomiting or give anything by mouth to an unconscious person. Drink promptly, a large quantity of milk, eggwhites, gelatin solution, or if these are not available, drink large quantities of water. Avoid alcohol. **IF IN EYES:** Hold eyelids open and flush with a steady, gentle stream of water for 15 minutes. Get medical attention. **IF ON SKIN:** Wash with plenty of soap and water. Get medical attention. **NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed. For medical emergency, contact local Poison Control center. In case of spill, contact Chemtrec at 1-800-424-9300.

**ADDITIONAL PRECAUTIONARY STATEMENTS: FOR INSTITUTIONAL HEALTHCARE USE ONLY. WEAR HAND AND EYE PROTECTION WHEN HANDLING POTENTIALLY INFECTIOUS WASTE. NOT FOR INTERNAL USE. FOR SINGLE USE ONLY.**



### Liquid Medical Waste Sanitizer

For The Safe Disposal Of Suction Canister Contents

**Active Ingredients:**

Iodine .....	3.50%
Phosphoric Acid .....	24.70%
<b>Inert Ingredients .....</b>	<b>71.80%</b>
<b>Total .....</b>	<b>100.00%</b>

**KEEP OUT OF REACH OF CHILDREN DANGER**

See side panels for additional precautionary statements and Statement of Practical Treatment.

EPA Reg #69251-1  
EPA Est #14167-IN-01  
Net Contents 136 ml (4.6 oz)

**DIRECTIONS FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Read entire label and Material Safety Data Sheet before using product.

Use one part Viodine per every eleven parts of whole blood contained in the suction canister. Viodine should be introduced 30 minutes before disposal. After procedure, empty Viodine into canister in one of two ways: 1) With the suction running, insert the suction wand or patient tubing into the Viodine bottle and drain; or, 2) open the accessory port and pour the contents of the Viodine bottle completely into the canister. Re-cap canister and mix thoroughly. Transport canister to hopper and dispose of contents according to institutional procedure. Dispose of empty canister according to institutional procedure.

**STORAGE AND DISPOSAL:** Do not store below 0°F. Do not contaminate water, food or feed by storage or disposal. Pesticide wastes are hazardous. Improper disposal of excess pesticides or rinsate is a violation of Federal Law. **Bottle Disposal:** Triple rinse.

**SANITIZER:** Viodine has been tested and is effective as a liquid medical waste sanitizer.

**MANUFACTURED FOR:**

ViAtro, Corp.  
Cleveland, Ohio U.S.A.  
(800) 337-8243

Catalog Reorder Number 2015



2015-11

ACCEPTED with COMMENTS in EPA Letter Dated:

JAN 11 1999

Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended by the FIFRA Act of 1972

69251-1

383



**PRECAUTIONARY STATEMENTS:** Hazards to Humans and Domestic Animals. **CORROSIVE:** Causes irreversible eye damage and skin irritation. Harmful if swallowed. Do not get in eyes, on skin or on clothing. Avoid contamination of food. Avoid breathing fumes or vapors. Wash hands thoroughly with soap and water after handling. Remove and wash contaminated clothing before reuse. **Environmental Hazards:** This pesticide is toxic to fish and aquatic organisms.

**STATEMENT OF PRACTICAL TREATMENT: IF SWALLOWED:** Call a doctor or get medical attention. Do not induce vomiting or give anything by mouth to an unconscious person. Drink promptly, a large quantity of milk, eggwhites, gelatin solution, or if these are not available, drink large quantities of water. Avoid alcohol. **IF IN EYES:** Hold eyelids open and flush with a steady, gentle stream of water for 15 minutes. Get medical attention. **IF ON SKIN:** Wash with plenty of soap and water. Get medical attention. **NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed. For medical emergency, contact local Poison Control Center. In case of spill, contact Chemtec at 1-800-424-9300.

**ADDITIONAL PRECAUTIONARY STATEMENTS:** FOR INSTITUTIONAL HEALTHCARE USE ONLY. WEAR HAND AND EYE PROTECTION WHEN HANDLING POTENTIALLY INFECTIOUS WASTE. NOT FOR INTERNAL USE. FOR SINGLE USE ONLY.



**Liquid Medical Waste Sanitizer**

**For The Safe Disposal Of Suction Canister Contents**

**Active Ingredients:**

Iodine .....	3.50%
Phosphoric Acid .....	24.70%
Inert Ingredients .....	71.80%
<b>Total .....</b>	<b>100.00%</b>

**KEEP OUT OF REACH OF CHILDREN DANGER**

See side panels for additional precautionary statements and Statement of Practical Treatment.

EPA Reg #69251-1  
EPA Est #14167-IN-01

Net Contents 136 ml (4.6 oz)

**DIRECTIONS FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Read entire label and Material Safety Data Sheet before using product.

Use one part Viodine per every eleven parts of whole blood contained in the suction canister. Viodine should be introduced 30 minutes before disposal. After procedure, empty Viodine into canister in one of two ways: 1) With the suction running, insert the suction wand or patient tubing into the Viodine bottle and drain; or, 2) open the accessory port and pour the contents of the Viodine bottle completely into the canister. Re-cap canister and mix thoroughly. Transport canister to hopper and dispose of contents according to institutional procedure. Dispose of empty canister according to institutional procedure.

**STORAGE AND DISPOSAL:** Do not store below 0°F. Do not contaminate water, food or feed by storage or disposal. Pesticide wastes are hazardous. Improper disposal of excess pesticides or rinsate is a violation of Federal Law. **Bottle Disposal:** Triple rinse.

**SANITIZER:** Viodine has been tested and is effective as a liquid medical waste sanitizer.

**MANUFACTURED FOR:**

**VIATRO, Corp.**  
Cleveland, Ohio U.S.A.  
(800) 337-8243

Catalog Reorder Number 2015



2015-11