UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

FEB 21 1997

Dr. Jospeh G. Martorano Medtrol, Inc. 7157 Austin Avenue Niles, IL 60714

Dear Dr. Martorano:

Subject: Super Chlor

EPA Registration No. 69687-1

Your Resubmission Dated November 25, 1996

This is in response to your resubmission of revised labeling for the subject product.

The labeling referred to above, submitted connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act, as amended is acceptable, subject to the following comment.


Change "First Aid" to read "Statement of Practical Treatment (First Aid)."

A stamped copy is enclosed for your records.

If you have any questions, please call Marianne Clark at (703) 305-7879.

Sincerely yours,

Robert Brennis Acting Product Manager (32) Antimicrobial Program Branch Registration Division (7505W)



EPA Form 1320-1A (1/90)

Printed on Recycled Paper

☆ U.S. GOVERNMENT PRINTING OFFICE 1995 - 519-457

SUPER-CHLOR™

For wiping environmental surfaces & patient care equipment.

Active Ingredient: Sodium hypochlorite 0.525% *Inert Ingredients*: 99.475% *Inert Ingredients do not include weight of towelette.

SUGGESTED AREAS OF USAGE: Medical, Dental and Laboratory Counters, Exam Tables, Carts, Point-of-Care Equipment, Telephones, Sink Tops, Toilet Seats See panel on reverse for: Directions for Use

KEEP OUT OF REACH OF CHILDREN WARNING

FIRST AID

•If swallowed, drink large amounts of water. Do not induce vomiting.

•If contact with eyes occurs, flush with water for 15 minutes. Get prompt medical attention.

If contact with skin occurs wash with plenty of soap and water.

FOR INSTITUTIONAL AND COMMERCIAL USE ONLY

•Do not reuse empty packet and towelette •Discard used towelette in designated waste container.

EPA Reg. No. 69687-1 EPA Est #069687-1L-001 Net Contents 0.6 fl. oz.

in IP. Lotter Dated

Under the Federal Inserticide,
Franchischer St. 100 Federal Reg. No. 69687-1
registere: Coder EZA Reg. No.

Directions for use:

- •It is a violation of Federal Law to use this product in a manner inconsistent with its labeling
- •Always use personal protective equipments
- •Open Super-Chlor Packet
- •Remove premoistened 8"x10" towelette
- ·Apply towelette and wipe desired surface to be decontaminated
- ·Allow to air day and discard (see front panel)

This product is not to be used as a terminal sterifant/high level distribution or any surface instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarilly penetrate the blood barriers or otherwise enter normally sterile areas of the body. This product may be used to precise not decontaminate critical or semi-critical medical devices prior to sterilization or high level distribution.

Super-Chlor is a 1:10 dilution of 5.25% sodium hypochlorite solution with distilled water to comply with recommended infection control protocol.

Complete product line information available by calling:

MEDIROL

7157 Austin Avenue, Niles, IL. 60714 PH: 847-647-6555 FAX: 847-647-6611

.001

ticides

69687-1