

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, D.C. 20460



United States
Environmental Protection
Agency

Office of Pesticide Programs

Garrison Dental Solutions
110 DeWitt Lane
Spring Lake, MI 49456

JUN 16 2008

Attention: Clarke Reberg, Director
Product Development

Subject: H2Pro Maintenance Treatment
EPA Registration No. 75829-1
Amendment Dated April 22, 2008

The amendment, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act (FIFRA), as amended, to request a product name change, is acceptable.

OLD PRODUCT NAME

NEW PRODUCT NAME

"H2Pro Bottle A Maintenance
Treatment"

"H2Pro Maintenance
Treatment"

A stamped copy of the "accepted" product labeling is enclosed for your records.

If you have any questions concerning this letter, please contact Martha Terry at (703) 308-6217.

Sincerely

Marshall Swindell
Product Manager 33
Regulatory Management Branch 1
Antimicrobials Division (7510P)

Enclosure

2 8 2

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

H₂Pro

PESTICIDE STORAGE

Store at room temperature out of direct sunlight. Do not contaminate water, food or feed by storage or disposal.

PESTICIDE DISPOSAL

Pesticide wastes are toxic. Improper disposal of excess pesticide is a violation of federal law. If these wastes cannot be disposed of by use according to label instructions, contact your state pesticide or environmental control agency, or hazardous waste representative at the nearest EPA regional office for guidance.

CONTAINER DISPOSAL

Securely wrap original container in several layers of newspaper and discard in trash.

Exp Date: 042209 Lot # 11308A



To reorder call 888.437.0032 or visit us at www.garrisondental.com

Maintenance Treatment

Dental Water Line Treatment

H₂Pro has been specially formulated and clinically proven to control odor causing bacteria and clean dental unit water lines.

Active Ingredient:	
Silver	0.0015%
Inert Ingredients:	99.9985%
Total	100.0000%

This product contains 0.0001245 lbs of silver per gallon.

CAUTION

Keep Out of Reach of Children

Manufactured by
Garrison Dental Solutions
Spring Lake, MI 49456

EPA Registration No. 75829-1
EPA Establishment No. 075829-MI-001
Net Contents: 4 US Ounces (118 ml)

DIRECTIONS FOR USE

For additional instructions see enclosed Clean Start Guide.

FOR PROFESSIONAL USE ONLY

A general-purpose cleaner for use in self contained dental unit water line systems.

MAINTENANCE TREATMENT

After initial dental unit water line treatment using H₂Pro Clean Start Bottles A plus B, pour entire contents of H₂Pro Maintenance Treatment (alternate formula of EPA Reg. No. 75829-1 Bottle A) into the 1 gallon dispensing container provided. Add 1 US gallon (3.785 liters) of distilled or deionized water (32:1 dilution). Dispense from container into self contained dental unit water bottle. Run all lines attached to the dental chair until blue coloration appears exiting the lines. Use continuously in dental water lines instead of distilled, deionized or tap water.

U.S. Patent No. 6,544,427
Additional Patents Pending

**ACCEPTED
with COMMENTS
EPA Letter Dated:**

JUN 16 2008

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No.

75829-1