

88804-1

2/27/2013

1/3



U.S. ENVIRONMENTAL PROTECTION AGENCY

Office of Pesticide Programs
Antimicrobials Division (7510C)
1200 Pennsylvania Avenue NW
Washington, D.C. 20460

EPA Reg. Number:
88804-1

Date of Issuance:
FEB 27 2013

Term of Issuance:
Conditional

Name of Pesticide Product:
Aqua-Tec

NOTICE OF PESTICIDE:

Registration
 Reregistration

(under FIFRA, as amended)

Name and Address of Registrant (include ZIP Code):

Swizz Style, Inc.
465 W. 3rd Street
Dover, OH 44622

Note: Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby registered/reregistered under the Federal Insecticide, Fungicide and Rodenticide Act.

Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA sec 3(c)(7)(A) provided that you:

1. Submit and/or cite all data required for registration of your product under FIFRA sec. 3(c)(5) when the Agency requires all registrants of similar products to submit such data; and submit acceptable responses required for re-registration of your product under FIFRA section 4.

2. Make the labeling changes listed below before you release the product for shipment:

a. Revise the "EPA Registration Number to read, "EPA Reg. No. 88804-1".

Signature of Approving Official:

Marshall Swindell
Product Manager Team-33
Regulatory Management Branch I
Antimicrobials Division (7510P)

Date:

FEB 27 2013

Page 2
EPA Registration No. 88804-1


b. Please note that the studies for the storage stability and corrosion characteristics must be submitted to the Agency within one (1) year from the date of this Notice of Registration.

4. Submit two (2) copies of your final printed labeling before distributing or selling the product bearing the revised labeling.

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the "accepted" label is enclosed for your records.

Sincerely,


Marshall Swindell
Product Manager 33
Regulatory Management Branch I
Antimicrobials Division (7510P)

Enclosure: (Stamped Label)

Front Panel

Protects your lab samples
With confidence

Swizz Style, Inc.

Aqua-Tec™

Antimicrobial Water Preservation Cell for Use in CO2 Incubators and Water Baths
To Prevent Mold, Mildew, Algae, and Non-Public Health Bacteria
Laboratory Use Only

PREVENTS water contamination
~~Salt and heavy metal~~
Protects up to 6 months
Designed for all laboratory water types

**ACCEPTED
with COMMENTS
EPA Letter Dated:
FEB 27 2013**

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended for the pesticide,
registered under EPA Reg. No.
88804-1

**Keep Out of the Reach of Children
CAUTION**

Active Ingredient:
Silver (as elemental) 62%
Other Ingredients..... 38%
Total..... 100%

Precautionary Statements and First Aid
See Instructions on Reverse

Net weight:

Back Panel

It is violation of Federal Law to use this product in a manner inconsistent with its label.

Directions for Use:

1. Place Aqua-Tec in water reservoir (pan or external bottle) of the CO2 incubator or water bath.
2. Use the suction cups to keep the Aqua-Tec in a fixed position.
3. Replace Aqua-Tec every 6 months to maintain effective protection against contamination. Use the peel off sticker as a reminder.

STORAGE: Do not store together with: Strong oxidants, strong acids (nitric acid, sulphuric acid), strong lye, halogens, halogen connections, concentrated hydrogen peroxide solutions.

DISPOSAL: Dispose of by wrapping the exhausted cartridge in newspaper and place in trash can. Ensure compliance with all local and state regulations.

Hazards to Human and Animals

CAUTION. Causes eye irritation. May be harmful if swallowed. Avoid contact with skin, eyes or clothing. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash before reuse.

First Aid

IF IN EYES	Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call poison control center or doctor for treatment advice
IF SWALLOWED	Immediately call a poison control center or doctor. Do not induce vomiting unless directed by a poison control center or doctor. Do not give any liquid to the person. Do not give anything by mouth to an unconscious person
IF ON CLOTHING or SKIN	Take off contaminated Clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call poison control center or doctor for treatment advice.

FOR MEDICAL ASSISTANCE IN THE US, CALL 1-800-222-1222

EPA Reg No. 88804-R
EPA Est. No. 88804 - XX - XXX

Swizz Style, Inc.
465 W. 3rd St.
Dover, OH 44622
1-877-663-7895

