

PCB MQLs Table (Unit: pg/l)

Congener	MQLs*	Congener	MQLs*	Congener	MQLs*	Congener	MQLs*
1	25	57	25	124	25	181	25
2	25	58	25	126	25	182/187	25
3	25	61/70	25	127	25	183	25
4/10	50	62	25	128/162	25	184	25
5/8	50	63	25	129	25	185	25
6	50	65	25	130	25	186	25
7/9	50	66/76	25	131	25	188	25
11	50	67	25	132/161	25	189	25
12/13	50	68	25	133/142	25	190	25
14	50	73	25	134/143	25	191	25
15	50	74	25	135	25	192	25
16/32	25	77	25	136	25	193	25
17	25	78	25	137	25	194	25
18	25	79	25	138/163/164	25	195	25
19	25	80	25	139/149	25	196/203	25
20/21/33	25	81	25	140	25	197	25
22	25	82	25	141	25	198	25
23	25	83	25	144	25	199	25
24/27	25	84/92	25	145	25	200	25
25	25	85/116	25	146/165	25	201	25
26	25	86	25	147	25	202	25
28	25	87/117/125	25	148	25	204	25
29	25	88/91	25	150	25	205	25
30	25	89	25	151	25	206	25
31	25	90/101	25	152	25	207	25
34	25	93	25	153	25	208	25
35	25	94	25	154	25	209	25
36	25	95/98/102	25	155	25		
37	25	96	25	156	25		
38	25	97	25	157	25		
39	25	99	25	158/160	25		
40	25	100	25	159	25		
41/64/71/72	25	103	25	166	25		
42/59	25	104	25	167	25		
43/49	25	105	25	168	25		
44	25	106/118	25	169	25		
45	25	107/109	25	170	25		
46	25	108/112	25	171	25		
47	25	110	25	172	25		
48/75	25	111/115	25	173	25		
50	25	113	25	174	25		
51	25	114	25	175	25		
52/69	25	119	25	176	25		
53	25	120	25	177	25		
54	25	121	25	178	25		
55	25	122	25	179	25		
56/60	25	123	25	180	25		

Note * If adjusted Reporting Limits (RL) are used to adjust MQLs due to laboratory's contemporary ambient background, such adjusted RL shall be updated no less than once per six months. If laboratory method blank, field blank or trip blank subtraction are used in calculation of sample analytical result, supporting document shall be submitted with the Semiannual Status Report.