

NDC 0178-0401-64

# MISSION® WHIRLPOOL SUPER-CONCENTRATE

A POVIDONE-IODINE GERMICIDE FOR THE  
DISINFECTION OF  
COMMERCIAL WHIRLPOOL BATHS

Active Ingredient: Polyvinylpyrrolidone-Iodine ..... 20%  
Inert Ingredients: ..... 80%  
(2% Titratable Iodine)

**CAUTION: KEEP OUT OF THE REACH OF CHILDREN**  
See other caution on the back panel.

Net contents: ONE-HALF GALLON

**ACCEPTED**  
43733-1  
OCT 14 1982  
Under Federal Insecticide  
Fungicide, and Rodenticide Act,  
as amended, for the pesticide  
registered under  
EPA Reg. No.

**MISSION®  
WHIRLPOOL  
SUPER-CONCENTRATE**

FOR THE DISINFECTION OF WATER IN  
WHIRLPOOL BATHS

**MISSION® WHIRLPOOL SUPER-CONCENTRATE** is designed for use in stationary and mobile whirlpools. **MISSION® WHIRLPOOL SUPER-CONCENTRATE** is very concentrated and disinfection is accomplished with only 1/2 fluid ounce per 20 gallons of water, providing 4 ppm of titratable iodine.

**MISSION® WHIRLPOOL SUPER-CONCENTRATE** produces minimum foam, will not stain skin or natural fabrics, and does not produce any objectionable odor.

**PRECAUTION:** Tests of thyroid function of some patients may be affected. Individuals with thyroid conditions should use this product under medical supervision.

**DIRECTIONS FOR USE:** Fill the tank with water. Measure and pour approximately 1 fl. oz. of **MISSION® WHIRLPOOL SUPER-CONCENTRATE** per 40 gallons of water. Drain and disinfect tank after each use with a suitable disinfectant.

**CAUTION:** Avoid contact of undiluted **SUPER-CONCENTRATE** with eyes. Should contact occur, wash thoroughly with water.

STORE IN A COOL PLACE

EPA EST. NO. 43733-TX-01  
EPA REG. NO. 43733-1

 **Mission**  
PHARMACEUTICAL COMPANY  
SAN ANTONIO, TEXAS 78209  
Manufacturer of Five Pharmaceuticals

BEST DOCUMENT AVAILABLE