

NOT REVIEWED  
in Accordance with FR Notice 02-2  
Based on Draft Labeling Date

**DIRECTIONS FOR USE:**

**COLD DISINFECTANT**

For disinfection of previously cleaned surgical instruments, barber tools and dental equipment, simply submerge instruments into solution containing 1 1/2 ounce per gallon of water for 10 minutes.

**RESTAURANT AND BAR RINSE**

Dishes, glassware, silverware, cooking utensils, wash with soap or synthetic detergent, rinse thoroughly and immerse in sanitizing solution containing 1 ounce to 4 gallons of water. No Terminal Potable Water Rinse Required.

**FOR ROOM FOGGING**

**DIRECTIONS:**

To minimize the danger of cross contamination from environmental surfaces apply solution as a fog as described below, before applying standard routine terminal cleaning and disinfecting practices.  
1. Remove all human, animal, and plant life from room. 2. Open closet doors and drawers. 3. Set up 34 in. revolving platform in center of room. 4. Mount Jet FOGMASTER TRI-JET MODEL 6208 or other fogging device delivering equivalent spray at 3 RPM. 5. Fill sprayer reservoir with 1.5 oz (o) 1 gallon of water to product 1150 ppm active ingredients. 6. Set sprayer mechanism to deliver 1 gallon of solution. 7. Fog for 15 minutes for an average 2100 cu. ft. room. For different room sizes, vary spray time proportionately to ensure complete wetting of exposed surfaces. 8. Wait 2 hours before entering room after treatment. 9. Rinse fogging equipment thoroughly with clear water following use.

**DANGER**

**KEEP OUT OF REACH OF CHILDREN.** Corrosive. Causes severe eye and skin damage. do not get in eyes, on skin or on clothing. Protect eyes and skin when handling concentrate. Harmful if swallowed. Avoid contamination of food.

**FIRST AID**

In case of contact, immediately flush eyes or skin with Plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink egg whites, gelatin solution, or if these are not available, drink large quantities of water. Call a physician.

**NOTE TO PHYSICIAN:**

Probable mucosal damage may contraindicate the use of gastric lavage.

# FORMULA #200 SANITIZER DISINFECTANT

**ACTIVE INGREDIENTS:**

n-Alkyl (60% C 14, 30% C 16, 5% C 12, 5% C 18) dimethyl benzyl ammonium chlorides.....	5%
n-Alkyl (68% C 12, 32% C 14) dimethyl ethylbenzyl ammonium chlorides.....	5%
<b>INERT INGREDIENTS .....</b>	<b>90%</b>
	<b>100%</b>

This product is a concentrate and should be diluted before using.

**NET CONTENTS:**

EPA Reg. No. 11712-24  
EPA EST NO. 11712-GA-01

## KEEP OUT OF REACH OF CHILDREN DANGER

See side panel for additional precautionary statements.  
Rinse empty container thoroughly with water before discarding  
or returning drum to drum reconditioner.

Manufactured By:  
**KOR-CHEM, INC.**  
5800 Bucknell Drive  
Atlanta, Georgia 30336

**BEST AVAILABLE COPY**

**USE DIRECTIONS:**

**HOUSEHOLD:**  
For sanitizing hard surface, base Pails, shower stalls with 1 ounce to 5 gallon. Salmonella tested areas.

**HOSPITALS:**  
For floors, Homes, Ins. Telephones of water. Ap will be effective. 10. For he cleaning st

**FOOD PROCESSING:**  
For sanitizing equipment of water to p level. No Po

**DAIRIES:**  
To sanitize d milk cans, fit Then apply gallons of w Rinse is req Board.

This hard offer virus