

Cross-contamination is of major housekeeping concern not only in hospitals, but in schools, institutions, and industry. Environ-3 is formulated for this problem area. It both cleans and disinfects effectively and is virucidal when used as directed. Its hard surface disinfecting action will reduce the hazard of cross-contamination. Environ-3 is an extremely effective one-step sanitizer for non-food contact surfaces. It cleans and sanitizes in waters up to 850 ppm hardness. This hard water tolerance is in excess of most potable waters and will insure complete sanitizing action when used as directed.

Two ounces of Environ-3 per gallon of water will deodorize surfaces in bathrooms, garbage storage areas and other areas where bacterial growth can cause malodors. Prevents odor caused by mold and mildew.

Precautionary Statements Hazard to humans and Domestic Animals:

WARNING:

Keep out of reach of children. May cause severe eye irritation or eye damage. Causes skin irritation. Do not get in eyes, on skin, or on clothing. Harmful if swallowed. Avoid contamination of food.

STORAGE AND DISPOSAL

THIS PRODUCT MUST BE KEPT UNDER LOCKED STORAGE SUFFICIENT TO MAKE IT INACCESSIBLE TO CHILDREN OR PERSONS UNFAMILIAR WITH ITS PROPER USE.

- DO NOT CONTAMINATE WATER, FOOD, OR FEED BY STORAGE OR DISPOSAL
- OPEN DUMPING IS PROHIBITED
- DO NOT REUSE EMPTY CONTAINER

PESTICIDE DISPOSAL

PESTICIDE, SPRAY MIXTURE OR RINSATE THAT CANNOT BE USED OR CHEMICALLY REPROCESSED SHOULD BE DISPOSED OF IN A LANDFILL APPROVED FOR PESTICIDES OR BURIED IN A SAFE PLACE AWAY FROM WATER SUPPLIES.

CONTAINER DISPOSAL

TRIPLE RINSE (OR EQUIVALENT) AND DISPOSE IN AN INCINERATOR OR LANDFILL APPROVED FOR PESTICIDE CONTAINERS, OR BURY IN A SAFE PLACE.

GENERAL

CONSULT FEDERAL, STATE OR LOCAL DISPOSAL AUTHORITIES FOR APPROVED ALTERNATIVE PROCEDURES SUCH AS LIMITED OPEN BURNING

HOLCOMB ENVIRON-3

Cleaner-Disinfectant-Deodorizer-Fungicide-Virucide* Sanitizer with Organic Soil Tolerance for Hospital, Institutional and Industrial Use.

Active Ingredients

Didecyl dimethyl ammonium chloride	4.5%
Tetrasodium ethylenediamine tetraacetate	1.9%
Sodium carbonate	1.0%
Sodium metasilicate, anhydrous	0.5%

Inert Ingredients

92.1%
100.0%

ACCEPTED
APR 16 1981
313-56

KEEP OUT OF REACH OF CHILDREN

WARNING:

Statement of Practical Treatment

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of water, egg whites, gelatin solution. Avoid alcohol. Call a physician immediately.

SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS

EPA Registration No. 323-56

EPA Establishment No. 00323-TX-01

HOLCOMB ENVIRON-3

Cleaner-Disinfectant-Deodorizer-Fungicide-Virucide* Sanitizer with Organic Soil Tolerance for Hospital, Institutional and Industrial Use.

Active Ingredients

Didecyl dimethyl ammonium chloride	4.5%
Tetrasodium ethylenediamine tetraacetate	1.9%
Sodium carbonate	1.0%
Sodium metasilicate, anhydrous	0.5%
Inert Ingredients	92.1%
	100.0%

KEEP OUT OF REACH OF CHILDREN

WARNING:

Statement of Practical Treatment

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of water, egg whites, gelatin solution. Avoid alcohol. Call a physician immediately.

SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS

EPA Registration No. 323-56

EPA Establishment No. C0323-TX-01

Directions for Use GENERAL CLASSIFICATION

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Environ-3 is a proven "one step" cleaner-disinfectant-sanitizer-fungicide-virucide in the presence of moderate amounts of organic soil.

Apply Environ-3 to walls, floors and other hard (inanimate) non-porous surfaces, such as tables, chairs, countertops, sinks and bedframes with a cloth, mop or mechanical spray device so as to thoroughly wet surfaces. For heavily soiled areas, a preliminary cleaning is required. Prepare a fresh solution daily or when use solution becomes visibly dirty.

Disinfection - to disinfect hard, non-porous hospital surfaces, add 2 1/2 oz. per gallon of water. Treated surfaces must remain wet for 10 minutes. At this use level, Environ-3 is also virucidal and fungicidal. At 1 1/4 oz. per gallon, Environ-3 will disinfect hard non-porous surfaces in schools, industry and non-medical institutions.

Sanitizing - to sanitize porous or non-porous, non-food contact surfaces, add 1 oz. per 2 1/4 gallons of water. Treated surfaces must remain wet for 30 seconds.

For hospitals and nursing homes, add 2.5 oz. per gallon of water. The broad spectrum activity of Environ-3 has been evaluated in the presence of 5% serum and found to be effective against the following organisms by the AOAC Use Dilution test:

<i>Pseudomonas aeruginosa</i>	<i>Klebsiella pneumoniae</i>
<i>Staphylococcus aureus</i>	<i>Enterobacter aerogenes</i>
<i>Salmonella choleraesuis</i>	<i>Salmonella schottmuelleri</i>
<i>Escherichia coli</i>	<i>Streptococcus faecalis</i>
<i>Streptococcus pyogenes</i>	<i>Shigella dysenteriae</i>
<i>Proteus mirabilis</i>	<i>Brevibacterium ammoniagenes</i>

AOAC Fungicidal Test - Environ-3 is an effective fungicide against *Trichophyton mentagrophytes* (the athlete's foot fungus) when used on surfaces in areas such as locker rooms, dressing rooms, shower and bath areas, exercise facilities, etc.

Virucidal Performance - *At 2 oz. per gallon use level, Environ-3 was evaluated in the presence of 10% serum and found to be effective against the following viruses: Herpes Simplex (a member of the virus family that causes infectious mononucleosis), Vaccinia (representative of the pox viruses), Adenovirus Type 5 (a causative virus in respiratory diseases) and influenza A2 as represented by the strains commonly called the Hong Kong Flu and the London Flu Virus, on inanimate environmental surfaces.

For schools, non-medical institutional and industrial uses - At 1-1/4 oz. per gallon of water, Environ-3 delivers excellent cleaning and germicidal effectiveness. It is effective against *Staphylococcus aureus*, *Salmonella choleraesuis* and *Escherichia coli*. The same AOAC tests used above to confirm suitable performance for hospitals were used.

Sanitizing non food contact surfaces (such as walls, floors, tables, etc.) - At 1 oz. / 2 1/4 gallons use level, Environ-3 is an effective sanitizer in waters of up to 850 ppm hardness against *Staphylococcus aureus* and *Klebsiella pneumoniae* in the presence of 1% serum on hard porous and non-porous environmental surfaces.

19521 (12/80) PRINTED IN U.S.A.

73