

6 UNDER THE FEDERAL FUNGICIDE AND FOR ECONOMIC ED UNDER NO. 334-20

**IONS** This product must al-  
r use. Never use full strength. See  
active with Soft, Medium, Hard,  
r.

**LEANING AND DISINFECTING:** Use  
r best results add the water to the

**ROOMS, NURSERIES, DELIVERY  
' ROOMS, MATERNITY WARDS—**  
spital Disinfectant and Deodorant  
surfaces such as Metal Beds and  
airs, Linen Carts, Hampers, Refuse  
rts, Clean-Up Carts, Mattress Cov-  
Ducts, Toilet Seats, Telephones,  
nals, Light Switches, Office Ma-  
Chutes, Door Knobs, and other  
the surface has been properly  
in reducing the hazard of cross

of operating rooms: Soak mat with

s Covers: Clean with 2% solution  
or 24 hours.

or utensils or instruments that may  
with patients: Use a 2% solution.

s, ward rooms, etc., should be wet  
in dry swept to prevent dust from  
germs. Use a 2% solution.

Use a 1% solution. Wash bathtubs,  
als, sinks, floors, and baseboards  
at least twice a day.

ie bowl, remove as much water as  
ing it over the trap.

**REFUSE PAILS**  
te by sprinkling or spraying with a  
his will aid in repelling and pre-  
ing of flies and check putrefaction  
often as necessary.

**ING AREAS**  
d products and packaging materials  
d from the room or carefully pro-  
2% solution for general cleaning  
After use, equipment and utensils  
be thoroughly rinsed with potable


o. 334-205 AA

MANUFACTURED BY



AGO, ILLINOIS 60609

**lysan**



# Hyte 13

tuberculocidal  
staphylocidal • fungicidal  
salmonellacidal  
pseudomonacidal

# Hospital DISINFECTANT

<b>ACTIVE INGREDIENTS</b> .....	<b>38.71%</b>
Isopropyl Alcohol, Sodium Dodecylbenzenesulfonate, Ortho-Benzyl-Para-Chlorophenol, Orthophenylphenol, Tetrasodium Ethylene Diamine Tetraacetate, Para-Tertiary-Amylphenol	
<b>INERT INGREDIENTS</b> .....	<b>61.29%</b>
Water, Propylene Glycol, Potassium Sulphite, Citrus Perfume.	
	<b>100.00%</b>

**WARNING: KEEP OUT OF REACH OF CHILDREN. See side panel for further precautionary statements**

**Net \_\_\_ Gal**

This **DISINFECTANT** kills strains of **STAPHYLOCOCCUS AUREUS EASILY**. This assures a pleasant, clean, surrounding. It is used in Hospitals, Nurseries, Operating Rooms, Corridors, and on Equipment. In Toilet Rooms, Public Conveyances, Asylums, Motels, Hotels, Institutions, Schools, to destroy many types of pathogenic organisms as well as many odors. Its high potency permits broad latitude in dilution and thus is economical to use.

Equal caution and care should be taken in **SCHOOLS, BUSES, and Other Public Conveyances AS WELL AS FACTORIES, OFFICES, HOMES, PUBLIC BUILDINGS, AUDITORIUMS, and various other INSTITUTIONS** by mopping or spraying with a 1/2% solution of **DISINFECTANT** on all floors, walls, seats, passageways, lockers, cabins, and pantries.

Effective on odors from any source in hospitals, nursing homes, ambulances, from postcolostomy, post-ileostomy dressing and bags, malignant diseases; in the home, kitchen, bathroom or the like, industrial plant processes, food processing and packaging, fishing equipment, food stores, burnt and smoke odors.

DILUTION CHART		
% SOLUTION	TO MAKE 1 GALLON	TO MAKE 1 QUART
1/2% Solution	20 cc or 2/3 fl. oz. Remainder Water	One teaspoonful Remainder Water
1% Solution	40 cc or 1-1/3 fl. oz. Remainder Water	Two teaspoonfuls Remainder Water
2% Solution	80 cc or 2-2/3 fl. oz. Remainder Water	Four teaspoonfuls Remainder Water

**RECOMMENDED USE**  
**USE A 1/2% SOLUTION** for all GENERAL DISINFECTION, FLOORS, WALLS, EQUIPMENT, SICK ROOM RECEPTACLES, RUBBER GOODS, INSTRUMENT STORAGE, and as ANTI-RUST AGENT IN HEAT STERILIZERS.  
**USE A 1% SOLUTION** for fungicidal action against Trichophyton interdigitale.  
**USE A 2% SOLUTION** for TUBERCULOSIS AND PSEUDOMONAS DISINFECTION.  
**USE A 2% SOLUTION** for COLD INSTRUMENT DISINFECTION: Soak for 15 minutes in covered container.

**WARNING:**  
HARMFUL IF SWALLOWED. AVOID CONTACT WITH SKIN. AVOID SPLASHING IN EYES. If in contact with eyes, wash with large amounts of flowing water. If irritation persists, get medical attention. If taken internally, drink plenty of olive oil or other vegetable oils. CALL A PHYSICIAN AT ONCE.  
Avoid contamination of food and feeds. Rinse empty container thoroughly with water and discard it.  
KEEP OUT OF REACH OF CHILDREN

## DIRECTIONS FOR USE

### GENERAL CLEANING AND DISINFECTING

Dilute 1 1/2 fl. oz. per gallon before use.

Start with clean equipment for all applications. When using a string or sponge mop or rag, the preferred procedure is to empty the rinse bucket to flush clean the applicator before redipping into the diluted solution. This technique minimizes contamination of the cleaner-disinfectant and permits more economical use of the product.

Most soils respond more easily to cleaning if the diluted solution is allowed to remain in contact with the surface for a few minutes before working or picking up.


### VINYL, ASPHALT TILE, TERRAZZO, MARBLE, LINOLEUM, CORK, WOOD, RUBBER AND ALL OTHER TYPES OF FLOORING

Use 1 1/2 oz. per gallon of warm water. For best results, the solution should be in contact with the surface for at least 5 minutes. The solution may be applied to the floor with a mop or sprayer.

**COLD INSTRUMENT DISINFECTION:** Soak for 15 minutes in covered container.

**HOSPITAL OPERATING ROOMS:** Use the "flooding technique." Apply the solution to the floor so that the floor is very wet. Allow to set for at least 5 minutes—then vacuum. This will clean and disinfect the floor.

**WALLS, FURNITURE, PLUMBING FIXTURES, BEDS** should be wiped with a sponge or rag dampened with a 1 1/2% solution. Allow to dry without rinsing.



UNIVERSITY OF ILLINOIS  
GENERAL INSECTICIDE  
FUNGICIDE AND RODENT CIDE ACT  
ECONOMIC POISON REGISTER  
EPA ORDER NO. 334-208

3 34-208

000000

# fyt

# HOSPITAL GERMICIDAL DETERGENT

tuberculocidal • staphylocidal  
salmonellacidal • fungicidal  
pseudomonacidal

ACTIVE INGREDIENTS: ..... 38.38%  
Isopropyl Alcohol, Sodium Dodecylbenzene-  
sulfonate, Ortho-benzyl-para-chlorophenol,  
Tetrasodium Ethylene Diamine Tetraacetate,  
Tripotassium Phosphate, Orthophenylphenol,  
Para-tertiary-amyphenol.  
INERT INGREDIENTS: ..... 61.62%  
Water, perfume.

This Hospital Germicidal Detergent combines the broad spectrum disinfectant activity of phenolics with the efficient cleaning ability of biodegradable alkyl benzene sulfonate and sequesterant. Pleasantly and lightly scented for general patient acceptance. Does not contain film-forming soaps or fatty acids and is effective in hot, cold, hard or soft water.

E.P.A. REG. NO. 334-208AA

**WARNING: KEEP OUT OF REACH OF CHILDREN.** See side panel for further precautionary statements.

NET: \_\_\_ GAL

### TERMINAL DISINFECTION AND CLEANING

Use a 1 1/2% solution to wash floors, walls, furniture, interiors of drawers and cupboards, lighting fixtures, air conditioners, ventilating fans, and all other stationary equipment. **DO NOT RINSE.**

### FUNGICIDAL ACTIVITY (Trichophyton

interdigitale) Use 1-1/2 oz. per gallon and follow directions for general cleaning and disinfecting.

**GARBAGE PAILS, CELLARS, REFUSE CANS, WASH TANKS:** Wash well with a 1 1/2% solution for general disinfection and cleaning.

### TUBERCULOSIS AND PSEUDOMONAS

**DISINFECTION** Use 4 fl. oz. per gallon and follow instructions for terminal disinfection and cleaning.

### FOOD PROCESSING AREAS

Before using, food products and packaging materials must be removed from the room or carefully protected. Use 1-1/2 fl. oz. per gallon for general cleaning and disinfection. After use, equipment and utensils in the area must be thoroughly rinsed with potable water.

### WARNING:

**HARMFUL IF SWALLOWED. AVOID CONTACT WITH SKIN. AVOID SPLASHING IN EYES.** If in contact with eyes, wash with large amounts of flowing water. If irritation persists, get medical attention. If taken internally, drink plenty of olive oil or other vegetable oils. **CALL A PHYSICIAN AT ONCE.**

Avoid contamination of food and feeds.

Rinse empty container thoroughly with water and discard it.

KEEP OUT OF REACH OF CHILDREN.

MANUFACTURED BY

**Kysan**  
CORPORATION

CHICAGO, ILLINOIS 60609