

**Disinfects Influenza A₂/Hong Kong Virus
on inanimate environmental surfaces.**

This Product Will Control the
Following Organisms on Environmental Surfaces: Tubercle Bacilli,
Salmonella Typhosa, Shigella Paradysenteriae, Salmonella Schott-
muelleri, Escherichia Coli, Staphylococcus Aureus, Streptococcus
Pyogenes, Mycobacterium Smezmatis, Penicillium Glaucum. Controls
mildew (mold or fungi). Completely inhibits growth on cloth. Mark-
edly inhibits growth on leather.

This Powerful Hospital Disinfectant and Deodorant acts on germ
laden surfaces such as metal beds and springs, wheel chairs, linen
carts, hampers, refuse cans, dressing carts, clean-up carts, air
ducts, toilet seats, telephones, wash basins, urinals, light switches,
office machinery, laundry chutes, door knobs, and other equipment,
after the surface has been properly cleaned. Will reduce the hazard
of cross infection by bacteria.

Non-Selective Odor Suppressant for space odor conditioning and
deodorizing surgical dressings. Molecular aerosol deodorant attacks
malodor at the source. Does not cause anosmia. Removes unpleas-
ant odors of decay, smoke, garlic, onions, food putrefaction, per-
spiration, stale air, tobacco, sick rooms, musty closets, fish.

Suppresses and Eliminates Odors by chemically neutralizing odor-
producing materials. Effective on odors from any source in nursing
homes, hospitals, ambulances, from post-colostomy, post-ileostomy
dressing and bags, malignant diseases. Acts instantly and remains
effective for a long period. Also for industrial plant processes, food
processing and packing, fishing equipment, food stores, smoke.

DIRECTIONS: SHAKE WELL BEFORE USING

FOR DISINFECTING: Hold dispenser upright about 6 to 8
inches from surface to be disinfected and spray 3 seconds
until surface is wet.

FOR DEODORIZING: Hold dispenser upright—point the open-
ing away from face. Spray into air towards the center of room
for about 4 seconds.

FOR MILDEW CONTROL: Clean surface. Spray surface until
wet. Repeat application every week to maintain control.

IMPORTANT: Do not use on polished wood, furniture or
rayon fabrics.

ACTIVE INGREDIENTS:	57.3356%
4-Chloro-2-Cyclopentylphenol	0.0256%
Ortho-phenylphenol	0.1400%
3, 4', 5-Tribromosalicylanilide	0.3600%
and other related Bromosalicylanilides including:	
3, 5-Dibromosalicylanilide	
4', 5-Dibromosalicylanilide	
2', 3, 4', 5-Tetrabromosalicylanilide	
Essential Oils	2.0000%
Alcohol	54.8100%
INERT INGREDIENTS	42.6644%

WARNING—Do not get in eyes. In case of contact with eyes,
flush at once with water for 15 minutes. Call a physician. Avoid

**BACTERICIDAL • PSEUDOMONICIDAL
STAPHYLOCIDAL • TUBERCULOCIDAL**

RESIST

DISINFECTANT

DEODORANT

**LEMON
SCENTED**

041173

334-296

ACCEPTED

June 1, 1973

UNDER THE FEDERAL INSECTICIDE
FUNGICIDE AND RODENTICIDE ACT
FOR ECONOMIC POISON REGISTER-
ED UNDER NO.
TO COMMENTS.