

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Disinfectant Pump Spray is designed specifically for use in homes, nursing homes, schools and hotels as a disinfectant and deodorant. It is intended for use as a disinfectant for telephones, bed springs, door knobs, linen carts, hampers, wheel chairs, tables, cabinets, light switches, garbage pails, and other clean, inanimate, non-porous hard surfaces. In addition, Disinfectant Pump Spray deodorizes those areas which generally are hard to keep fresh smelling, such as garbage storage areas, empty garbage bins and cans, basements, restrooms, shower stalls, and other areas which are prone to odors caused by microorganisms. Disinfectant Pump Spray can also be used as a room air deodorizer and refreshener.

SURFACE DISINFECTION AND DEODORIZATION - For disinfection, first remove soil deposits from surface to be treated. Hold container six inches from surface and spray until completely wet, allowing the surface to remain wet for 10 minutes.

AIR DEODORIZER AND AIR REFRESHENER - Spray upwards to center of room for three (3) seconds.

MILDEWSTAT - To control mold and mildew on pre-cleaned, hard, non-porous surfaces spray surface to be treated making sure to wet completely. Let air dry. Repeat application at weekly intervals or when mildew growth appears.

BACTERICIDAL ACTIVITY - When used as directed, Disinfectant Pump Spray exhibits effective disinfectant activity against the organisms: *Staphylococcus aureus*, *Salmonella choleraesuis* and *Pseudomonas aeruginosa* and meets all requirements for hospital use.

FUNGICIDAL ACTIVITY - This product is fungicidal against the pathogenic fungi, *Trichophyton mentagrophytes* (Athlete's Foot Fungus) when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, or other clean, non-porous, hard surfaces commonly contacted by bare feet.

KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS.)

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids use disposable latex gloves, gowns, masks, or eye covering.

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of Disinfectant Pump Spray.

CONTACT TIME: Allow surface to remain wet for 10 minutes.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

Stepan

Stepan Company
Northfield, IL 60093

DISINFECTANT PUMP SPRAY

DISINFECTANT AIR REFRESHENER	DEODORIZER
ACTIVE INGREDIENTS	
n-Alkyl (60% C ₁₄ , 30% C ₁₆ , 5% C ₁₂ , 5% C ₁₈) dimethyl benzyl ammonium chlorides	
	0.14%
n-Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzyl ammonium chlorides	
	0.14%
Isopropanol	75.71%
INERT INGREDIENTS	24.01%

Total	100.0%

EPA REG. NO. 1839-82
EPA EST. NO. 1839-IL-01

KEEP OUT OF REACH OF CHILDREN

WARNING

Do not use or store near heat or open flame.

See side panels for additional precautionary statements

NET CONTENTS

PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS

WARNING

KEEP OUT OF REACH OF CHILDREN. Causes eye and skin irritation. Do not get in eyes, on skin or on clothing. Harmful if swallowed. Avoid contamination of food.

STATEMENT OF PRACTICAL TREATMENT

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse.

If swallowed, drink egg whites or gelatin solution, or if these are not available, drink large quantities of water. Call a physician immediately.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage.

STORAGE AND DISPOSAL

DO NOT CONTAMINATE WATER, FOOD, OR FEED BY STORAGE OR DISPOSAL.

CONTAINER DISPOSAL - Do not reuse container (bottle, can, or jar). Rinse thoroughly before discarding in trash.

STORAGE - Store in a dry place no lower in temperature than 50°F or higher than 120°F.

