

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

CD 3.2 detergent/disinfectant has been designed specifically for hospitals, nursing homes, schools, food processing plants, and other institutions where housekeeping is of prime importance.

DISINFECTION — To disinfect hard inanimate surfaces (such as walls, floors, table tops), apply solution with mop, cloth, sponge or mechanical sprayer so as to wet all surfaces thoroughly. Allow to remain wet for 10 minutes and then let air dry.

DILUTION RATES: General Disinfection-add 2 ounces CD 3.2 per gallon of water. Hospital Disinfection-add 3 ounces CD 3.2 per gallons of water.

To disinfect toilet bowls, flush toilet, add 2 ounces CD 3.2 detergent/disinfectant directly to the bowl water. Swab the bowl completely using a scrub brush or toilet mop, making sure to get under the rim. Let stand for 10 minutes and flush.

For heavily soiled areas, a pre-cleaning step is required. Prepare a fresh solution for each use.

MILDEWSTAT: To control mold and mildew on hard non porous surfaces (such as floors, walls, table tops) clean surfaces then add 2 ounces CD 3.2 per gallon of water. Apply with a cloth, mop or sponge to surfaces to be treated making sure to wet thoroughly. Let air dry. Prepare a fresh solution for each use. Repeat application at weekly intervals or when mildew growth appears.

ACCEPTED

NOV 16 1988

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 1839-103

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER

KEEP OUT OF REACH OF CHILDREN. Causes severe eye and skin damage. Do not get in eyes, on skin or on clothing. Wear goggles or face shield and rubber gloves when handling. Harmful or fatal if swallowed. Avoid contamination of food.

STATEMENT OF PRACTICAL TREATMENT

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse.

If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol. Call a physician immediately.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock as well as oxygen and measures to support breathing manually or mechanically may be needed.



CD 3.2 (D & F) Detergent/Disinfectant

CLEANER
DETERGENT
DEODORIZER

DISINFECTANT
FUNGICIDE (against
Pathogenic fungi)

MILDEWSTAT (on hard inanimate
surfaces)

ACTIVE INGREDIENTS

n Alkyl (30% C14, 30% C16, 5% C12,
5% C18) dimethyl benzyl ammonium
chlorides

1.60%

n Alkyl (68% C12, 32% C14)
dimethyl ethylbenzyl ammonium
chlorides

1.60%

INERT INGREDIENTS

96.80%

Total 100.00%

EPA REG. NO. 1839-103

EPA EST. NO. 1839/L-01

**KEEP OUT OF REACH OF CHILDREN
DANGER**

see side panels for additional precautionary statements

Stepan S

Stepan Company
Northfield, IL 60093

NET CONTENTS

CD 3.2 is a cleaner/disinfectant designed for general use in schools, food processing plants and light industrial and hospital applications where housekeeping is of prime importance in controlling the hazard of cross contamination. When used as directed, CD 3.2 is formulated to disinfect inanimate hard surfaces (such as walls, floors, sink tops, toilet bowls, tables, chairs and telephones). For larger areas such as patient care facilities and restrooms, the product is designed for both general cleaning and disinfecting.

In addition, CD 3.2 deodorizes those areas which generally are hard to keep fresh smelling, such as garbage storage areas, empty garbage bins and cans, toilet bowls and any other areas which are prone to odors caused by microorganisms.

DISINFECTION — At 2 ounces per gallon dilution, CD 3.2 exhibits effective disinfectant activity against the organisms *Salmonella choleraesuis* and *Staphylococcus aureus*.

At this level the product is also fungicidal against the Pathogenic fungi *Trichophyton mentagrophytes*.

At 3 ounces per gallon of dilution, CD 3.2 is an effective hospital disinfectant and demonstrates disinfectant activity against *Pseudomonas aeruginosa*.

EFFICACY TESTS HAVE DEMONSTRATED THAT CD 3.2 IS AN EFFECTIVE BACTERICIDE, AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOILING, 10% BLOOD SERUM.

MOLD AND MILDEWSTAT — CD 3.2 controls mold and mildew and the odors they cause on hard non porous inanimate surfaces at 2 ounces per gallon dilution.

STORAGE AND DISPOSAL

DO NOT CONTAMINATE WATER, FOOD, OR FEED BY STORAGE AND DISPOSAL.

STORAGE: Store in a dry place no lower in temperature than 50°F or higher than 120°F.

CONTAINER DISPOSAL — Do not reuse empty container. Triple rinse empty container with water. Return metal drum then offer for reconditioning or puncture and dispose of in a sanitary landfill or by other procedures approved by State and local authorities. Plastic containers may be disposed of in a sanitary landfill, incinerated or, if allowed by local authorities, by burning. If burned, stay out of smoke.

PESTICIDE DISPOSAL — Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.