

Head strong
 concentrated acid
 disinfectant toilet bowl cleaner
 designed for institutional
 and hospital use
 for
 Hospitals-Industry-Homes
 Disinfects-Cleans-Deodorizes
 Removes Stains
 Kills Many Germs in Toilet Bowls
 Including Those Found in Hospitals
 (eg, Pseudomonas aeruginosa & Staphylococcus aureus)

ACTIVE INGREDIENTS:

Octyl decyl dimethyl ammonium chloride.....	1.250%
Diocetyl dimethyl ammonium chloride.....	0.625%
Didecyl dimethyl ammonium chloride.....	0.625%
Hydrogen Chloride.....	17.500%
Essential Oils.....	0.200%
INERT INGREDIENTS.....	79.800%
	<u>100.000%</u>

Contains no phosphate compounds



POISON

DANGER

KEEP OUT OF REACH OF CHILDREN

See back panel for first aid statement and additional precautions.

E.P.A. Reg. No. 3525-61

E.P.A. Est. no. 3525-Nj-1..

Manufactured for:
 UTIKEM PRODUCTS
 Div. of Qualco, Inc.
 225 Passaic Street
 Passaic, NJ 07055

HEAD-STRONG is a concentrated acid disinfectant bathroom toilet bowl cleaner designed for institutional and hospital use. This disinfectant toilet bowl cleaner is especially formulated to kill many germs commonly found in hospital, industrial and household environments.

HEAD-STRONG is designed for use in toilet bowls and urinals.

HEAD-STRONG contains a high level of acid together with an acid compatible detergent system. This unusual combination facilitates the easy removal of resistant stains and permits its disinfectant action.

When used as directed, HEAD-STRONG cleans, disinfects and deodorizes in one easy step.

USE-DIRECTIONS

Flush toilet. Add 1 fluid oz. (2 tablespoons) directly to water in toilet bowl (3/4 gallon). Brush to cover all interior surfaces of the bowl including under the rim. Allow disinfectant solution to remain at least 10 minutes. Flush toilet.

For badly stained bowls, repeat above treatment using 2 to 4 fluid ounces of HEAD-STRONG. For use in urinals, add 1 fluid oz. (2 tablespoons) to 3/4 gallon of water. Apply by brush to all urinal surfaces. Allow disinfectant solution to remain at least 10 minutes. Flush.

At 1 fluid oz. per 3/4 gallon of water, HEAD-STRONG is effective against a broad spectrum of microorganisms commonly found in toilet bowls and urinals. Effectiveness against Pseudomonas aeruginosa, Salmonella choleraesuis, Staphylococcus aureus, Shigella dysenteriae, Streptococcus faecalis and Escherichia coli has been confirmed by the AOAC Use-Dilution procedure.

Reviewed. Registration to be in accordance with
 notification of... for Phase 2
 Status.

Special Note:

The cationic detergent system of HEAD-STRONG inhibits the acid corrosion of metal parts found around toilet bowls and urinals.

DANGER

Corrosive. Keep Out of Reach of Children. May be fatal or harmful if swallowed. Do not breathe vapor or fumes. Produces chemical burns. Do not get in eyes, on skin or on clothing.

FIRST AID

If swallowed, call a physician immediately. Drink a teaspoon or more of milk of magnesia or a cup of soapy water followed with milk or raw egg whites.

Eye Contact: Wash with water for 15 minutes. Get prompt medical attention.

Skin Contact: Wipe off the acid gently and immediately flood the surface with water, using soap freely, then cover with moist magnesia or baking soda.

Rinse empty container with water before discarding it.

