

HOSPITAL SAN-O-SIX

Cleaner • Disinfectant • Deodorizer • Fungicide

Designed for Hospital and Institutional Use

ACTIVE INGREDIENTS: 7.52%
 Alkyl (C14, C16, C18) dimethylbenzyl ammonium chlorides, 2.26% Alkyl (C12, C14, C16, C18) dimethyl ethylbenzyl ammonium chlorides, 2.26% Sodium carbonate, 3.00%
INERT INGREDIENTS: 92.48%

WARNING

KEEP OUT OF REACH OF CHILDREN

Causes eye and skin irritation. Do not get in eyes, on skin or in clothing. Harmful if swallowed. Avoid contamination of food.
FIRST AID: In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, do not induce vomiting. Drink plenty of fluid and call a physician immediately.

EPA Reg. No. 3635-119

NET CONTENTS _____ GALLONS

ACCEPTED
 5-12-72
 UNDER THE FEDERAL REGULATION OF THE ENVIRONMENT
 FUNGICIDE AND DISINFECTANT
 FOR RESIDENTIAL USE
 3635-119



DESCRIPTION

HOSPITAL SAN O SIX is a cleaner, disinfectant, deodorizer and fungicide designed for institutional use in FOOD PROCESSING PLANTS, HOSPITALS, NURSING HOMES, and SCHOOLS where housekeeping will be done on the interior hard surface disinfection using HOSPITAL SAN O SIX is a product that is safe when you need it in water. It is used and tested on a variety of surfaces to aluminum, tin, plastic, rubber, asphalt, and on any surface on which water can be used.

ADVAANTAGES - DILUTION INFORMATION - TESTS

Net weight per gallon	11.26
Net weight per quart	2.82
Risk of eye irritation (PPE)	1:100
.....
.....

DIRECTIONS

.....

 Mop or spray

- B For General Disinfection. To disinfect athlete's foot fungus on surfaces.
 1. Prepare a solution of 1 oz. HOSPITAL SAN O SIX per gallon of water. Mop or spray this onto floors and other surfaces normally contacted by the feet of your team, showers, and around swimming pools.
- C For Sanitizing and Disinfecting Food Processing Equipment.
 1. Clean and sanitize a Equipment as recommended in the manual. Mop or spray a solution of 1 oz. HOSPITAL SAN O SIX per gallon of water on the surface for 15 minutes contact time. Rinse with potable water.
 2. For general disinfection. Follow same procedure but use a solution of 2 oz. HOSPITAL SAN O SIX per gallon of water. Use Agitation where possible. Rinse with potable water.
 3. For equipment or other food contact items use 1/2 oz. HOSPITAL SAN O SIX per gallon of water. Mop or spray solution on equipment. One and one-half hours contact time. Rinse with water.

D72

RINSE EMPTY CONTAINER THOROUGHLY AND DISCARD

MANUFACTURED BY

