

WARNING: KEEP OUT OF REACH OF CHILDREN

See back panel for additional precautions

to disinfect surfaces:

1. Remove gross filth and heavy soil deposits.
2. Hold can upright 10" to 12" from surface.
3. Spray for 3-5 seconds until surface is covered with mist.
4. If desired, wipe dry with cloth.

to deodorize and refreshen air:

Spray as required.

to control mold and mildew on hard surfaces:

1. Clean surface.
2. Spray for 3-5 seconds until surface is covered with mist.
3. Repeat application every 7 days.

WARNING: FLAMMABLE

DO NOT USE NEAR FIRE OR FLAME.

Harmful if swallowed. Don't spray in face. In case of contact, flush eyes immediately with plenty of water. Avoid breathing spray mist. Keep out of reach of children.

Avoid contact with food and food utensils. If sprayed on floor, wipe dry.

Contents under pressure. Do not set on stove or radiator or keep where temperature will exceed 120°F, as container may burst. Do not puncture or throw in fire.

ACTIVE INGREDIENTS:

Ethyl alcohol 24.238%, Isopropyl alcohol 24.237%, n-Alkyl (60% C₁₂, 30% C₁₄, 5% C₁₆, 5% C₁₈) dimethyl benzyl ammonium chloride 0.065%, n-Alkyl (50% C₁₂, 30% C₁₄, 17% C₁₆, 3% C₁₈) dimethyl ethylbenzyl ammonium chloride 0.065%.

INERT INGREDIENTS:

51.395%

EPA Reg. No. 4822-120 AA

BY **Johnson WAX**
ANOTHER **END BAC** PRODUCT
FOR HOSPITAL DISINFECTION

Johnson
END
BAC
II
PRESSURIZED

disinfectant
spray

**kills staph,
strep, many
other surface
germs**

tuberculocidal

**bacteriostatic
residual action
when not rinsed
or washed off**

**eliminates
odors**

■ Kills staph, strep, T. interdigitale (athlete's foot fungus), many other dangerous surface germs

■ Destroys and prevents mold and mildew.

■ Destroys many odor-causing bacteria at their source

Johnson END-BAC II Pressurized disinfectant spray is specially formulated with modern germicide ingredients used in hospitals. Spray regularly around urinals, toilets, shower stalls, other rest room and bathroom fixtures, garbage cans, lockers, in hospital rooms. Disinfects toilet seats, door knobs, and other surfaces against staph, strep and many other dangerous surface germs.

Use to sanitize golf shoes and to disinfect athletic equipment. Refreshens air by eliminating cooking, smoke and rest room odors. This formulation has passed the official A.O.A.C. test against *Staphylococcus aureus*, *Salmonella choleraesuis*, *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis* and *Trichophyton interdigitale*.

Virucidal: Disinfects Influenza A₂/Hong Kong, Parainfluenza Type 1, Adenovirus Type 2 and Herpes simplex viruses carried on inanimate environmental surfaces.

Johnson, KIT-IN-A-CAN and END-BAC are trademarks.

© S. C. JOHNSON & SON, INC.
RACINE, WIS. 53403, U.S.A.



NET WT 17 OZ (1 LB 1 OZ)

A PRODUCT OF **Johnson WAX** SERVICE PRODUCTS DIVISION

42867