

TECH QUATERNARY

a disinfectant for Hospitals, Farms,
Schools and Food Industry
Deodorizes • Sanitizes • Disinfects

Dilute before using—Sanitizing solution—1 oz. TECH QUATERNARY to 4 gallons water will produce 200 ppm active ingredients.

Disinfectant solution—1 oz. TECH QUATERNARY to 2 gallons water will produce 400 ppm active ingredients.

CAUTION

KEEP OUT OF REACH OF CHILDREN

See side panel for other cautions.

U.S.D.A. Reg. No. 5664-22

Net Contents

TECH QUATERNARY

Prior to using TECH QUATERNARY as a disinfectant, thoroughly clean all surfaces with a suitable detergent. Do not mix with soap or ammonia detergent.

Directions

To sanitize hard surfaces, eating and drinking utensils, food processing areas and equipment, use one oz. TECH QUATERNARY in four (4) gallons water (200 ppm active ingredients) and follow with a rinse of clean potable water.

To disinfect and deodorize hard surfaces, such as metal beds, bedpans, food lockers, garbage cans, lockers and lavatories, use one (1) oz. TECH QUATERNARY in two (2) gallons water (400 ppm active ingredients) and follow with a clean potable water rinse.

For porous surfaces such as floors and counter tops, use for sanitizing, one (1) oz. TECH QUATERNARY to two (2) gallons water. For disinfecting, use one (1) oz. TECH QUATERNARY to one (1) gallon water. Allow solution to remain in contact with surface for 10 minutes before rinsing with clean potable water.

CAUTION

Harmful if swallowed. Avoid contact with skin and eyes. In case of contact with eyes, flush with plenty of water for at least 15 minutes. If irritation persists, get medical attention. Avoid contamination of food.

Active Ingredients

N-Alkyl (30% C₁₂, 7% C₁₄, 3% C₁₆)

Dimethyl Ethylbenzyl Ammonium Chloride 5%


N-Alkyl (60% C₁₂, 30% C₁₄, 5% C₁₆, 5% C₁₈)

Dimethyl Benzyl Ammonium Chloride 5%

Inert Ingredients

90%

100%



Manufactured exclusively by

TECHNICAL MAINTENANCE PRODUCTS, INC.

2211 N. American St., Phila., Pa. 19133

2011 Teall Ave., Syracuse, N. Y. 13206