

KANCEL: An effective new sanitizer for space air sanitation and positive deodorization. Specifically designed for use in hospitals, restaurants, industrial plants, homes and areas where mal-odors are an annoyance and reduction of air-borne bacteria, fungi and viruses are desirable.

FOR USE WITH SYNCRO-MIST ACTUATOR: KANCEL (odorless) air sanitizer is designed especially to be used with the Syncro-Mist automatic and periodic electrical actuator. The Syncro-Mist actuator can be set to activate the metered valve on a periodic schedule to treat large or small areas. Usually the effectiveness of KANCEL will be noticeable within one hour after the machine is activated, and if less than 24 hours control is required, the machine should be set to start one hour before control is required and run one hour after control is no longer needed.

Keep in mind that such conditions as air change rate, shape of space treated, humidity and temperature all affect results, therefore follow instructions supplied with each Syncro-Mist unit to achieve maximum results. As with any method of control, improved performance is achieved when good sanitary practices and other supplementary methods are used.

The 12 oz. can of KANCEL contains 3400 shots and will operate for a total of 850 hours when actuated at rate of four times per hour on a 24 hour basis.

Caution: Non-toxic when used as directed. Do not store above 120°F. Do not puncture or incinerate container.

CONTENTS: Net Weight 12 ounces

A Product of
**Syncro-Mist
Controls inc.**
New York 16, N.Y.

SYNCRO-MIST



(ODORLESS)

Kancel

**REDUCES RISK OF
AIR-BORNE INFECTION**

Deodorizer

Sanitizer

Reduces Air-Borne Bacteria

Active Ingre

n-Alkyl (50

17% C₁₆,

benzyl an

chlorides

Triethylene

Propylene

Isopropyl a

Inert Ingredi

Pa
26
17
Of
Pe

Sanitizer for space air sanitation and
 specially designed for use in hospitals,
 homes and areas where mal-odors
 and air-borne bacteria, fungi and

ACTUATOR: KANCEL (odorless) air
 is to be used with the Syncro-Mist
 mechanical actuator. The Syncro-Mist actu-
 ator has the metered valve on a periodic
 basis in all areas. Usually the effectiveness
 is within one hour after the machine
 is started. If 24 hours control is required, the
 actuator should be started one hour before control is re-
 quired. After control is no longer needed.
 Conditions as air change rate, shape of
 room, temperature all affect results, there-
 fore, included with each Syncro-Mist unit to
 be used with any method of control, im-
 proved when good sanitary practices
 are used.
 Contains 3400 shots and will operate
 when actuated at rate of four times per

as directed. Do not store above
 normal temperature container.
 Weight 12 ounces

Product of
Syncro-Mist
Products Inc.
 New York 16, N.Y.

SYNCRU-MIST



(ODORLESS)
Kancel
REDUCES RISK OF
AIR-BORNE INFECTION

Deodorizer

Sanitizer

Reduces Air-Borne Bacteria

Active Ingredients:

- n-Alkyl (50% C₁₂, 30% C₁₄,
 17% C₁₆, 3% C₁₈) dimethyl
 benzyl ammonium
 chlorides 1.8%
- Triethylene glycol 5.0%
- Propylene glycol 5.0%
- Isopropyl alcohol 12.0%

Inert Ingredients: 76.2%

Total 100.00%

Patent Numbers
 2613108
 177170 Des.
 Other Applications
 Pending

ACCEPTED WITH COMMENTS

KANCEL: An effective new sanitizer for space air sanitation and positive deodorization. Specifically designed for use in hospitals, restaurants, industrial plants, homes and areas where mal-odors are an annoyance and reduction of air-borne bacteria, fungi and viruses are desirable.

FOR USE WITH SYNCRO-MIST ACTUATOR: KANCEL (odorless) air sanitizer is designed especially to be used with the Syncro-Mist automatic and periodic electrical actuator. The Syncro-Mist actuator can be set to activate the metered valve on a periodic schedule to treat large or small areas. Usually the effectiveness of KANCEL will be noticeable within one hour after the machine is activated, and if less than 24 hours control is required, the machine should be set to start one hour before control is required and run one hour after control is no longer needed.

Keep in mind that such conditions as air change rate, shape of space treated, humidity and temperature all affect results, therefore follow instructions supplied with each Syncro-Mist unit to achieve maximum results. As with any method of control, improved performance is achieved when good sanitary practices and other supplementary methods are used.

The 12 oz. can of KANCEL contains 3400 shots and will operate for a total of 850 hours when actuated at rate of four times per hour on a 24 hour basis.

Caution: Non-toxic when used as directed. Do not store above 120°F. Do not puncture or incinerate container.

CONTENTS: Net Weight 12 ounces


A Product of
**Syncro-Mist
Controls inc.**
New York 16, N.Y.

SYNCRO-MIST

(ODORLESS)

Kancel

**REDUCES RISK OF
AIR-BORNE INFECTION**



Active Ingredients:

n-Alkyl (50% C₁₂, 30% C₁₄,
17% C₁₆, 3% C₁₈) dimethyl
benzyl ammonium
chlorides 1.8%
Triethylene glycol 5.0%
Propylene glycol 5.0%
Isopropyl alcohol 12.0%

Inert Ingredients: 76.2%

Total 100.00%

Deodorizer

Sanitizer

Reduces Air-Borne Bacteria

Patent Numbers
2613108
177170 Des.
Other Applications
Pending

ACCEPTED WITH COMMENTS