

# MULTI-PURPOSE HOSPITAL DISINFECTANT

# SBT-24

**SBT-24 CLEANS AND DISINFECTS IN ONE OPERATION • KILLS A WIDE VARIETY OF GRAM POSITIVE AND GRAM NEGATIVE BACTERIAL PATHOGENS, SUCH AS STAPHYLOCOCCI, STREPTOCOCCI, COLIFORMS, PROTEUS AND SALMONELLA • DESTROYS MOST DERMATOPHYTES AND MILDEW ORGANISMS  
SBT-24 IS TUBERCULOCIDAL-IN-VIVO CHALLENGE CONFIRMATION TEST (GUINEA PIGS)**

## SBT-24 TREATED SURFACES

continue to combat new invasions of pathogens. One treatment aids in protecting against repeat exposure to bacteria. Treated surfaces become self-sanitizing for prolonged periods.

Long lasting residual action of SBT-24 provides protection against odors.

**USE SBT-24** — in addition to regular cleaning and disinfecting procedures, 2 ounces per gallon of water to fog isolation rooms, laundry chutes, operating rooms and nurseries. Set fogger on floor or table, vacate room until fogging is completed. SBT-24 may be nebulized through oxygen tents, without damage to rubber or metal parts.

Disinfects all hard surfaces: floors, dressing containers, walls, tables, operating room equipment.

Disinfects and cleans in one operation. ADD 2 ounces SBT-24 per gallon water or per gallon diluted detergent.

Disinfects instruments; catheters, rubber sheeting, tubing. IMMERSE ten minutes in solution 2 ounces SBT-24 per gallon water.

Active Ingredients: n-Propyl alcohol, Potassium ortho-benzyl parachlorophenate, 3,4',5 Tribromosalicylanilide and 4',5 Dibromosalicylanilide Ethylenediamine tetraacetic acid and Potassium soaps of tallow fatty acids

Inert Ingredients: 67.65% Water

Store away from heat

## KEEP OUT OF REACH OF CHILDREN

LEVER BROTHERS COMPANY, INDUSTRIAL DIVISION, NEW YORK 22, N.Y.

828-016 270

1 GALLON

U.S.D.A. Reg. No. 6018-7

ACCEPTED WITH COMMENTS