

SIGN OF THE PROFESSIONAL



GlideRinse III

Disinfectant - Sanitizer
Fungicide - Deodorizer
with Organic Soil Tolerance
for Hospital, Institutional, Industrial and School Use

Active Ingredients:

Allyl (C ₁₂ , 50%, C ₁₂ , 40%, C ₁₂ , 10%)	
Dimethyl benzyl ammonium chloride	3.000%
Octyl decyl dimethyl ammonium chloride	2.250%
Didecyl dimethyl ammonium chloride	1.125%
Dioctyl dimethyl ammonium chloride	1.125%
Inert Ingredients	92.500%
Total	100.0%

KEEP OUT OF REACH OF CHILDREN

DANGER

Statement of Practical Treatment

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse.

If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution, or if these are not available, drink large quantities of water. Avoid alcohol. Call a physician immediately.

NOTE TO PESTICIDE: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression, and convulsion may be needed.

SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS

E.P.A. Reg. No. 6109-18

E.P.A. Est. No. 6109-IL-1, 6109-IL-2, 6109-PA-1

FOR SALE TO, USE, AND STORAGE BY SERVICE PERSONNEL
NET FIVE U.S. GALLON (18.9 LITERS)

ACCEPTED

JUL 28 1982

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 6109-18

Precautionary Statements

DANGER Hazard to Humans and Domestic Animals
Keep out of reach of children. Corrosive. Causes eye damage and skin irritation. Do not get in eyes, on skin, or on clothing. Wear goggles or face shield and rubber gloves when handling this product. Harmful if swallowed. Do not breathe spray mist. Avoid contamination of food.

Directions for Use

GENERAL CLASSIFICATION

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
A-1: GlideRinse III with a cloth, mop or mechanical spray device. When applied with a mechanical spray device, surface must be sprayed until thoroughly wetted. Treated surfaces must remain wet for 10 minutes. Fresh solution should be prepared daily or when the use solution becomes visibly dirty.

Disinfection in Hospitals, Nursing Homes and Other Health Care Institutions: For disinfecting floors, walls, countertops, bathing areas, lavatories, restrooms, tables, chairs, garbage pails and other hard non-porous surfaces.

A-2: 1.0 oz. of GlideRinse III per gallon of water (8 ml/l). Apply to previously cleaned hard surfaces. At this use level, GlideRinse III is effective against *Pseudomonas aeruginosa*, *Staphylococcus aureus* and *Salmonella choleraesuis* in the presence of 5% blood serum when evaluated by 0% ADAC Use-Dilution test.

Disinfection in Institutions, Industry and Schools: Add 1/2 oz. of GlideRinse III per gallon of water (1 ml/l). At 0.5 oz. to 2.25 gallons of water (12.5-ml/US l) use level, the fungicidal effectiveness of GlideRinse III in the presence of 5% blood serum against *Trichophyton mentagrophytes* has been shown utilizing the ADAC Fungicidal Test.

STORAGE AND DISPOSAL

- Keep in cool area.
- Do not contaminate water, food, or feed by storage or disposal.
- Open sampling is prohibited.
- Do not reuse empty container.

PESTICIDE RESIDUE

Pesticides or residues that cannot be used or chemically reprocessed should be disposed of in a landfill approved for pesticides or buried in a safe place away from water supplies.

CONTAINER DISPOSAL

Input rinse (or equivalent) and dispose in an incinerator or landfill approved for pesticide containers, or bury in a safe place.

GENERAL

Consult Federal, state or local disposal authorities for approved alternative procedures such as incineration, burning.

GLIDERINSE III

GlideRinse III is a multiuse disinfectant, sanitizer, fungicide, and deodorant which is effective in the presence of organic soil. GlideRinse III is designed for use in hospitals, schools, and other institutions where hard surface disinfection can aid in controlling cross contamination.

Sanitizing of Food Processing Equipment and Other Hard Surfaces in Food Contact Locations: For sanitizing food processing equipment, dairy equipment, food utensils, dishes, silverware, glasses, sink tops, countertops, refrigerated storage and display equipment and other hard non-porous surfaces. No potable water rinse is required.

Wash and rinse all articles thoroughly, then apply a solution of 1 oz. GlideRinse III in 4 gallons of water (150 ppm active). Surfaces should remain wet for at least one minute followed by adequate draining and air drying. Fresh solution should be prepared daily or when use solution becomes visibly dirty. For mechanical application, use solution may not be reused for sanitizing applications.

Apply to sink tops, countertops, refrigerated storage and display equipment and other stationary hard surfaces by cloth or brush or mechanical spray device. No potable water rinse is required.

Dishes, silverware, glasses, cooking utensils and other similar size food processing equipment can be sanitized by immersion in a 1 oz./4 gallon solution of GlideRinse III. No potable water rinse is required.

At 1 oz./4 gallons GlideRinse III fulfills the criteria of Appendix F of the Grade "A" Pasteurized Milk Ordinances 1978 Recommendations of the U.S. Public Health Services in waters up to 800 ppm of hardness calculated as CaCO₃ when evaluated by the AOAC Germicidal and Detergent Sanitizer Method against *Escherichia coli* and *Saphylococcus aureus*.

TO REMEMBER:
23205

Manufactured by:
ServiceMaster Industries Inc.
Downers Grove, Illinois 60515, U.S.A.

SERVICEMASTER

SERVICEMASTER INDUSTRIES INC.

300 Warrenville Road
Downers Grove, Illinois 60515
312) 964-1300

BEST DOCUMENT AVAILABLE