

Label corrected
 in person 4/2/70
 superseded label
 Submitted 3/11/70

Under a proper condition of use, Sani-Det is highly effective for all the common bacteria, viruses, and fungi. Its effectiveness is particularly high in the presence of organic matter. Sani-Det is also effective against a wide range of other organisms, which tend to develop in the presence of organic matter. Sani-Det is particularly recommended for use on mattresses, plastic mattress covers, pillows, and room equipment which is being moved from room to room.

Sani-Det should be used EXHAUSTED in a well-ventilated area. The suggested spraying distance is between one and three feet from the surface. This distance will produce a noticeable condensation without excessive soaking. Sani-Det is highly effective against mold and mildew. It will remain so for only a few minutes necessary for absorption and evaporation of the solvent. Sani-Det is a sanitizer on porous surfaces such as mattresses, pillows, covers, etc., sprayed from a distance of one to three feet until moist to touch, avoid excessive soaking. This will dry in a few minutes without harmful effects to the fabric. Sani-Det is not to be used for terminal disinfection.

Cameraman's Note
 Poor Copy

(See Insert 1 Label)

(See Insert 2 Label)

CONTENTS

INGREDIENTS:

Insert 1

- 1. Red Water, Sodium Borate, Triethylene Glycol,
- 2. Ethylene Glycol, Mineral Oil, Acetylphenylpolyethoxy-
- 3. U. S. P. Gentian Color

TOTAL: 78.355%

[Handwritten signature]

TESTS:

Insert 2

- 1. ...
- 2. ...
- 3. ...

ACCEPTED WITH COMMENTS