

LONZA TBC 301 is a disinfectant bathroom bowl cleaner, designed for hospital use. This disinfectant toilet bowl cleaner is especially formulated to kill many germs commonly found in hospital, industrial and institutional environments. LONZA TBC 301 is designed for use in toilet bowls and urinals.

LONZA TBC 301 contains an acid compatible detergent system that removes resistant stains and increases its disinfectant action.

When used as directed, LONZA TBC 301 cleans, disinfects and deodorizes in one easy step.

Precautionary Statements

Hazardous to Humans and Domestic Animals

DANGER

Corrosive: Keep Out of Reach of Children. Causes eye and skin damage (or skin irritation); produces chemical burns. Do not get in eyes, on skin, or on clothing. Wear goggles or face shield and rubber gloves when handling. Do not breathe vapor or fumes. May be harmful or fatal if swallowed.

Storage and Disposal

- Do not contaminate water, food, or feed by storage disposal.
- Do not store on side.
- Avoid creasing or impacting of side walls.

PESTICIDE DISPOSAL

Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray or mixture of rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL

Triple rinse (or equivalent). Then offer for recycling or reconditioning or puncture and dispose of in a sanitary landfill, or incineration, or, if allowed by state and local authorities by burning. If burned, stay out of smoke.

Meta Container: Triple rinse (or equivalent), then offer for recycling or reconditioning, or dispose of in a sanitary landfill, or by other procedures approved by state and local authorities. (If container is 1 gallon or less, use this container disposal statement.)

Do not reuse empty container (bottle, can, bucket). Wrap container and put in trash.

LONZA TBC 301

High Acid Disinfectant Toilet Bowl Cleaner
With Organic Soil Tolerance
For Hospital-Institutional-Industrial Use

Based on a Synergistic Combination of
Quaternary Ammonium Compounds
And an Acid Compatible Detergent System
Disinfects-Cleans-Deodorizes
Removes Stains and Hard Water Deposits

Active Ingredients

Octyl decyl dimethyl ammonium chloride	0.2250%
Dioctyl dimethyl ammonium chloride	0.0900%
Didecyl dimethyl ammonium chloride	0.1350%
Alkyl (C ₁₄ , 50%; C ₁₂ , 40%; C ₁₆ , 10%) dimethyl benzyl ammonium chloride	0.3000%
Hydrogen Chloride	9.6000%

Inert Ingredients

89.6500%
100.0000%

TPD - 8.6916/gal.
KEEP OUT OF REACH OF CHILDREN

DANGER

STATEMENT OF PRACTICAL TREATMENT

If swallowed, call a physician immediately. Drink a teaspoonful or more of milk of magnesia or a cup of soapy water followed with milk or raw egg whites.

Eye Contact: Wash with water for 15 minutes. Get prompt medical attention.

Skin Contact: Wipe off the acid gently, immediately flood the surface with water, using soap freely, then cover with moist magnesia or baking soda.

SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS

EPA Registration No. 6836-85

EPA Establishment No. 6836-IL-1

Net Contents

MANUFACTURED BY:

LONZA INC., 17-17 Route 208, Fair Lawn, N.J. 07410

For Sale For Use And Storage By
Maintenance Personnel Only

Reg # 6836-85 PM-32

DIRECTIONS FOR USE 191

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Flush toilet. Add 3 fluid oz. (6 tablespoons) directly to water in toilet bowl (¼ gallon). Brush to cover all interior surfaces of the bowl including under the rim. Allow disinfectant solution to remain at least 10 minutes. Flush toilet.


For heavily soiled areas, a preliminary cleaning may be required.

To clean badly stained bowls, repeat above treatment using 5 to 6 fluid ounces of LONZA TBC 301.

For use in urinals, add 3 fluid oz. (6 tablespoons) to ¼ gallon of water. Apply by brush to all urinal surfaces. Treated surfaces must remain wet for 10 minutes. Flush.

The broad spectrum efficacy of LONZA TBC 301 in toilet bowls and urinals has been evaluated in the presence of 5% serum and found to be effective against the following organisms by the AOAC Use-Dilution test:

Pseudomonas aeruginosa
Staphylococcus aureus
Salmonella choleraesuis
Escherichia coli
Salmonella typhosa
Serratia marcescens
Shigella dysenteria
Streptococcus faecalis



BEST AVAILABLE COPY