


# Chem-Spray



# STAF HOSPITAL SPRAY DISINFECTANT

### DIRECTIONS FOR USE

For use in Hospitals, Nursing Homes, Institutions, Offices, Schools, Motels and Hotels.

Kills Staphylococcus aureus, Salmonella choleraesuis and Pseudomonas aeruginosa.

For use on telephones, door knobs, night switches, bed springs, diaper pails, toilet seats, sink tops, bed pans, garbage pails, hampers, linen carts and desk tops.

To disinfect, hold can six (6) inches from surfaces and spray until completely wet.

To deodorize and air refreshen, point spray upwards to center of room and spray for three (3) seconds.

**WARNING:** Avoid getting spray concentrate in eyes. In case of contact, flush with water. If irritation persists, get medical attention.

**WARNING:** Flammable! Contents under pressure. Do not puncture. Do not use or store near heat or open flame. Exposure to temperatures above 130°F. may cause bursting. Never throw container into fire or incinerator.


### Active Ingredients:

n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides	0.1%
n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride	0.1%
Isopropanol	53.0%
Essential Oils	0.5%
Inert Ingredients:	46.3%
	100.0%

Contains Sodium Nitrite as a corrosion inhibitor

E.P.A. Reg. No. 7056-67

E.P.A. EST. 7056 IX 1



**ELIMINATES ODORS  
KILLS MOST GERMS**

**WARNING: KEEP OUT OF REACH OF CHILDREN.  
CONTENTS UNDER PRESSURE.  
SEE BACK PANEL FOR ADDITIONAL CAUTIONS.**

## NET WT. 1 lb. 2 oz. AVDP

Manufactured by  
**CHEM SPRAY AEROSOLS, INC.**  
HOUSTON, TEXAS 77088